

NL4000 BodyTom[®] Elite User Manual

1-NL4000-060 Revision 19



NeuroLogica Corporation

14 Electronics Avenue, Danvers, MA 01923 USA

Under the copyright laws, the documentation may not be copied, photocopied, reproduced, translated, or reduced to any electronic medium or machine-readable form in whole or in part without the prior written consent of NeuroLogica Corporation.

Copyright © applied in 2022

Printed in 2017

Contents

Intended use of the system	24
Clinical benefit	24
Consumer information	24
Proprietary rights	24
Legal disclaimer	25
Contact information	25
Damage in transportation	26
User requirements	26
Essential Performance	27
About this user manual	27
Identified symbols and system classifications	28
Conventions used in this user manual	28
Applicable versions of BodyTom Elite	30
 Chapter 1 Compliance and Safety Requirements	31
IEC classification and symbols	31
Environmental specifications	34
Considerations when preparing gantry for use	35
Site specification	36
Site specification for enclosed CT room	37
Hazardous substances	38
Part numbers and product-marking plates	38
Class 1 Type B medical devices	41
Focal spot	42
Anode input power	42
Filtration	43
Source to Detector distance (SID)	43
Compliance statement	43
EMI/EMC terms	44
Electromagnetic Compatibility	44
Electromagnetic Interference (EMI)	45
Susceptibility	45
Countermeasures against EMC related issues	50
Use recommendations	54
Installation recommendations	55
Cable shielding and grounding	55
Adjacent components and equipment	55

Static magnetic field limits	55
Electrostatic discharge environment and recommendations	55
Facility IT-NETWORK.....	56
Hazard information	56
Laser safety	58
Scanner mobility safety	60
Electrical safety	61
Mechanical safety	63
Radiation safety	64
Fire and explosion safety.....	66
EMERGENCY STOP button.....	66
Battery safety and information.....	67
Maintenance and service	71
Cybersecurity.....	71
Contraindication(s).....	72
Personnel privileges and terminology	72
Qualified operator.....	72
Operator of record	72
Scanning privileges	72
Protocol privileges	72
Administrative privileges	73
Clinical operation.....	73
Clinical scanning.....	73
Clinical protocol	73
Kernel	73
Chapter 2 System Overview	74
BodyTom Elite system	74
Overview of the scanner control panels.....	75
Controls on the left end of the scanner.....	76
Controls on the right end of the scanner	78
Identifying operator control panel buttons.....	80
Overview of the pendant.....	81
Overview of the workstation	84
Workstation types	85
Workstation remote power controls	86
Workstation UPS	87
Workstation considerations before use	87
Keyboard and mouse	88

Workstation power cord.....	88
Parts that potentially come into contact with the patient	88
Chapter 3 Basic Scanner Operations	90
Powering on and off the BodyTom Elite system.....	90
Checking a connection between the workstation and the scanner	91
Moving and transporting the scanner	92
Drive direction of scanner	95
Safety bumper system.....	95
Positioning the scanner before a scan.....	96
Positioning the patient using the laser lights	97
Operating the E-STOP button	100
Restoring the system from E-STOP	100
Chapter 4 Basic Workstation Operations	101
Understanding the types of users	101
Using the workstation	101
Identifying the workstation's remote power display	101
Identifying the microphone, speaker, and controls	103
Powering the workstation.....	104
Logging in to the workstation	105
Locking and unlocking the workstation.....	108
Navigating around the workstation's main screen.....	110
Brief overview of the main menu.....	110
Brief overview of the File menu	110
Brief overview of the Tools menu	112
Brief overview of the Customize menu	113
Getting Help from the Help menu.....	114
Getting to know the status bar.....	116
The workstation tabs.....	119
Chapter 5 System and User Configuration and Setup.....	121
Setting user accounts	122
Setting or updating the user's information	124
Unlocking the user	125
Deleting the user.....	126
Modifying the order of the user in the accounts list.....	127
Assigning general settings	129

Managing DICOM servers	136
Assigning a server as a store or worklist server	137
Modifying a server	140
Echoing a server	141
Deleting a server	141
Moving a server up and down the server list	142
Saving DICOM servers to a PACS list	143
Selecting PACS options	144
Assigning DICOM settings	147
Assigning audio configuration	152
Finding and listening to audio files	152
Recording and saving an audio file	153
Transmitting an audio file	155
Deleting an audio file	155
Assigning dose report	156
Applying dose configuration	157
Setting Dose Check	158
Assigning Dose Configuration to a patient protocol	159
Updating saved dose	162
Deleting a saved dose limit	164
Applying Windowing Presets	165
Editing kernel presets	165
Setting Window Presets	168
Editing a window preset	170
Deleting a preset	172
Setting up the Audit Trail Viewer	172
Setting image orientation	175
Setting Filter Kernels	178
User configuration	180
Updating your user account	180
Selecting a room for the BodyTom Elite	188
Chapter 6 Protocol Manager	189
Creating a new protocol	190
Using Build From to create a new protocol	199
Editing an Existing Protocol	203
Copying and pasting protocols	208
Deleting a protocol	210
Adding breathing instructions to your protocol	212

Importing protocols from a storage device	215
Exporting protocols to a storage device	218
Changing the order of protocols in the list	220
 Chapter 7 Daily Calibration and Quality Assurance	222
The QA phantom overview	225
Starting Quality Assurance	227
Ensuring good image quality	231
Identifying filtration accuracy	232
Using Axial plane to determine image resolution	233
Using slice width	234
Measuring noise	234
Measuring low contrast	234
Finding uniformity	234
Identifying CT contrast scale	235
Using performance phantoms	235
Measuring dose	235
The BodyTom Elite dose information (21 CFR 1020.33 c)	236
Body CTDI _w phantom	238
Head CTDI _w phantom	238
The BodyTom Elite dose in air	239
QA measurements	239
Tube accuracy	239
ACR testing procedure	239
Identifying high-contrast resolution	241
Noise, uniformity, and mean CT number of water	241
Uniformity and mean CT number	242
Geometric efficiency in the Z axis direction notification	243
Half-value layer	244
Allowable variations	244
Scatter radiation	244
Dose profile/Geometric Efficiency	250
 Chapter 8 Patient Registration	252
Navigating the Patient Registration screen	252
Registering the patient	253
Querying patient information	253
Storing patients in the Stored Results list	256
Manually registering a patient	257

Viewing patient information	260
Deleting patients from the Stored Result list	260
Chapter 9 Patient Scanning	262
Identifying Scan Types	263
Performing a scan	264
Repeating an image	274
Extending a scan	275
Scanning with special features	277
Using the step-and-shoot option	277
Performing a scan with Automatic Exposure Control	279
Performing a CT angiography scan with Bolus Tracking	282
Performing Test Bolus	286
Performing a CT Perfusion Scan	289
Calculating and creating perfusion maps	291
Using the Interventional Package	293
Examining the scanned image with tools	299
Using tools on the Acquisition tab	299
Chapter 10 Patient Browser	301
Navigating the Patient Browser	302
Identifying symbols on Patient Browser	302
Using the vertical and horizontal scroll bars on Patient Browser	303
Locking a study	303
Marking a series to read	304
Using the preview window	305
Archiving patient series	306
Archiving to PACS	306
Archiving to Media	308
Archiving to Navigation	311
Import	312
Importing from PACS	312
Importing from media	314
Delete	315
Registering a patient from Patient Browser	316
Building dose from Patient Browser	317
Merge	318
Using Show Info to view, update, and move a series	321
Modifying a series scanned under the wrong patient	323

Loading a series into view	325
Appending a series	325
Chapter 11 Viewing Images	327
Using keyboard shortcuts	329
Setting window width and center	330
Viewing images in 2D	331
Comparing images	332
Comparing a scout and a scan	334
Using the ROI	335
Using Layout and Rotate in 2D view.....	336
Applying a grid to your images in 2D.....	337
Viewing images in MPR	338
Understanding and using slab	339
Creating the slab	339
Viewing images in 3D	344
Chapter 12 Post Reconstruction	346
Reconstruction Overview	346
Metal artifact reduction	346
Noise reduction	346
Windmill Correction	346
Performing Post Reconstruction	348
Resending images from the scanner to the workstation	352
Chapter 13 Accessories and Options.....	353
Using the Universal Transfer Board.....	353
Chapter 14 Cleaning and Storing the System and Workstation Specifications	357
Cleaning the scanner and workstation	357
Cleaning the outside of the scanner and workstation.....	358
Maintenance of the workstation	359
Storing the system.....	359
Storing the scanner and workstation.....	359
Storing the QA phantom	360
Workstation specifications	360
Understanding the symbols and product-marking plate	361
Locating the product-marking plate on the workstation	363

Listing of replacement parts for workstation	363
Product Safety and Electromagnetic Comparability	364
Appendix A Glossary	365
Appendix B Listing of All Buttons, Tools, and Icons	375
Status bar icons	375
System state orbs	377
Workstation buttons	378
Viewing tools	383
Appendix C Sample of Reference Protocols Provided	389
Appendix D Automatic Exposure Control	390
1 Introduction:.....	390
2 Image Noise:.....	390
3 AEC working instructions:	392
3.1 AEC input parameters:	392
3.2 The scan parameters:	393
3.3 Notes	393
3.4 Sample protocols:	394
4 AEC algorithm description:.....	394
5 The Noise measurements:.....	397
Appendix E Rotating Anode X-Ray Tube	398
Appendix F Error Code	406
Appendix G Revision History	421

List of Figures

Figure 1: Current cover (left) and former covers (middle & right).....	30
Figure 2: Product-marking plate on scanner	38
Figure 3: Scanner dimensions including drivebar	40
Figure 4: Identifying the scanner's safety label(s) - foot-crush-hazard label(s)	58
Figure 5: Laser aperture's direction	59
Figure 6: Dangerous-to-patient/operator safety-warning label location (left) and label (close-up, right).....	65
Figure 7: BodyTom Elite E-STOP locations (right and left)	67
Figure 8: Close-up of the scanner control panel and the E-STOP button.....	67
Figure 9: Scanner battery capacity icon	68
Figure 10: Workstation battery capacity icon	69
Figure 11: Scanner X-ray tube capacity icon.....	70
Figure 12: BodyTom Elite system configuration	75
Figure 13: Left end of the scanner	76
Figure 14: Right end of the scanner.....	78
Figure 15: Operator control panel buttons and indicators	80
Figure 16: BodyTom Elite remote-control pendant	83
Figure 17: Scanner's positional display	83
Figure 18: Workstation with leaded-glass shield (optional installation).....	84
Figure 19: Workstation safe distance location (two views).....	86
Figure 20: The remote power display.....	87
Figure 21: Workstation keyboard and mouse	88
Figure 22: AC cord and storage on scanner (120V left plug in or 240V right plug out).....	90
Figure 23: Scanner's power cord receptacle for 120VAC and 240/VAC.....	91
Figure 24: Scanner hardwired to the workstation with an ethernet cable to data-access ports	91
Figure 25: Transport button on the operator control panel	92
Figure 26: Rocker-Switch-Lift UP (top) and DOWN (bottom) button.....	93
Figure 27: Drive bar front.....	93
Figure 28: Drive bar.....	93
Figure 29: Scanner drive direction (right side view).....	95
Figure 30: Bumper system	95
Figure 31: Patient centered in bore (height positioning)	96
Figure 32: Phantom positioned in center of FOV.....	98
Figure 33: Pendant use for positioning lasers upon patient	98
Figure 34: Positioning lasers upon patient	99
Figure 35: Positional display	99
Figure 36: BodyTom Elite E-STOP locations (right and left).....	100
Figure 37: E-STOP button on the scanner control panel on both the left and right sides of the scanner	100
Figure 38: Workstation remote power display.....	102

Figure 39: Microphone, speaker, and controls	103
Figure 40: Remote power display on workstation	104
Figure 41: Dell workstation power on button	105
Figure 42: Shutdown computer popup	105
Figure 43: User ID dropdown box	106
Figure 44: User ID dropdown list	106
Figure 45: Password text box	106
Figure 46: Patient Registration tab	107
Figure 47: User ID, current date, and time	107
Figure 48: System Lock button	108
Figure 49: Lock/Unlock System popup to lock the workstation	108
Figure 50: Unlock button	109
Figure 51: Lock/Unlock System popup to unlock the workstation	109
Figure 52: Main menu	110
Figure 53: File menu	111
Figure 54: File > Log Off	111
Figure 55: Login popup	111
Figure 56: File dropdown menu	111
Figure 57: Restart Application or Restart Computer popup	112
Figure 58: File dropdown menu	112
Figure 59: Tools dropdown menu	113
Figure 60: Customize dropdown menu	113
Figure 61: Help dropdown menu	114
Figure 62: Support Connection browser window	114
Figure 63: About Us popup	115
Figure 64: Scanner and workstation status bar	116
Figure 65: Workstation tabs to perform a patient examination	119
Figure 66: User Accounts tab	122
Figure 67: User account fields filled in	123
Figure 68: Save aborted popup message - Password requirements	124
Figure 69: List of users	124
Figure 70: Update Aborted popup message - Password requirements	125
Figure 71: List of users not selected	126
Figure 72: List of all available users	127
Figure 73: Down arrow	128
Figure 74: Up arrow	128
Figure 75: Save button for list order	129
Figure 76: General Settings tab	130
Figure 77: General Settings > Hospital Setup subtab	131
Figure 78: General Settings > Application Setup subtab	132

Figure 79: General Settings > Scanner Setup subtab	134
Figure 80: General Settings > Remote Support Setup subtab	135
Figure 81: DICOM Servers tab	136
Figure 82: DICOM Servers tabs	137
Figure 83: DICOM Servers > Servers tabs	138
Figure 84: Action Succeeded popup message - Server saved.....	140
Figure 85: Action Succeeded popup message - Server updated	140
Figure 86: Echo Successful and Echo Failed popup	141
Figure 87: Up and Down arrows to move up and down server list.....	142
Figure 88: Save button	143
Figure 89: DICOM Servers > PACS List tab	144
Figure 90: PACS List Saved popup message - PACS saved.....	144
Figure 91: DICOM Servers > Options tab	145
Figure 92: Time (increase and decrease time) arrows	146
Figure 93: PACS List Saved popup	147
Figure 94: DICOM Settings tabs (six)	148
Figure 95: DICOM Settings > HIS/RIS Query	149
Figure 96: DICOM Settings > MPPS	149
Figure 97: DICOM Settings > Patient Module	150
Figure 98: DICOM Settings > Study Module	150
Figure 99: DICOM Settings > Series Module.....	151
Figure 100: DICOM Settings > Image Module	151
Figure 101: Audio Configuration tab	152
Figure 102: Audio files list.....	153
Figure 103: Audio files list.....	154
Figure 104: New audio file.....	154
Figure 105: Audio files transmitted to save to the scanner	155
Figure 106: Dose Report tab.....	156
Figure 107: Generated dose report.....	157
Figure 108: Dose Configuration > Dose Check	158
Figure 109: Save Successful popup - Dose Check successfully saved.....	159
Figure 110: Dose Configuration > Dose Configuration for adult and pediatric	160
Figure 111: Anatomical orbs.....	160
Figure 112: Pediatric Dose Configuration Parameters.....	161
Figure 113: Invalid Parameter popup message - Dose setting kV already exists	161
Figure 114: Save Successful popup message - Maximum dose saved.....	162
Figure 115: Saved Doses List.....	163
Figure 116: Save Successful popup message - Maximum dose saved.....	163
Figure 117: Dose Configuration > Dose Check tab	164
Figure 118: Save popup message - Maximum dose saved	165

Figure 119: Windowing Preset tab	166
Figure 120: Windowing Presets > Kernel Presets tab	166
Figure 121: Sharpness dropdown	167
Figure 122: Action Succeeded popup message - Preset saved.....	168
Figure 123: Window Presets tab.....	168
Figure 124: Window Presets > Name.....	169
Figure 125: Window Presets > Width	169
Figure 126: Window Presets > Center	170
Figure 127: Action Succeeded popup message - Preset saved.....	170
Figure 128: Listing update	171
Figure 129: Action Succeeded popup message - Preset saved.....	171
Figure 130: Action Succeeded popup message - Preset deleted	172
Figure 131: Audit Trail Viewer tab.....	173
Figure 132: Adding a date or a date span.....	173
Figure 133: Audit Trail Viewer > Audit Type dropdown	174
Figure 134: Audit Trail Viewer > User ID dropdown	174
Figure 135: Audit results	175
Figure 136: Image Orientation tab.....	176
Figure 137: Image Orientation > New Flip Orientation dropdown	177
Figure 138: Settings Saved popup message - Image orientation settings saved	177
Figure 139: Filter Kernels tab	178
Figure 140: Selected Axial kernel	179
Figure 141: Selected Helical kernel	179
Figure 142: Last Name, First Name, Password, and Verify Password fields	180
Figure 143: Update Succeeded popup message - Account updated	181
Figure 144: Column Settings dialog box with HIS/RIS Query option	182
Figure 145: Column Settings with a selected query in HIS/RIS	182
Figure 146: Column Settings with HIS/RIS Query option using Up and Dwn buttons	183
Figure 147: Make Default option	183
Figure 148: Column Settings with Patient Browser option	184
Figure 149: Column Settings with Patient Browser Series option - using Up and Dwn buttons.....	185
Figure 150: Make Default option	185
Figure 151: Scan Dosage Report tab	186
Figure 152: Date, Protocol, and mA Range filled.....	187
Figure 153: Scan Dosage Report results.....	187
Figure 154: Available rooms before moving the scanner	188
Figure 155: Protocol Manager for adult and pediatric	190
Figure 156: Adult and pediatric anatomical orbs, with Chest orb selected	191
Figure 157: Adult and pediatric protocol lists	192
Figure 158: New Protocol dialog box	192

Figure 159: Patient position handles	193
Figure 160: New Series dialog box.....	194
Figure 161: New Reconstruction popup.....	197
Figure 162: Edit Series dialog box	198
Figure 163: Save New Protocol.....	199
Figure 164: Close Button	199
Figure 165: Anatomical orbs.....	200
Figure 166: Build from protocol selected.....	201
Figure 167: Build From button	201
Figure 168: New Protocol dialog box.....	202
Figure 169: Build from save.....	203
Figure 170: Build from close.....	203
Figure 171: Edit protocol orbs.....	204
Figure 172: Edit protocol selected	205
Figure 173: Edit button.....	205
Figure 174: Edit Protocol dialog box	206
Figure 175: Edit series update button.....	206
Figure 176: Edit protocol update button	207
Figure 177: Edit protocol close button.....	207
Figure 178: Protocol Manager for Adult and Pediatric	208
Figure 179: Anatomical orbs, in this case the chest orb	209
Figure 180: Copy right-click floating menu	209
Figure 181: Paste right click floating menu.....	210
Figure 182: Protocol Manager with a protocol selected	211
Figure 183: Delete Confirmation popup message - Yes or No to delete selection	211
Figure 184: Edit button.....	212
Figure 185: Edit Protocol dialog box	213
Figure 186: Add breathing edit button.....	213
Figure 187: Edit Series dialog box	214
Figure 188: Use Breathe Indicator Audio option	214
Figure 189: Breathe Indicator Audio Files popup.....	215
Figure 190: Import button	216
Figure 191: Select File popup	216
Figure 192: Select file	217
Figure 193: Import button active in Select File when file(s) selected.....	217
Figure 194: Protocols Imported popup message - Protocols imported.....	218
Figure 195: Export button.....	218
Figure 196: Select Directory popup	219
Figure 197: Protocols Exported popup message - Protocols exported	219
Figure 198: Changing protocol order with Up and Down (arrow) buttons.....	220

Figure 199: Protocol Save button	221
Figure 200: Perform Daily Cal popup	223
Figure 201: Perform Daily Cal popup with count down	224
Figure 202: Perform Daily Cal summary popup	224
Figure 203: Air freshness icon changes as the air quality drops from green to yellow to red	225
Figure 204: QA phantom	225
Figure 205: Phantom holder	226
Figure 206: Phantom on the phantom holder	228
Figure 207: Place QA phantom	228
Figure 208: Proper QA stand positioning	228
Figure 209: QA phantom positioning	229
Figure 210: Laser button	229
Figure 211: Quality Assurance popup	230
Figure 212: System Ready to Scan popup message - System is ready to begin	230
Figure 213: QA results of QA image	230
Figure 214: Phantom image	231
Figure 215: Locked QA results shown in Patient Browser	231
Figure 216: Results of QA image after the QA test	232
Figure 217: MTF	241
Figure 218: Catphan 515 using 120kV, 300mA, 1 rotation, and 5mm slice	243
Figure 219: Geometric efficiency confirmation	243
Figure 220: Scatter plot (120kV, 100mA in μ Rad)	246
Figure 221: Scatter plot (120kV, 100mA in μ Gy)	247
Figure 222: Scatter plot (140kV, 300mA in μ Rad)	247
Figure 223: Scatter plot (140kV, 300mA in μ Gy)	248
Figure 224: Scatter measurements (X–Y plane)	248
Figure 225: Scatter measurements (Y–Z plane)	249
Figure 226: Dose profile for 8 rows	250
Figure 227: Dose profile for 32 rows	251
Figure 228: Activated Patient Registration tab	252
Figure 229: Patient Registration tab	253
Figure 230: Query Information dialog box	254
Figure 231: Edit Value popup for name	255
Figure 232: Patient Registration Query Results table	255
Figure 233: Patient Registration Stored Results table	256
Figure 234: Patient Registration tab	257
Figure 235: Exam Information dialog box	257
Figure 236: Patient ID field	258
Figure 237: Patient data filled in	259
Figure 238: Expand link in context and close up	259

Figure 239: Exam Information popup	259
Figure 240: Patient Registration tab	260
Figure 241: Patient Registration tab	260
Figure 242: Active Acquisition tab.....	262
Figure 243: What appears on Acquisition.....	263
Figure 244: Exam Planner for Adult and Pediatric	265
Figure 245: Anatomical orbs, with the Chest orb selected	266
Figure 246: Protocol selected and Edit button active	266
Figure 247: Edit Protocol dialog box	267
Figure 248: Edit Series dialog box	268
Figure 249: Update button.....	268
Figure 250: Is Scanner Properly Positioned? popup	269
Figure 251: System Ready to Scan.....	269
Figure 252: Scanner control panel - START button	270
Figure 253: Countdown popup	270
Figure 254: Scanner control panel - CANCEL button.....	270
Figure 255: Scouts and FOV button.....	271
Figure 256: Continue button	272
Figure 257: Pending Scanner Movement popup message	272
Figure 258: System Ready to Scan popup message - System is ready to begin scan.....	272
Figure 259: Scanner control panel - START button	273
Figure 260: Perform Reconstructions popup message - To perform post reconstructions.....	273
Figure 261: Dose report	274
Figure 262: Protocol Viewer - Start Position and End Position.....	274
Figure 263: Repeat Protocol popup	275
Figure 264: Protocol Viewer dialog box.....	276
Figure 265: Extend Protocol popup.....	276
Figure 266: Extend Protocol popup.....	277
Figure 267: Step & Shoot option in the Edit Series dialog box.....	278
Figure 268: System Ready to Scan popup	278
Figure 269: Step & Shoot popup	279
Figure 270: Edit Series dialog box with AEC options selected	280
Figure 271: Toggle Graph button	281
Figure 272: Graphs on the scout(s).....	281
Figure 273: AEC modulation graph	282
Figure 274: Minimum mA and maximum mA; noise level	282
Figure 275: Bolus Tracking option	283
Figure 276: Active Scan Region - Bolus Reference or Helical CTA.....	284
Figure 277: Scout line (blue) and Reference line (green)	284
Figure 278: Bolus ROI tool	285

Figure 279: ROI on the Reference scan	285
Figure 280: Scan triggers when bolus enters reference point/ROI.....	286
Figure 281: Scan at peak enhancement.....	286
Figure 282: Test Bolus option.....	287
Figure 283: Active scan region	287
Figure 284: Scan planning lines	288
Figure 285: Bolus ROI	288
Figure 286: ROI on the Reference scan	288
Figure 287: Bolus timing graph	289
Figure 288: Edit Series CTP Scan Time	290
Figure 289: Brain Perfusion Image	291
Figure 290: CTP tools	291
Figure 291: Perfusion maps	292
Figure 292: Arterial Venous Flow	293
Figure 293: Interventional Tab	293
Figure 294: Interventional Tab - Patient exam details	294
Figure 295: Scan Tree, Current Scan, and Updated Dose Gauge	295
Figure 296: Viewing Tools, Windowing, and Zoom options	295
Figure 297: Interventional Workflow	296
Figure 298: Interventional Workflow - Protocol Information dialogue box.....	296
Figure 299: Interventional Workflow - Scan acquired.....	297
Figure 300: Modify protocol parameters.....	297
Figure 301: Initiate Scans - Interventional protocol	298
Figure 302: Move the Scanner	298
Figure 303: Repeat Scans - Interventional protocol.....	298
Figure 304: Active Patient Browser tab	301
Figure 305: Patient Browser sections.....	302
Figure 306: Patient Browser locked, read, PACS and Stored (archived), and media symbols	303
Figure 307: Patient Browser horizontal and vertical scroll bars	303
Figure 308: Floating menu - Lock	304
Figure 309: A locked series.....	304
Figure 310: Floating menu - Mark.....	305
Figure 311: Preview Button	305
Figure 312: The series appears in the preview window	305
Figure 313: Archive Destination popup.....	306
Figure 314: Archive To Server popup.....	306
Figure 315: Store/Print Queue dialog box	307
Figure 316: Archive Destination popup.....	309
Figure 317: Archive to Media popup	309
Figure 318: Archive Button active	310

Figure 319: Archiving complete.....	311
Figure 320: Archive to Server popup	312
Figure 321: Import Location popup.....	312
Figure 322: Import from PACS dialog box	313
Figure 323: PACS Query Information dialog box.....	313
Figure 324: Import PACS dialog box with active Get Series button	314
Figure 325: Import Location popup.....	314
Figure 326: Import from Media popup	315
Figure 327: Active Import button.....	315
Figure 328: Confirm Deletion popup	316
Figure 329: Patient browser register button.....	316
Figure 330: Create New Study popup	316
Figure 331: Build dose button.....	317
Figure 332: Please Wait popup.....	317
Figure 333: Dose Build Failed popup	318
Figure 334: Two series selected to merge	318
Figure 335: Active Merge button.....	319
Figure 336: Image viewer without image(s) selected	319
Figure 337: Selecting images to move to image viewer or bottom filmstrip	320
Figure 338: Finished filmstrip	320
Figure 339: Merge Series Saved popup.....	321
Figure 340: View/Update Information dialog box.....	321
Figure 341: View/Update Information dialog box.....	323
Figure 342: Move Series popup.....	323
Figure 343: Floating menu - Append Images	325
Figure 344: Please Wait popup.....	326
Figure 345: (Appended) series created	326
Figure 346: Active Viewing tab	327
Figure 347: Windowing preset dropdown list	330
Figure 348: Windowing Width and Center text boxes, and the Apply button	331
Figure 349: Right click menu	331
Figure 350: 2D tools	332
Figure 351: Floating menu - Mark For Compare	332
Figure 352: Floating menu - Compare with Selected Series.....	333
Figure 353: Compared series	333
Figure 354: Link button	333
Figure 355: Using the Compare button	334
Figure 356: Scout and scan selected to compare	334
Figure 357: Comparing a scout (1) and a scan (2).....	335
Figure 358: ROI	336

Figure 359: Layout (viewing tools).....	336
Figure 360: Rotate dropdown.....	337
Figure 361: Rotate (viewing tools).....	337
Figure 362: Grid (mm)	337
Figure 363: MPR tools	338
Figure 364: Image reformat selections.....	338
Figure 365: Tilt tool.....	339
Figure 366: Image formats	340
Figure 367: Enable Slab option	340
Figure 368: Enable Slab option under Secondary Series.....	340
Figure 369: Cyan Line and cyan circle to drag for FOV	340
Figure 370: Slab Thickness and Slab Spacing text boxes	341
Figure 371: Small yellow boxes to manually drag for desired slab thickness	341
Figure 372: Slab Rendering Options dropdown.....	341
Figure 373: Tilt tool.....	342
Figure 374: Tilt white circle	343
Figure 375: Capture Complete popup message - Series saved.....	343
Figure 376: MPR images in Patient Browser	343
Figure 377: 3D tools.....	344
Figure 378: 3D Color Presets.....	344
Figure 379: 3D Render modes.....	345
Figure 380: 3D Orientation options	345
Figure 379: Active Post Reconstruction tab	347
Figure 380: Post Reconstruction areas	348
Figure 381: Post Reconstruction study and series tables.....	348
Figure 382: Post Reconstruction viewing image area	349
Figure 383: FOV resizing boxes	349
Figure 384: Sharpness on the Reconstruction Parameters tab	350
Figure 385: Reconstruction Parameters Sharpness dropdown	350
Figure 386: Reconstruction Parameters Slice Thickness/Spacing dropdown	350
Figure 387: # of Expected Images.....	350
Figure 388: Noise Reduction on the Options tab for a Helical scan	351
Figure 389: Perform Windmill Correction and/or Noise Reduction on the Options tab for a Helical scan	351
Figure 390: Metal artifact removal	351
Figure 391: Please wait while the system performs data reconstruction message.....	351
Figure 392: Resend button	352
Figure 393: Universal transfer board and stiffeners	354
Figure 394: Four types of mattress stiffeners.....	355
Figure 395: Mattress stiffener in place.....	355
Figure 396: Universal transfer board properly positioned on the bed on a mattress stiffener.....	356

Figure 397: Universal transfer board with safety strap installed 356

Figure 398: BodyTom Elite castor wheels 359

Figure 399: Product-marking plate on side of the workstation 363

List of Tables

Table 1: Conventions used in this user manual	28
Table 2: Applicable symbols	32
Table 3: Operating environment	34
Table 4: System operating parameters and specifications	34
Table 5: System operating parameters.....	36
Table 6: Battery operating parameters	37
Table 7: Site specification	37
Table 8: Hazardous substances	38
Table 9: Core-system-component part numbers and product-marking plate locations	38
Table 10: Core-system component dimensions	40
Table 11: Workstation dimensions	40
Table 12: Filtration.....	43
Table 13: Acronyms and abbreviations.....	46
Table 14: Emission declaration for BodyTom Elite systems.....	47
Table 15: EMC Immunity declaration for BodyTom Elite systems.....	48
Table 16: Recommended separation distances	51
Table 17: Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment	52
Table 18: Left end of the scanner	76
Table 19: Right-end of the scanner.....	78
Table 20: Operator control panel buttons and indicators.....	80
Table 21: Pendant buttons	81
Table 22: Workstation power-control buttons	102
Table 23: Speaker control buttons.....	103
Table 24: Status bar icons.....	116
Table 25: System state orbs	118
Table 26: System configuration tabs	121
Table 27: Protocol Manager command buttons	189
Table 28: Scan protocols used by the QA	227
Table 29: Modulation Transfer Function (MTF) direction	233
Table 30: QA results.....	233
Table 31: Load factors	235
Table 32: Body CTDI _w (mGy/100mAs)	238
Table 33: Head CTDI _w (mGy/100mAs)	238
Table 34: Normalized CTDI of body phantom	238
Table 35: Normalized head CTDI	238
Table 36: CTDI air (mGy/100mAs)	239
Table 37: Mean and standard deviation of CTDI air	239
Table 38: Tube accuracy	239

Table 39: The CT number and linearity of the different inserts in the ACR phantoms	240
Table 40: The NeuroLogica head and abdomen ACR scan protocols	240
Table 41: The cutoffs of some of the common reconstruction kernels.....	241
Table 42: Uniformity and Mean CT Number using Water Phantom.....	242
Table 43: Uniformity and Mean CT Numbers using Water Phantom.....	242
Table 44: Half-value layer	244
Table 45: Scatter measurements (X–Y plane) (μ Rad/100 mAs)	249
Table 46: Scatter measurements (Y–Z plane) (μ Rad/100 mAs)	249
Table 47: The geometric efficiency of the two different collimations of the BodyTom Elite	251
Table 48: Patient Registration buttons	252
Table 49: Acquisition buttons	262
Table 50: Bolus tracking parameters and tools.....	282
Table 51 CTP Tools.....	292
Table 52: Image tools	299
Table 53: Command buttons	301
Table 54: Store and Print Queue buttons	307
Table 55: 2D, MPR, 3D, and CTP image tools.....	327
Table 56: Arrow key navigation.....	330
Table 57: Reconstruction tools.....	347
Table 58: Universal Transfer Board weight-bearing restrictions.....	354
Table 59: Workstation specifications	360
Table 60: Symbols and product-marking plate	361
Table 61: Status bar icons	375
Table 62: System state orbs.....	377
Table 63: BodyTom Elite workstation buttons.....	378
Table 64: Viewing tools	383
Table 65: Pendant buttons	387
Table 66: Sample of BodyTom Elite adult protocols and important estimates.....	389
Table 67: Sample of BodyTom Elite pediatric protocols and important estimates	389
Table 68: The measured noise at 120 kV.....	397
Table 69: Error code list	406
Table 70: Workstation (only) revision history before combined with BodyTom Elite User Manual (below).....	421
Table 71: Revision history	422

Intended use of the system

The NL4000 BodyTom Elite CT system is intended to be used for x-ray computed tomography applications for anatomy that can be imaged in the 85cm aperture. The CT system is intended to be used for both pediatric and adult imaging and as such has preset dose settings based upon weight and age.

The CT images can be obtained either with or without contrast. BodyTom Elite CT systems can be used for low dose lung cancer screening. The screening must be performed in compliance with the approved and established protocols as defined by professional medical societies.

* Please refer to clinical literature, including the results of the National Lung Screening Trial (N Engl J Med 2011; 365:395-409) and subsequent literature for further information.

Clinical benefit

Computed Tomography (CT) provides real time imaging of bone, soft tissue and blood vessels that can provide detailed information to diagnose, plan treatment for, and evaluate many conditions in adults and children. Additionally, the detailed images provided by CT scans may eliminate the need for exploratory surgery.

Consumer information

Proprietary rights

NeuroLogica® and BodyTom® Elite are registered trademark of NeuroLogica Corporation, a subsidiary of Samsung Electronics Co., Ltd., in the United States, other countries, or both. Catphan® is a registered trademark of Phantom Laboratory, Inc. ACR Appropriateness Criteria® is a registered trademark of the American College of Radiology. Image Gently® is a registered trademark of Society for Pediatric Radiology. Teflon® is a registered trademark of E.I. DuPont and Company. TB Quat™ is a trademark of ABC Compounding Co. Wex-cide™ is a trademark of Wexford Labs, Inc., product number Wexcide128.

Legal disclaimer

This user manual is intended as a guide for material supplied by NeuroLogica Corp. It provides the operator with necessary information to carry out specific procedures and maintain NeuroLogica produced equipment. Use this manual in conjunction with instruction and training supplied by qualified NeuroLogica personnel.

Any information or descriptions contained in this manual may not be reproduced and released to any of the general public or used in conjunction with other professional instruction without written consent of NeuroLogica Corp., USA – a subsidiary of Samsung. Direct any written inquiries to the appropriate address found in the section “Contact information” on page 25.

Unauthorized copying of this user manual may not only infringe copyright but also reduce the ability of NeuroLogica Corp. to provide accurate and up-to-date information to users: limited and restricted operators and administrators.

This user manual, though complete and accurate, may not provide answers to undocumented changes or unexpected results that could occur from system anomalies.

Contact information

Keep user information readily available to contact **Customer Service** about general assistance or reporting on serious incidents (should they occur).

In the case of a serious incident or adverse event, please notify NeuroLogica at the below contact information and establishments local competent authorities.

To provide any comments, suggestions, or corrections to this user manual, write to and include chapter title and page number:

NeuroLogica Corporation	
Customer Service	14 Electronics Avenue, Danvers, MA 01923 USA
USA and Canada	1-888-564-8561
International	1-978-564-8561
Email	support@neurologica.com

If you have any questions about faults or errors on the system, battery or charge issues, or mechanical issues with the scanner, contact a **Technical Representative**.

If you have questions about the clinical use of your system, building protocols, creating MPRs, imaging artifacts, creating a clinical workflow or process, logging in or access issues, and general usage of the system, contact **Customer Service**.

		
Winckels Medical Devices Expertise	Europe Bergerweg 18 6085 AT Horn The Netherlands	Tel: +31 (0)475 582285 Fax: +31 (0)475 582278
Australian Sponsor	Level 8/15 Talavera Road PO Box 646 North Ryde NSW 2113 Australia	Mobile: +61 (0)412 563 016 Tel: +61 (0)2 8114 1535 Fax: +61 (0)2 8114 1599 Customer Service: 1-800 060 168
Brazilian Authorized Distributor	2 nd floor, 515 Rua Peixoto Gomide Sao Paulo Brazil	Tel: 55-11-3371-1500

Damage in transportation

Closely examine all packages at the time of delivery. If you see damage, notate **“damage in shipment”** on all copies of the freight bill **before** you accept or sign for delivery (by the facility receiving agent).

Whether damage is noted immediately or concealed (noticed after delivery), damage **must** be reported to carrier **immediately** upon discovery, or within 14 days after receipt, and content and containers held for inspection by carrier.

Keep in mind – the transportation company **will not** pay a claim for damage if an inspection is not requested within the 14-day period.

User requirements

The equipment can *only* be operated by users who have received professional medical education and training, such as radiologic technologists, physicians, radiologists, and other medical specialists.

Users should be trained professionals who are certified to operate such systems **before** scanning or diagnosing patients. This training must include medical and x-ray education, as well as NeuroLogica applications training.

Everyone that uses this equipment must read, understand, and follow all instructions, precautions, and warnings.

Keep this user manual near the equipment. It is important to review the procedures and safety precautions periodically.

Essential Performance

The NL4000 BodyTom Elite has the following essential performance factors mitigated by design:

- Over Radiation protection
- Re-scan prevention
- Stray Radiation exposure prevention
- Diagnostic performance

About this user manual

The instructions in this user manual describe how to use the NeuroLogica NL4000 BodyTom Elite Computed Tomography (CT) system, manufactured by NeuroLogica Corp.

The manual applies to both old and new colors of BodyTom Elite.

This user manual **does not** provide medical explanations but does suggest potential applications for some of the software features. This user manual describes potential safety problems and how to avoid them.

Anyone who operates this system should have received training **before** attempting to scan or diagnose patients, to include medical and x-ray education, as well as NeuroLogica applications training.

This manual is made available in electronic format to the customer as part of each product delivery.

For electronic manuals, please go to: Forms.samsungneurologica.com.

Click on “Downloads” and choose “CT Manuals”. The site will ask for the serial number of your product and a password. The password can be provided from Field Service.

Translation of this manual is available for any country that does not allow for English labeling. Please reach out to NeuroLogica directly if translation is required.

Identified symbols and system classifications

The specifications and details of this user manual may change to improve the product or to enhance its performance.

Throughout this user manual, a yellow triangle with a black border and exclamation point is used to draw attention to those conditions or situations that fit one or more of the following criteria, which are definitions from ANSI Z535.5:



DANGER Indicates a hazardous situation, which if not avoided **will** result in death or serious injury.



WARNING Indicates a hazardous situation, which if not avoided **could** result in death or serious injury.



CAUTION Indicates a hazardous situation, which if not avoided **could** result in minor or moderate injury.

Conventions used in this user manual

Table 1: Conventions used in this user manual

Convention	Use
Commands to perform actions	To perform a string of commands, this user manual will present them as follows: Customize > System. This means click Customize and then click System .
Bold	When content refers to commands, windows, screens, dialog boxes, popups, tabs, buttons, options, keyboard keys, statuses, and modes, these items appear in bold for faster identification, especially in a procedure.
<i>Italic</i>	Identifies a word that is emphasized for your attention.
Numbered steps	Numbered paragraphs represent sequential steps that require you to take the action in the sequence provided – unless otherwise instructed. Procedures that are numerical mean that the sequence is important to follow. You may perform some procedures out of the recommended sequence; however, the results may vary.
Note	The appearance of a note is as such:

Convention	Use
	Note Indicates additional information to help you operate this product.
Hyperlink (an electronic cross-reference)	<p>A cross reference appears in the electronic (.pdf) user manual as a hyperlink. To retrieve an electronic copy of this user manual (in .pdf), click Help > User Manual from the workstation.</p> <p>A hyperlink is a quick way to go to another area of the user manual with a simple click. Hyperlinks appear like this: “Understanding the types of users” on page 101. In this case, hover the mouse pointer over the gray hyperlink text. Hold the Ctrl key on your keypad and click the mouse button. After you click the hyperlink, the hyperlink takes you to the referenced area in the user manual.</p>
Click vs right-click	In this user manual, click means to press the left mouse button. This user manual never says ‘left click’ as it is assumed that is the traditional way to click; however, it does point out when to right click the mouse button.

Understanding the use of “you” in this user manual

Unless specifically noted, the implied “you”, in this user manual, is the user/operator. It is assumed users/operators are certified and medically trained personnel, qualified to use these systems.

The following identifies those actions each user is permitted to perform:

Administrator	Full access to the system and its configurations. Can create protocols, User ID’s, and passwords, as well as access all functions of the system.
Limited operator	Modified access to the system. Users with Limited access can modify protocols during scanning but cannot create and save protocols; has no access to system configurations.
Restricted operator	Users with Restricted access can scan with the system but are unable to make any changes to protocol parameters while scanning, they also have no access to system configurations.

Active and inactive objects

When a menu command, option, button, tab, field, is gray, the item is not active or enabled. When an item is gray, it can mean additional or required tasks must be completed first or you do not have permission to access that option. An active menu command, option, button, tab, and field means you can use the item to perform an action. Active items are green and/or highlighted.

Applicable versions of BodyTom Elite

This user manual is applicable to all versions of BodyTom Elite, despite the color variations.



Figure 1: Current cover (left) and former covers (middle & right)

Chapter 1 Compliance and Safety Requirements

It is important that you are aware of and familiar with compliance and safety requirements to ensure you, the patient, and the systems are safe at **all** times.

IEC classification and symbols



In accordance with International Safety Standard IEC 60601-1, the BodyTom Elite CT scanner is classified as Type B equipment; Class 1 equipment, internally powered equipment, and continuous connection to the supply mains in standby state and for specified loading.

Type B equipment provides an adequate degree of protection against shock, regarding:

- Allowable leakage current.
- Reliability of the protective earth connection.
- (In accordance with the International Safety Standard IEC 60601-1), the manufacturer is not responsible for any consequences caused by the unauthorized modification of this equipment.
- Earth leakage current.



WARNING Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.



















WARNING Equipment is not suitable for use with oxygen or oxygen enriched atmospheres.















BodyTom Elite Computer Tomography systems comply with Class I- Type B equipment as defined in IEC 60601-1 standard.

Mode of operation is a continuous connection to the supply mains in standby state and for specified loading conditions.

The BodyTom Elite CT scanner is patient-environment equipment.

Table 2: Applicable symbols

Symbol	Description
	Alternating current
	Protective earth (ground)
	Functional Earth
	Caution: consult accompanying documents
	Caution: risk of electrical shock
	Electrostatic sensitive devices
	Type B equipment
	X-ray warning
	X-ray source assembly emitting
	Non-ionizing radiation
	Warning: laser in use
	Warning: Laser Radiation Do Not Stare into Beam Class 2 Laser Product
	Laser Output and Standards Information Label
	Warning: FDA Laser Information
	Warning: high temperature
	Emergency switch
	Crush warning

Symbol	Description
	Foot/toe crush warning when lowering machine
	System up
	System down
	Temperature limits
	Keep away from rain for packaging
	Humidity limit for packaging
	Warning: battery charging
	Fuse usage
	Refer to instruction in user manual/booklet
	Medical Device Symbol
	Legal Manufacturer Symbol
	CE Mark or Conformité Européenne ; number below CE represent Notified Body number
	Intertek ETL (Edison Testing Laboratories) Mark
	European Authorized Representative Symbol



WARNING This x-ray unit may be dangerous to patient and operator unless safe exposure factors and operating instructions are observed.

Environmental specifications



CAUTION The specified environment must be constantly maintained: 24-hours a day, seven days a week.

Table 3: Operating environment



Operating	
Ambient temperature	15° C to 35° C (59° F to 95° F)
Relative humidity	20% to 80% (non-condensing)
Altitude	0-3048m (0-10,000 ft.)
Storage	
Temperature	-25° C to 70° C (-13° F to 158° F)
Relative humidity	20% to 85% (non-condensing)
Transport	
Temperature	-20° C to 60° C (-4° F to 140° F)
Relative humidity	20% to 85% (non-condensing)
Powering system	
Time period prior to powering the system	24 hours ¹
Floor	
Flatness	<+/- 0.120in. (3mm) per 10ft. (3.048m)
Recommended minimum scan area	10ft. x 15ft. (3.048m x 4.572m)

Note For good image quality, the recommended practice is to keep the system free from vibration and to maintain the flatness specification noted.

Table 4: System operating parameters and specifications

Phase	120V~ Single	240V~ Single
Voltage range	100-120V~±10%	208-240V~ ±10%
	Nominal voltage for U.S. 120, 240	
Circuit protection	20 amps	30 amps
Facility outlet	NEMA 5-20R	NEMA 6-30R
	Outline of outlet:	

¹ If the system is in a facility outside the noted operating temperature, it is recommended to allow the noted time for the system to acclimate to the environment.

Phase	120V~ Single	240V~ Single
		
Frequency	50 or 60Hz	50 or 60Hz
Battery capacity	Fully charged, 12 hours typical	
Typical usage	110-120V~ 60Hz	230-240V~ 50Hz
Wiring	125V, 2 pole, 3 wire grounding	250V, 2 pole, 3 wire grounding
Main power supply's apparent resistance	0.105Ω	
Heat dissipation (when system is not idle)	55kW	

Considerations when preparing gantry for use



CAUTION Check for obstructions before moving and system setup.



CAUTION Monitor scanner motion to prevent collision with surrounding environment and foreign objects.



CAUTION Press the red **EMERGENCY STOP** button immediately in case of abnormal or unexpected motion.



WARNING Verify scanner is on its centipedes (fully down position) prior to positioning patient at scanner entrance.



WARNING Make sure all extremities are not under scanner while lowering or raising it.



WARNING In the case of a single pendant failure, the additional pendant is available for use to prevent loss of system function.



WARNING Always keep patient in view. Ensure that the patient can be seen when the operator is near the scanner control panel and **EMERGENCY STOP** button. Never leave patient unattended when the patient is in the gantry.

- NeuroLogica advises complying with local regulations and/or site recommendations as specified by the facility physicist or certified representative for the following:

- Use mobile x-ray protective-shielding devices. Technologists should be at correct location and consider wearing personal radiation protective equipment.
- A radiation safety plan in the working area boundaries, to optionally include as needed, mobile x-ray protective shielding devices. Otherwise assign a larger, working area to avoid radiation to the public. Effective dose for people outside the working area should be less than 0.25mSv annually (equals to 5 uSv weekly). The air kerma rate 0.3 meters away from the working area will be smaller than 2.5 uGy/h. Have monitoring and personal dose management for occupational exposure and related public health care personnel.
- There should be a working plan before scanning. The plan should include CT condition, time, location, working area, scanning plan, and site-clearing method; clearly state the responsibilities of working, protection, and management personnel. Keep a good record of the whole process.
- Restrict the working control and monitor area. Place obvious warning signs at the control-area boundaries to prevent unauthorized personnel from entering. Installation of a working status indication light is recommended.
- In accordance with the safety plan, self-monitor during the scanning process. A certified radiation representative should monitor the working area and take measurements immediately if abnormal circumstances are detected. Additionally, this should be reported to the local environmental administrative and health departments. There should be a public notice at the working area, to include the nature of work, time, location, control area, name of the working department, person in charge of the project, contact telephone number, radiation report telephone number.

Site specification

Note For good image quality, the recommended practice is to keep the system free from vibration and to maintain the flatness specification noted.

Table 5: System operating parameters

Operating voltage	100-240 VAC~
Operating frequency	50Hz-60Hz
Apparent resistance of supply mains at 120VAC	0.3 ohms
Operating current at 120VAC	13 amps
Heat dissipation	1672 watts



CAUTION For domestic purposes, scanner can be powered using either 120V~ or 240V~. If the scanner is using 120V~, the facility's circuit must be capable of providing 20 amps (single phase). If the scanner is using 240V~, the facility's circuit must be capable of providing 30 amps (single phase). If other devices are connected to the same circuit, the facility's circuit breaker may trip and, therefore, prevent the scanner from being ready when needed.

Table 6: Battery operating parameters

Operating voltage	480 to 585 VDC
Output current _(peak)	100 amps

Site specification for enclosed CT room

Table 7: Site specification

Issue	Comment
Receiving area	Secured
Packing material and waste	Near availability of a trash receptacle for dunnage
Room dimensions for use	12ft. x 15ft. room with a finished level floor; recommended the room be well lit
Power availability	120VAC/20amp wall outlets (2x)
Floor flatness	$\leq \pm 0.120$ in. (3mm) per 10ft.
Floor strength	Site must be able to support product weight

Note Not all beds are compatible with this system. Please contact **Customer Service** for assistance.

- NeuroLogica advises complying with local regulations and/or site recommendations as specified by the facility physicist or certified representative for the following:
- There should be enough space inside the CT room. The area should not be less than 30m². Any side should not be less than 4m. Leaving any unnecessary items inside the room is prohibited.
- Carpet and soft material cannot be used on the floor. There should not be obstacles on the floor. Ensure flat surface area no less than 12ft. x 15ft. level degree $\leq \pm 3$ mm per 10ft.
- Appropriate protection measures should be taken to protect staff and to ensure annual-dose-rate is less than 0.25mSv (equals to 5 uSv weekly).

Hazardous substances

Table 8: Hazardous substances

Substance/material	≅ Weight/system
Lead	7.69 kg (17.0lbs.)
Cadmium	0.036kg (0.079lbs.)
Mercury	0kg (0lbs.)
Hexavalent chromium	0kg (0lbs.)
PolyBrominated Biphenyls (PBB)	<0.46kg (1lb.)
PolyBrominated Diphenyl Ethers (PBDE)	<0.46kg (1lb.)

Part numbers and product-marking plates

Table 9: Core-system-component part numbers and product-marking plate locations

Component	Part number	Product-marking plate locations
BodyTom Elite gantry	0-NL4000-001 0-NL4000-002 10-00345-0001	Near the main input plug or on the side of the system. See Figure 2 below.
BodyTom Elite workstation	40-00157-000	On the back of the workstation.
QA phantom	10-00268-001	On the back of the phantom.

Note The applicable components making up the BodyTom Elite CT scanner is identified with the nameplate statement “This product complies with radiation performance standards, 21 CFR sub-chapter J.”




Figure 2: Product-marking plate on scanner

Table 10: Core-system component dimensions

Component / mode		Size L x W x H
BodyTom Elite	Scan	256.54cm x 103.1cm x 199.6cm
		101in. x 40.58in. x 78.57in.
	Transport	256.54cm x 103.1cm x 205.7cm
		101in. x 40.58in. x 81.00in.
	Bore	85cm
		33in.

Weight: 3510 lbs 1592 kg




Figure 3: Scanner dimensions including drivebar

Table 11: Workstation dimensions

Component	Size (inches) L x W x H	Size (centimeters) L x W x H	Weight (lbs)	Weight (kg)
BodyTom Elite system Workstation Cart	26.3 x 24.4 x 79.8	66.8 x 62.0 x 202.7	207	94

Class 1 Type B medical devices

This equipment generates, uses, and can radiate radio-frequency energy. The equipment may cause radio-frequency interference to other medical and non-medical devices and to radio communications. To provide reasonable protection against such interference, this product complies with emission limits for Class 1 medical devices as stated in EN 60601-1-2.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment is found to cause interference (which can be determined by switching the equipment on and off), the user should attempt to correct the problem using one or more of the following measures:

- Re-orient or relocate the affected device(s).
- Increase the separating space between the equipment and the affected device.
- Power the equipment from a source different from that of the affected device.
- Consult the point of purchase or the service representative for further suggestions.

NeuroLogica Corp. is not responsible for any interference caused either by the use of interconnect cables other than those recommended or by unauthorized changes or modifications to this equipment. Unauthorized changes or modifications could void the user's authority to operate the equipment.

To comply with the regulations applicable to an electromagnetic interface for a **Group 1 Class A** medical device, note the following:

- All interconnect cables to peripheral devices **must be** shielded and properly grounded.
- Use of cables not properly shielded and grounded may result in the equipment causing radio-frequency interference in violation of the European Union's Medical Device Directive and FCC regulations.



CAUTION Ensure there is no potentially detrimental interaction of system's irradiation with a patient's active and implantable medical devices and/or body-worn and active medical devices.



CAUTION Do not use devices that intrinsically transmit radio waves, such as a cellular phone, radio transceiver, mobile radio transmitter, radio-controlled toy, and so on. Use of these devices near this equipment could cause this equipment to malfunction. Keep power of these devices turned off when near this equipment. Medical staff in charge of this equipment is required to instruct technologists, patients, and other people who may be around this equipment to fully comply with the above regulation.

- Medical staff in charge of this equipment are required to instruct technologists, patients, and other people who may be around this equipment to fully comply with the above regulations.

Focal spot

Nominal size is: ~1.2 x 1.4mm

Size limit is: ~1.2 to 1.7mm width and ~1.4 to 1.9mm length.

Testing standard for reference is as follows:

- IEC 60336:2005

Focal spot centering is within 1mm of center of bolt pattern. Maximum motion due to gravity in X, Y, and Z axis is 0.1mm.

Maximum motion from anode rotation is 0.1mm.

Maximum motion from anode heating in X axis is 0.1mm. Maximum motion from anode heating in Z axis is 0.3mm.

Anode input power

The maximum, anode, cooling rate is 8,750W (12,250 HU/sec).

The maximum, anode, heat dissipation is 3,400W (4,760 HU/sec).

The nominal, anode, input power is 42kW.

Continuous anode input power when applied at the nominal, x-ray, tube-voltage is 150kV, 23mA.

Filtration

Table 12: Filtration

Tube Voltage (kV)	100	120	140
Half-value layer (aluminum equivalent)	6mm	7mm	8mm
Filters consist of 0.0014in. [0.036mm] of copper and 0.086in. [2.18mm] of aluminum, along with a variable thickness bowtie filter made from Teflon®.			
X-ray tube's total filtration of irremovable layers is 5.7mm of equivalent aluminum.			

Source to Detector distance (SID)

The SID value is 1041.9mm.

Compliance statement

Note All editions and years of revisions for standards noted in this chapter are static as of **Revision 15** of this user manual (and Revision 02 of the **BodyTom Workstation User Manual**), which has been incorporated into this user manual to include the scanner and the workstation.

The BodyTom Elite system complies with the regulatory requirements of the following:

- AAMI ES60601-1 Issue: 2005 Version - Medical Electrical Equipment, Part 1: General Requirements for Basic Safety and Essential Performance; Amendment No. 2: 2010/05/30.
- CAN/CSA-C22.2 No. 601.1-M90(R2005) Issued: 2003/11/01 Medical Electrical Equipment – Part 1: General Requirements for Safety; General Instruction No. 1: 1990, Supplement 1: 1994, Amendment 2: 1998, General Instruction No. 2: 2003.
- CENELEC EN 60601-1 2nd Edition, Medical Electrical Equipment - Part 1: General Requirements for Safety, includes Amendment A1:1993 and A2:1995.
- CENELEC EN 60601-1 3rd Edition, Medical Electrical Equipment - Part 1: General Requirements for Safety.
- CSA C22.2#60601-1 Issued: 2008/02/01 Ed 3 Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance.
- IEC 60601-1 Issued: 2005/01/01 Ed 3 Medical Electrical Equipment Part 1: General Requirements for Basic Safety and Essential Performance.

- IEC 60601-1-1 (2000) 2nd Edition: Medical Electrical Equipment, Part 1: General Requirements for Safety - Collateral Standard: Safety Requirements for Medical Electrical Systems.
- IEC 60601-1-2:2014 Ed4.0 Medical Electrical Equipment - Part 1-2: General Requirements for Basic Safety and Essential Performance - Collateral Standard: Electromagnetic Disturbances
- IEC 60601-1-3 (2008), Medical Electrical Equipment Section 1-3: General Requirements for Safety. Collateral Standard: General Requirements for Radiation Protection in diagnostic X-ray Equipment.
- IEC 60601-1-4 (2005), 3rd Edition Consolidated Edition, Medical Electrical Equipment Part 1-4: General Requirements for Safety – Collateral Standard: Programmable Electrical Medical Systems.
- IEC 60601-1-6 Issued: 2008/12/08 Ed 2 Medical Electrical Equipment - Part 1-6: General Requirements for Safety. Collateral standard: Usability.
- IEC 60601-2-28 (1993) Medical Electrical Equipment, Part 2: Particular Requirements for the Safety of X-ray Source Assemblies and X-ray Tube Assemblies for Medical Diagnosis.
- IEC 60601-2-32 (1994) Part 2-32, Particular Requirements for Safety sections 2.32 Specification for Associated Equipment of X-ray Equipment.
- IEC 60601-2-44 (2009) Medical Electrical Equipment -Part 2-44: Particular Requirements for the Safety of X-ray Equipment for Computed Tomography.
- IEC 60825-1:2007 Safety of Laser Products - Part 1: Equipment Classification, and Requirements 2nd Ed.
- International Electrotechnical Commission (IEC) International Standards Organization, when applicable.
- Intertek Testing Service (ITS), an independent testing laboratory.
- Medical Device Good Manufacturing Practice Manual issued by the FDA (Food and Drug Administration). Department of Health, USA.
- NeuroLogica Corporation is ISO 13485:2016 and MDSAP certified.

EMI/EMC terms

Electromagnetic Compatibility

Electromagnetic Compatibility (EMC) is the branch of electrical sciences that studies the unintentional generation, propagation, and reception of electromagnetic energy with reference to the unwanted effects **(Electromagnetic Interference (EMI))** that such energy may induce. The goal of EMC is the correct operation, in the same electromagnetic environment, of different equipment, which use electromagnetic phenomena and the avoidance of any interference effects.

To achieve this, EMC pursues two different kinds of issues. Emission issues are related to the unwanted generation of electromagnetic energy, to the countermeasures that should be taken to reduce such generation, and to avoid the escape of any remaining energies into the external environment. Susceptibility or immunity issues, in contrast, refer to the correct operation of electrical equipment in the presence of unplanned electromagnetic disturbances.

Interference, or noise, mitigation, and hence EMC is achieved by addressing both emission and susceptibility issues, that is., quieting the sources of interference, making the coupling path between source and victim less efficient, and making the potential victim systems less vulnerable.

Electromagnetic Interference (EMI)

Electromagnetic Interference (EMI), also called **Radio Frequency Interference (RFI)** is an unwanted disturbance that affects an electrical circuit due to electromagnetic radiation emitted from an external source. The disturbance may interrupt, obstruct, or otherwise degrade or limit the effective performance of the circuit. The source may be any object, artificial or natural, that carries rapidly changing electrical currents, such as an electrical circuit, the sun, or the northern lights.

Susceptibility

Susceptibility is the capability of an electronic system to respond to unwanted electrical energy.

EMI/EMC compliance

This equipment complies with IEC 60601-1-2 EMC standard for medical devices.

The BodyTom Elite system is suitable to be used in an electromagnetic environment, as per the limits and recommendations described in the tables hereafter:

- Emission Compliance level and limits (see Table 14 on page 78).
- Immunity Compliance level and recommendations to maintain equipment clinical utility.

Note This system complies with the above-mentioned EMC standard when used with supplied cables. If different cable lengths are required, contact a qualified service representative for advice.

Table 13: Acronyms and abbreviations

Acronym and abbreviation	Definition
AEC	Automatic Exposure Control
CBF	Cerebral Blood Flow
CBV	Cerebral Blood Volume
CT	Computed Tomography
CTA	CT Angiography
CTP	CT Perfusion
CTDI_{vol}	Volume Computed Tomography Dose Index
CTDI_w	Weighted average Computed Tomography Dose Index
DICOM	Digital Imaging Communication in Medicine
DLP	Dose Length Product (DLP)
DHCP	Dynamic Host Control Protocol
EMC	Electromagnetic Compatibility
EMI	Electromagnetic Interference
FOV	Field Of View
HIS	Hospital Information System
HU	Hounsfield Unit
IBC	Iterative Bone Correction
MAR	Metal Artifact Reduction
MIP	Maximum Intensity Projection
MPPS	Modality Performed Procedure Step
MPR	Multi-Planar Reformation, sometimes referred to as Multi-Planar Reconstruction
MTT	Mean Transit Time
PACS	Picture, Archiving, and Communication System
QA	Quality Assurance
RIS	Radiology Information System
RSO	Radiation Safety Officer
RFI	Radio Frequency Interference
SCP	Service Class Provider
SCU	Service Class User



WARNING Medical, electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to EMC information provided in accompanying documents.



CAUTION Portable and mobile RF communications equipment can affect medical electrical equipment.



CAUTION Do not use or stack the equipment or system with other equipment and if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

Note The EMC tables and other guidelines included in this user manual provide information to the user essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use; to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment.


Note The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Table 14: Emission declaration for BodyTom Elite systems

BodyTom Elite system is intended for use in electromagnetic environment specified below. The user of the BodyTom Elite system should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment guide
RF emissions CISPR 11	Group 1	BodyTom Elite systems use RF energy only for internal function. Therefore, RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	BodyTom Elite systems are predominantly intended for use in non-domestic environments, and not directly connected to the Public Mains Network. BodyTom Elite systems are predominantly intended for use (for example, in hospitals) with an appropriate power supply (see operation manual) and recommended shielding for portable use.
Harmonic emissions, IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions, IEC 61000-3-3	Complies	

Table 15: EMC Immunity declaration for BodyTom Elite systems

BodyTom Elite systems are intended for use in the electromagnetic environment specified below. The customer or user of an BodyTom Elite system should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2kV, ±4kV, ±8kV, 15kV air	±8kV contact ±2kV, ±4kV, ±8kV, 15kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV line-line ±2kV line-ground	±1kV line-line ± 2kV line-ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 40% UT; 5 cycles 70% UT; 25 cycles 0% UT; 250 cycles	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 40% UT; 5 cycles 70% UT; 25 cycles 0% UT; 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of a BodyTom Elite system requires continued operation during power interruptions, it is recommended that the BodyTom Elite system be powered from its internal batteries.
Immunity test	IEC 60601-1-2 Test Level	Compliance level	Electromagnetic environment guidance.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m, 50Hz or 60Hz	30 A/m, 50Hz or 60Hz	Power-frequency magnetic-fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Clause 8.10	Table 14	Per Table 14	IEC 60601-1-2:2014

BodyTom Elite systems are intended for use in the electromagnetic environment specified below. The customer or user of an BodyTom Elite system should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC 61000-4-6	3 Vrms 150kHz to 80MHz	V1 = 3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of an BodyTom Elite system, including cables, than recommended separation distance calculated from the equation appropriate for transmitter frequency. Recommended separation distance: See Table 16.
Radiated RF IEC 61000-4-3 (alternative method: IEC 61000-4-21)	3 Vrms 80MHz to 2,5GHz	E1 = 3 V/m	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is maximum power rating in watts and D is recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than compliance levels (V1 and E1). Interference may occur in vicinity of equipment marked with the following symbol:</p> 

Note: The wireless receiver operates within the following bands.

2.412 to 2.462 GHz (11 channels)

5.180 to 5.240 GHz (4 channels)

5.260 to 5.320 GHz (4 channels)

5.500 to 5.700 GHz (8 channels, excluding 5.600 to 5.640 GHz)

5.745 to 5.825 GHz (5 channels)

The preferred frequency band is 5.189 to 5.240 GHz at 40MHz bandwidth.

The wireless transmitter operates within the following frequency bands and power.

802.11b:

Typ. 26 ± 1.5 dBm @ 1 Mbps, Typ. 26 ± 1.5 dBm @ 2 Mbps

Typ. 26 ± 1.5 dBm @ 5.5 Mbps, Typ. 25 ± 1.5 dBm @ 11 Mbps

802.11g:

Typ. 23 ± 1.5 dBm @ 6 to 24 Mbps, Typ. 22 ± 1.5 dBm @ 36 Mbps

Typ. 20 ± 1.5 dBm @ 48 Mbps, Typ. 19 ± 1.5 dBm @ 54 Mbps

802.11n (2.4 GHz):

Typ. 23 ± 1.5 dBm @ MCS0/8 20 MHz,

Typ. 18 ± 1.5 dBm @ MCS7/15 20 MHz

Typ. 23 ± 1.5 dBm @ MCS0/8 40 MHz,

Typ. 17 ± 1.5 dBm @ MCS7/15 40 MHz

802.11a:

Typ. 23 ± 1.5 dBm @ 6 to 24 Mbps, Typ. 21 ± 1.5 dBm @ 36 Mbps

Typ. 20 ± 1.5 dBm @ 48 Mbps, Typ. 18 ± 1.5 dBm @ 54 Mbps

802.11n (5 GHz):

Typ. 23 ± 1.5 dBm @ MCS0/8 20 MHz,

Typ. 18 ± 1.5 dBm @ MCS7/15 20 MHz

Typ. 23 ± 1.5 dBm @ MCS0/8 40 MHz,

Typ. 18 ± 1.5 dBm @ MCS7/15 40 MHz

The device includes 4 dBi gain antennas

Countermeasures against EMC related issues

Generally, it is very difficult to grapple with issues related to EMC. It may take a variable amount of time and cost to identify issues causing interference.

General countermeasures of electromagnetic interference with other equipment:

- Electromagnetic interference may be alleviated by positioning other equipment far from system.
- Electromagnetic interference may be mitigated by changing relative location (installation angle) between system and other equipment.
- Electromagnetic interference may be eased by changing wiring locations of power/signal cables of other equipment.
- Electromagnetic influence may be reduced by altering the power-supply path of other equipment.
- Electromagnetic environment specified (see Table 15 on page 48).

Table 16: Recommended separation distances

Recommended separation distances between portable and mobile RF communications equipment and the BodyTom Elite system			
BodyTom Elite system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the BodyTom Elite system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BodyTom Elite system as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output Power (P) if transmitter Watts (W)	150kHz to 80MHz Separation distance meters ²	80MHz to 800MHz Separation distance meters ²	800MHz to 2,5GHz Separation distance meters ²
0.01	.12	.12	.23
0.1	.38	.38	.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

² Separation distance according to frequency of transmitter (m)

Recommended separation distances between portable and mobile RF communications equipment and the BodyTom Elite system			
BodyTom Elite system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the BodyTom Elite system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BodyTom Elite system as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output Power (P) if transmitter Watts (W)	150kHz to 80MHz Separation distance meters ²	80MHz to 800MHz Separation distance meters ²	800MHz to 2,5GHz Separation distance meters ²
For transmitters rated at a maximum output power not listed above, the separation distance is estimated using the equation in the corresponding column, where P is the maximum output power rating of the transmitters in Watts (W) according to the transmitter manufacturer.			
Note At 80MHz and 800MHz, separation distance for higher frequency range applies. Note These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

Table 17: Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test Frequency (MHz)	Band ^{a)}	Service ^{a)}	Modulation ^{b)}	Max Power (W)	Distance (m)	Immunity Test Level (V/m)
385	380-390	Tetra 400	Pulse Modulation ^{b)} 18Hz	1.8	0.3	27
450	430-470	GMRS 460,	FM ^{c)}	2	0.3	9

Test Frequency (MHz)	Band ^{a)}	Service ^{a)}	Modulation ^{b)}	Max Power (W)	Distance (m)	Immunity Test Level (V/m)
		FRS 460	±5kHz deviation 1 kHz sine			
710	704-787	LTE Band 13,17	Pulse Modulation ^{b)} 217Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900 TETRA 800, iDEN 820, CMDA 850, LTE Band 5	Pulse Modulation ^{b)} 18Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CMDA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation ^{b)} 217Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN,	Pulse Modulation ^{b)} 217Hz	2	0.3	28

Test Frequency (MHz)	Band ^{a)}	Service ^{a)}	Modulation ^{b)}	Max Power (W)	Distance (m)	Immunity Test Level (V/m)
		802.11 b/g/n, RFID 2450, LTE Band 7				
5240	5100-5800	WLAN 802.11 a/n	Pulse Modulation ^{b)} 217Hz	0.2	0.3	9
5500						
5785						

Use recommendations

This product complies with IEC 60601-1-2 standard for medical devices and with radio frequency emission requirements per CISPR11 Group 1 Class A standard limits. The BodyTom Elite system is predominantly intended for use in hospitals.

Do not use devices that intentionally transmit RF signals (cellular phones, transceivers, or radio-controlled products) in the vicinity of this equipment as it may cause performance outside the published specifications. Keep the power to these types of devices turned off when near this equipment.

Adhering to the distance separation (recommended in Table 16 on page 83) between 150kHz and 2.5GHz, will reduce disturbances recorded at the image level, but may not eliminate all disturbances; however, when installed and operated as specified herein, the system will maintain its essential performance by continuing to safely acquire controlled, radiological, x-ray exposures in a mobile radiography environment. For example, a 1W mobile phone (800MHz to 2.5GHz carrier frequency) is put 2.3 meters apart from the BodyTom Elite system (to avoid image interference risks).

The use of accessories, transducers, and cables, other than those specified, may result in degraded, electromagnetic compatibility of the BodyTom Elite system.

The medical staff in charge of this equipment is required to instruct technologists, patients, and other people who may be around this equipment to comply fully with the above equipment requirements.

Installation recommendations

This system complies with above-mentioned EMC standard when used with supplied cables. To minimize interference risks, the following requirements apply.

Cable shielding and grounding

All interconnect cables to peripheral devices must be shielded and properly grounded. Use of cables not properly shielded and grounded may result in the equipment causing radio-frequency interference.

Adjacent components and equipment

BodyTom Elite system should not be used adjacent to or stacked with other equipment; if adjacent or stacked use is necessary, the BodyTom Elite system should be tested and verified to make sure normal operation in the configuration in which it is used. Consult NeuroLogica and Facility **Technical Support** staff regarding device/system configurations.

Static magnetic field limits

To avoid interference on the BodyTom Elite system, static-field limits from the surrounding environment are specified. The static field is specified less than <1 Gauss around the unit.

Electrostatic discharge environment and recommendations

- To reduce electrostatic-discharge interference, install a charge-dissipative floor-material to avoid electrostatic charge-buildup.
- The relative humidity must be at least 30 percent.
- The dissipative material must be connected to the system ground-reference.

Facility IT-NETWORK

The BodyTom Elite system utilizes the IT-NETWORK for the customer (as applicable) to communicate with the modality worklist and for supplemental, image-storage space. As part of the installation, the applicable IT-NETWORK is reviewed to create the appropriate setup to the system. Setup is done to ensure no potential concerns arise with the system.



CAUTION It is possible that the IT-NETWORK connection from the system could result in previously unidentified issue(s) to the respective population. Should this occur, contact Customer Service right away to identify, analyze, evaluate, and resolve the issue(s).



CAUTION It is possible that any changes to the IT-NETWORK made by the facility could introduce a new issue where Customer Service needs to be contacted to address the concern, right away.

Changes to the IT-NETWORK include – but are not limited to the following:

- Changes in network configuration
- Connection to additional items
- Disconnection to items
- Updating equipment
- Upgrading equipment.

Hazard information

Review this material before using the system and observe basic, common-sense safety rules when operating this scanner.

General safety considerations and statements

Review the following before using the system (Scanner and Workstation (as applicable)) to observe basic, common-sense safety rules when operating the scanner:

- Become familiar with the functional hardware to help recognize serious problems.
- Do not use scanner or workstation if it appears damaged or fails.
- Wait for qualified personnel to correct any problem.

Note The scanner is provided with a video-camera monitoring-system to help navigate the unit while being transported within a facility.



WARNING Modification of this equipment is **not** allowed.



WARNING Equipment maintenance of non-medical electrical equipment should not be performed in the patient environment.



CAUTION All non-medical electrical equipment will comply with relevant IEC and ISO safety standards.



CAUTION Federal law restricts the use of this device without a prescription by a physician.



CAUTION Always store and/or use unit in a well-ventilated area. Keep air pollution to a minimum. Always keep floor clean.



CAUTION Do not touch parts of non-medical electrical equipment in patient environment and patient simultaneously.



CAUTION For disposal of any material emanating from the system; follow local regulations.



CAUTION This system was designed for use by individuals trained in CT system operation. The user should be familiar with this user manual before scanning patients.



CAUTION It is the user's responsibility to make sure that after installation or subsequent modification, the system complies with the requirements of collateral standard IEC 60601-1.



WARNING Installation of this product is performed in accordance with Installation Manual (1-NL4000-059). All installation processes and qualified personnel are outlined in that document.



WARNING Proper disposal of batteries is required to ensure compliance with environmental safety guidelines. Contact authorized NeuroLogica representative for instructions.



WARNING Observe safety-exposure factors and operating procedures to protect patient from physical harm during contact with this x-ray scanner.



WARNING Observe safety requirements to prevent excessive dose exposure to patient and/or operator.



CAUTION Improper system (including workstation) usage could endanger patients and/or users and void the warranty if not operated correctly.



CAUTION Should the workstation encounter a computer related virus, be sure to contact Technical Support for assistance with removing said virus from the equipment.



CAUTION Radiation dose exposure to patients should not exceed maximum of 1Gy CTDI.



CAUTION For proper disposal of material at equipment's end-of-useful life; contact NeuroLogica for instructions.



WARNING Equipment in which protection against electric shock relies on basic insulation *only*, should not be used in this system.



WARNING If the system fails to move due to loss of power, the patient can be easily removed from the scanner by moving the patient bed.

Four, foot crush hazard labels are affixed to the scanner in four places, above the four soft bumpers. The following shows a safety label:



Figure 4: Identifying the scanner's safety label(s) – foot-crush-hazard label(s)

Laser safety

There are four lasers used with the BodyTom Elite system as indicated in Figure 5 on page 59: 1 laser (Coronal) at position **1**, 1 laser (Sagittal) at position **2**, and 1 external and 1 internal laser (Transverse or Axial) at position **3**. Laser **2** is mounted internally (fixed) to the disk assembly, which spins within the system's bore. Therefore, the lasers output light is always aimed at and rotating within the bore itself.



WARNING Viewing the laser output with certain optical instruments (for example, eye loupes, magnifiers, and microscopes) within 100mm may pose an eye hazard.



WARNING Complies with 21 CFR 1040.10 and 1040.11 except for conformance with IEC 60825-1 Ed. 3., as described in Laser Notice No. 56, dated May 8, 2019.

This statement will be a label placed onto the system by the Driver Side.

Laser parameters

- Lasers **1** and **3** (see Figure 5):
 - Wavelength = 650nm
 - Output Power = 1mW
- Laser **2** (see Figure 5):
 - Wavelength = 650nm
 - Output Power = 4mW




Figure 5: Laser aperture's direction



CAUTION Instruct the patient to close his/her eyes before turning **ON** the alignment light.



CAUTION Closely monitor infants and infirm patients to prevent them from accidentally staring into the beam.



CAUTION Class 2 laser radiation when open. **Do not stare into the beam** or view directly with optical instruments.



CAUTION Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.



CAUTION The warning label (below, for “laser in use”) is located on the front of the scanner cover and inside the scanner to identify the presence of a laser.



Scanner mobility safety



CAUTION Due to the mobility of the system, an external interlock is not available; however, a prescribed scan can, at any time, be terminated from the scanner or the workstation. When the user activates the scan, a 10-second, countdown-clock, scan delay (adjustable to 99 seconds) triggers. This countdown allows the user time to perform needed tasks before the scan begins.



CAUTION To prevent healthcare provider injuries, a single healthcare professional should not move the scanner and workstation. Although one person can drive the BodyTom Elite when moving the scanner about the facility, NeuroLogica recommends two people move the scanner (lengthwise, only) to ensure no collisions occur when maneuvering through tight hallways and around corners. Be especially cautious when moving the system about an inclined floor.



WARNING To prevent involuntary movement, do not position scanner on an incline while in **Transport** mode.



WARNING Do **not** move the system right or left if transport on an incline becomes necessary. Always keep the system in a straight motion.



WARNING Contact **Technical Support** for assistance when movement is required on an incline.

Note Be sure there are no obstacles in front of the scanner while you move the scanner.

If the system needs to be moved over a threshold it is critical that the scanner be oriented so that it is driven in the forward or reverse direction. The scanner does not have the capability of moving laterally over thresholds.



CAUTION Check to ensure proper clearance is provided to allow removal of patient from scanner in case of a power failure. This is accomplished by moving patient's support (after unlocking wheel-locks) away from scanner.



CAUTION To prevent patient entrapment or entanglement with accompanying equipment, slowly move scanner away from patient by using control panel switches or pendant controls while observing patient.



CAUTION The scanner is equipped with a video camera to help the operator prevent collisions when transporting system to different locations that could otherwise result in personal injury or facility damage.



CAUTION Do not station or operate the system on an uneven floor. The flatness requirement is 0.12in. over 10ft. or 3mm over 305cm.



CAUTION Prior to transporting the scanner, verify that power cord is unplugged from wall to avoid damage to cord and outlet, and avoid tripping. Verify that the ethernet cable is unplugged from the workstation to avoid damage to the cable and connector.

Floor level

For proper operation, the system must be operated on an even, level, hard surface.

Carpeting

Do not use the system on a carpeted floor. Poor image quality could result due to unevenness of the floor.

Electrical safety



WARNING The system's external AC power cord should be checked prior to use to verify there are no exposed wires or damaged insulation/prongs. Damaged prongs could result in sparking and fire. In case of such damage, contact Customer Service, immediately.



WARNING To separate the device from AC power simply disconnect the power cord from the wall and turn off the main breakers, located on the same side of the scanner as the power cord.



WARNING Access to the main breaker is critical for safety. Do not position the scanner so that the access to the breaker is diminished.



CAUTION Check to ensure the AC outlet is working properly before plugging in the system's AC power cord. NeuroLogica recommends using a dedicated outlet for powering the BodyTom Elite system, *only*.



WARNING To prevent electrical shock, do not connect items that are not specified as part of the system, including the workstation.



WARNING To prevent electrical shock, do not remove the covers from the equipment. The covers protect the user and the patient from moving parts or electrical shock. Hazardous voltages are present within this equipment. The covers provide protection from radiation exposure given off from the x-ray tube. The covers also protect the equipment.



WARNING An electrical shock hazard: no user should apply serviceable parts; refer to qualified service personnel for any service.



WARNING Always electrically isolate this equipment from the main electrical supply before cleaning and disinfecting it to prevent short-circuiting or possible electrical shock.



WARNING Never position the mobile system and/or workstation in a manner that prohibits access to unplugging it or prohibits pressing the **EMERGENCY STOP** button.



WARNING To minimize shock hazard, the system chassis must be connected to an electrical ground. The system is grounded through the ground conductor of the supplied, three-conductor power cord. The power cord must be plugged into a three-conductor electrical outlet receptacle. Do not alter the ground connection.



WARNING Avoid all contact with any electrical conductor as follows:

- Allow only qualified personnel who know the proper procedures and use the proper tools to install, adjust, repair, or modify the equipment.
- Only use this equipment in rooms or areas that comply with all applicable laws (or regulations having the force of law) concerning electrical safety for this type of equipment.
- Always electrically isolate this equipment from the main electrical supply before cleaning and disinfecting it.
- The detachable cord is the disconnecting device, which is used to remove mains power from the wall socket.
- The system is internally powered.



WARNING For Class 1 equipment (for example, the workstation, AC power cord, and so on) using an alternate internal source: a warning to use the alternate source if the integrity of the protective earth conductor is in doubt.



WARNING Do not position the system so that it is difficult to access the AC power cord.



CAUTION Protect the system power cord against mechanical damage.

Where the integrity of an external, protective conductor, in the installation or its arrangement, is in doubt, equipment is operated from its internal, electrical power source.

Parts of non-medical electrical equipment in the patient environment that, after removal of covers, connectors, and so on, without the use of a tool, may be contacted by the operator during routine maintenance and calibration, will operate at a voltage not exceeding 25VAC or 60VDC or peak value supplied from a source that is separated from the supply mains in accordance with one of the methods described in IEC 60601-1.



CAUTION All systems within the patient environment will provide the same level of safety as medical equipment complying with IEC 60601-1.



CAUTION To help prevent tripping hazards, use care in the arranging of any cords (for example, AC cord, ethernet cable, and so on) when connecting to the system/workstation.



CAUTION To prevent damaging electrical outlet cords, check to ensure they have been removed and properly stored before transporting the scanner.



WARNING The BodyTom Elite CT scanner contains high-voltage circuits for generating x-rays. Only trained and qualified personnel should be permitted access to the internal parts of this equipment.



CAUTION Use the **Universal Serial Bus (USB)** terminal located near the **EMERGENCY-STOP (E-STOP)** button for archiving to USB, *only*. **Do not** use the USB terminal located near the **E-STOP** for connecting any other device to equipment.



CAUTION For proper disposal of material at the end of the useful life of the equipment, contact NeuroLogica for instructions.

Mechanical safety



WARNING In case of unwanted movement or motion, press the **EMERGENCY STOP (E-STOP)** button.



WARNING Physically assist all patients on and off the bed and into position on the scan board. Adjust the bed to the specified height for patient loading and unloading (see “Positioning the patient” on page 97).



WARNING Adjust scanning platform to specified height for patient loading and unloading; see “Positioning the patient” on page 97.



WARNING When positioning the scanning platform, be careful when moving the scan table to avoid having it hit the scanner covers.



WARNING Position any lines (IVs and so on) attached to the patient so the lines cannot catch on the scanner during scanner travel.



CAUTION Prevent pinching or crushing of the patient’s extremities. Keep patient’s hands on the side of his/her body. Always watch the patient and equipment carefully during scanner movement.



CAUTION To prevent pinching or crushing of the operator’s feet/toes, be sure extremities are not positioned under the scanner when it is being lowered from **Transport** mode to **Scan** mode.



WARNING Equipment maintenance of non-medical electrical equipment should not be performed in the patient environment.



WARNING Maintenance checks and all service must be performed by NeuroLogica trained service personnel.



CAUTION Ask patient to scoot up into universal scan board or manually aid them into position.



CAUTION When the scan board is in place, be especially careful when moving the bed to avoid driving it into the gantry covers.










CAUTION Periodically check all accessories for damage and remove them from service if damaged or cracked.

Radiation safety

Two **Dangerous to patient and operator** labels are affixed to both sides of the scanner, just above the operator controls. See Figure 6 on page 65.



Figure 6: Dangerous-to-patient/operator safety-warning label location (left) and label (close-up, right)

-  **WARNING** Improperly used x-ray equipment may result in unwanted radiation exposure. Read and understand the instructions in this user manual before attempting to operate this equipment.
-  **CAUTION** Use technique factors prescribed by the radiologist or diagnostician. Use a dose that produces the best diagnostic results with the least x-ray exposure.
-  **CAUTION** All persons authorized to use the equipment must understand the dangers posed by excessive x-ray exposure. NeuroLogica recommends use of protective materials and devices.
-  **WARNING** Everyone having anything to do with x-ray must take adequate steps to insure protection against injury.
-  **CAUTION** The use of this device requires its users to receive proper training in accordance with local and national laws.
-  **CAUTION** **Never** perform calibration with patients in the scanner or while personnel are present in the vicinity of the scanner to prevent exposure to unwanted radiation.
-  **CAUTION** Amber indicator lights (on the top of the scanner) illuminate during x-ray exposure.



CAUTION Ensure that there is no potential for detrimental interaction of the system's irradiation with a patient's active implantable medical devices and/or body-worn, active, medical devices.

X-rays can only be produced during the following conditions:

- The scanner is in the **Scan** mode position.
- The workstation is connected.
- The **START** button is activated when the patient is registered, the protocol is selected, the **Begin** button is clicked, and the protocol is prepared. The **START** button on the scanner's control panel illuminates when the scanner is ready to begin. See Figure 8 on page 67 to identify the **START** button.

Fire and explosion safety



DANGER This equipment is not suitable for use in presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.



DANGER This equipment is not suitable for use in the presence of oxygen or oxygen-enriched atmosphere.

- Fire regulation for the type of medical area being used should be fully applied, observed, and enforced. Fire extinguishers should be provided for both electrical and non-electrical fires.
- All operators of the BodyTom Elite scanner should be fully aware of and trained in the use of fire extinguishers and the firefighting equipment, and in local fire procedures.



WARNING Only use extinguishers on electrical or chemical fires that are specifically labeled for those purposes. Using water or other liquids on an electrical fire can lead to fatal or other serious injury.

If it is safe to do so, attempt to isolate the equipment from electrical and other supplies before attempting to fight a fire. This will reduce the risk of electrical shocks.

EMERGENCY STOP button



CAUTION Check the **EMERGENCY-STOP (E-STOP)** button at least one time a month to ensure proper function.



CAUTION Every user should take a few minutes to locate the **E-STOP** before scanning the first patient.



CAUTION In case of emergencies, stop scanner movement immediately by pressing the **E-STOP** red push-button located on the scanner, below the control panel.




Figure 7: BodyTom Elite E-STOP locations (right and left)




Figure 8: Close-up of the scanner control panel and the E-STOP button

Note When the **E-STOP** button is activated, the moving gantry may overrun by less than 10mm.

Battery safety and information

The **System battery capacity** icon shows an indication of the scanner's battery capacity, which is identical to the indicator on the scanner. The user should always check the indicator on the scanner to verify the batteries' status; there are 145 lithium-ion batteries in the scanner; 144 are used for scanning and the remaining battery is used for moving the scanner while in transport mode.



CAUTION The system unit contains batteries and will always be charging when plugged into AC mains.



CAUTION In case of battery leakage, do not handle the batteries themselves nor continue to operate the system. Contact NeuroLogica for service. See “Contact information” on page 25.

Battery replacement and disposal



WARNING Battery replacement is to be performed by authorized and trained NeuroLogica service personnel, to ensure proper disposal of hazardous material.



WARNING Dispose batteries in accordance with federal, state, and local regulations.



WARNING Do not incinerate batteries.



WARNING Contact an authorized NeuroLogica representative for appropriate, product-disposal instructions.

Scanner battery capacity

The range is 0 to 100. Battery voltage and current are used to calculate the system’s battery capacity based on charging state.



Figure 9: Scanner battery capacity icon

Run time operation

During normal, run-time operation, the battery capacity is being calculated one time, per second. The **battery capacity** indicator is updated on the scanner’s **display screen** as well.

Note The **battery capacity** indicator on the scanner’s **display screen** is displayed in 5% increments when above 10%. When 10% and below, the displayed capacity is in 1% increments. That is, above 10% capacity, the displayed value is always rounded to the closest multiple of 5 (for example, 93 gets rounded to 95, 42 gets rounded to 40, 47 gets rounded to 45, and so on).

State changes

After each periodic update, capacity is checked to make sure it does not fall below certain thresholds, as follows:

- Low voltage alarm state

When **battery capacity** goes below 25%, a periodic alarm will sound. It will remain in this state until the battery capacity has gone back up to 27% or higher. The **Start** and **Cancel** buttons when scanning will not illuminate if the battery is too low.

- Low voltage lock-out state

When the battery capacity goes below 1%, the scanner screen-display buttons are disabled and starting a scan is prohibited; for example, the ability to move the scanner and certain protocol buttons are disabled. It will remain in this state until the battery capacity has gone back up to 2% or higher. The low voltage alarm will continue to be active in this state.

Predictive scanning

Before each scan, battery usage for that scan is predicted based on the selected load factors (for example, kV, mA, scan time) and is compared against the available battery capacity. If there is not enough battery capacity to scan, a popup appears on the workstation screen. The user can cancel the scan at that time or continue the scan, with the understanding that the scan may abort due to a low-power fault.

Under voltage protection

When the **system battery voltage** drops below the low-voltage cutout-level while unplugged, a system power-down sequence is initiated.

Workstation

The **Workstation battery capacity** icon shows an indication of the workstation's battery capacity. On the workstation's main screen, place the mouse pointer over the battery icon to see the capacity of the battery, ranging from 0 to 100%. The user should always check the screen to verify the status of the batteries.



Figure 10: Workstation battery capacity icon

Note: If the scanner's **display screen** is black, the system is not charging and/or the batteries are permanently damaged. A service call is required.



CAUTION In newer workstations, the workstation reports battery capacity on the workstation remote power display (under the monitor). The battery system is designed to be replaced by authorized and trained NeuroLogica service personnel, *only*.



CAUTION The workstation will not report the proper battery capacity and status if a network connection is not made.



CAUTION The system can only be charged from a correctly rated wall outlet. A rating information plate is located on the product system label (lower backside panel or lower left side panel, see Figure 2).



CAUTION The system (including the workstation) should be plugged in at all times, when not in **Transport** mode, being transported or in scanning use, to help maintain battery life and proper system operation. Failure to do so could result in permanent battery damage, which will require a service technician to repair.



CAUTION The system may not complete a scan when below 25% battery capacity while unplugged.



CAUTION If the system is unplugged and battery capacity reports to be 0%, permanent battery damage can occur.



CAUTION The power cord selection must not be less than 110v/12A (USA) and 220v/7.5A (EU and Asia), made of 2.08mm (diameter) copper wire in accordance with local power supply cable standards.

Note Medical grade power cords should be used at all times.

Scanner X-ray tube capacity

The percentage of the tube capacity required for a scan = $((\text{kV} \times \text{mA} \times \text{scan time(s)}) / 180000) \times 100\%$. Approximately 0.11% capacity is regained each second during cooling.



Figure 11: Scanner X-ray tube capacity icon

Maintenance and service



WARNING Equipment maintenance of non-medical electrical equipment should not be performed in the patient environment.



WARNING Maintenance checks and all service must be performed by NeuroLogica trained service personnel. Service personnel use Service manual (1-NL4000-062) to effectively perform needed service and preventive maintenance and inspection of the system. See “Contact information” on page 25 for NeuroLogica’s contact information.



WARNING The only calibration performed by the user on this system is called **daily calibration** and is described in detail later in this user manual. All other calibration needs that arise must be performed by trained service personnel at NeuroLogica Corp. See “Contact information” on page 25 for NeuroLogica’s contact information.



CAUTION Service personnel must complete training at NeuroLogica Corp. for the system and its accessories prior to conducting any service activities.

Users are not to perform service or maintenance on the system at any time. This includes battery maintenance.

Note NeuroLogica recommends that a quarterly preventive maintenance be conducted by NeuroLogica’s service personnel/trained facility bioengineer.

Instructions for replacing serviceable parts are identified in the Service Manual (1-NL4000-062).

Cybersecurity



WARNING Upon detection of a cybersecurity threat to the system or workstation, do the following:

- Immediately contact Technical Support.
- Discontinue use of system (enabling the **EMERGENCY-STOP** if needed).
- Remove any ethernet and/or wireless connection that has been made with the facilities’ IT-network.

Continued use of the system can occur after Technical Support has assessed the situation and provided the go-ahead to do so.

Note NeuroLogica Corp. recommends the customer facility utilize an IT-network that provides sufficient means of cybersecurity control to help maintain the requirements of HIPAA.

Contraindication(s)

There are no contraindications associated with CT scanning.

Personnel privileges and terminology

Qualified operator

The operator as determined by the healthcare facility and assigned by a user with administrative privileges – who by their education, certification, experience, and training, are sufficiently qualified to competently perform clinical scans on the CT system. See “Understanding the types of users” on page 101 for a description of the types of users.

Operator of record

The operator of record is an operator or health care professional currently logged onto the CT system with a unique User ID and password identifier.

Scanning privileges

Scanning privileges are granted to a qualified operator, assigned by a user with administrative privileges, to conduct clinical scans on the CT system. This privilege level allows use of all clinical protocols to scan the patient.

Protocol privileges

Protocol privileges are granted to a qualified healthcare professional, as determined by the healthcare facility, and assigned to users with administrative privileges, who by their education, certification, experience, and training, is sufficiently qualified to competently save or modify clinical protocols on the CT system. A healthcare professional with protocol privileges does not necessarily need scanning privileges on the CT system.

Administrative privileges

Administrative privileges are granted to qualified healthcare professionals as determined by the healthcare facility who by their education, certification, experience, and training, is sufficiently qualified to competently assign, maintain, and oversee the assignments of personnel to scanning privileges and/or protocol privileges on the CT system they administer. In addition, qualified healthcare professionals are authorized and qualified to pull system logs associated with this standard for **Quality Assurance** review. Healthcare professionals with administrative privileges do not necessarily need scanning privileges or protocol privileges on the CT system.

Clinical operation

CT system operation that involves scanning patients and/or creating or editing protocols intended for use on patients.

Clinical scanning

CT system operation that involves scanning of patients.

Clinical protocol

A protocol on the system intended for use on patients.

Kernel

The kernel or filter is defined as the reconstruction algorithm or mathematical equation used for convolution of the attenuation profiles and reconstruction of the CT images. The choice of the kernel determines the noise level and the contrast resolution of the reconstructed images.

Chapter 2 System Overview

The BodyTom Elite CT system lets you scan patients in a room or ward, an Emergency Room (ER), Operating Room (OR), Radiology, and procedure rooms.

This chapter provides a brief introduction to the BodyTom Elite system. You will learn how to use the BodyTom Elite system (to see an illustration, see Figure 12 on page 75) – in subsequent chapters.

Note Whether you turn on the scanner or the workstation first (see Chapter 4) does not matter; however, it is advised to power up the BodyTom Elite system hardware first, to allow time for the scanner to warm up.

BodyTom Elite system

The BodyTom Elite is a 32-slice, mobile, battery-operated CT scanner and software system with **Axial**, **Helical**, and **Dynamic**, capabilities.

The BodyTom Elite core system consists of the scanner, the workstation, and the phantom holder. Consider the following:

- The scanner and workstation communicate using a wireless connection. They communicate using an ethernet connection, if necessary.
- The BodyTom Elite workstation is a computer with custom software that allows the user to employ pre-defined, system protocols or devise unique protocols for performing patient studies. It also allows the user to update patient information and store images. The viewing portion of the BodyTom Elite workstation allows the user to view images in more detail and includes tools to help facilitate diagnosis by a physician.
- The maximum scout length is 2000mm.
- The scanner can scan a slice-thickness of 1.25mm, 2.5mm, 5.0mm, and 10mm in **Axial** mode.
 - In **Axial** mode, the BodyTom Elite scans 10mm of anatomy with each rotation.
 - The maximum scan-range in **Axial** mode is 900mm.
- The scanner can scan a slice-thickness of 1.25 x 0.625, 1.25 x 1.25, 2.5 x 1.25, 2.5 x 2.5, 5.0 x 2.5, and 5.0 x 5.0 in **Helical** mode.
 - In **Helical** mode, the BodyTom Elite scans 32mm of anatomy at a pitch of 0.8.

- The maximum scan-range in **Helical** mode is 2000mm.
- In **Dynamic** mode, the BodyTom Elite scans 40mm of anatomy.
- The scanner is compatible with surgical navigation, **HIS**, **RIS**, and **PACS**.



- 1 QA stand and phantom
- 2 BodyTom Elite or gantry
- 3 BodyTom Elite workstation
- 4 Bed with patient near bore

Figure 12: BodyTom Elite system configuration

For information on cleaning and storing the scanner, see “Cleaning the scanner and workstation” on page 357 and/or “Storing the system” on page 359.

Overview of the scanner control panels

Control panels appear on the right and left side of the scanner. The **scanner control panel** allows the operator to start and stop a scan. All motion and x-ray generation can also be quickly stopped using the **EMERGENCY STOP** button. After the patient is registered and you select a protocol, you must start the scan from the **scanner control panel**, which is located under the scanner’s display screen. The **operator control panel** lets you power on the scanner, lift, and lower the scanner, and lock scanner functions.



Note To determine where personnel should stand during a scan, consult with the hospital physicist. See “Scatter radiation” on page 244.




Controls on the left end of the scanner



Figure 13: Left end of the scanner

Table 18: Left end of the scanner

Scanner control	Description
	Video camera shows what is in front of the scanner during transport.
	The scanner's display screen shows the status of the scanner and workstation (see Table 24 on page 116 for a list and description of each status symbol).

Scanner control	Description
	START and CANCEL buttons (to the left) and the EMERGENCY-STOP (E-STOP) button (to the right) on the scanner control panel.
	Ethernet port to connect to the workstation, and four additional data-access ports on the data interface panel.
	120VAC/20A and 240VAC/30A plugs.
	The pendant ; see “Overview of the pendant” on page 81 for more information.





Controls on the right end of the scanner



Figure 14: Right end of the scanner

Table 19: Right-end of the scanner

Scanner control	Description
	<p>The scanner's display screen shows what the video camera captures and where the scanner is moving during transport.</p> <p>It also shows the status of the scanner and workstation (see Table 24 on page 116 for a list and description of each status symbol).</p>






Scanner control	Description
	START and CANCEL buttons (to the left) and the EMERGENCY-STOP (E-STOP) button (to the right) on the scanner control panel.
	Drive bar and enable bar (arrow).
	Key lock , the Power-ON and Power-OFF buttons for scanner, and the Rocker-Switch-Lift Up and Down buttons on the operator control panel.
	The pendant ; see “Overview of the pendant” on page 81 for more information.


Identifying operator control panel buttons



Figure 15: Operator control panel buttons and indicators

Table 20: Operator control panel buttons and indicators




Operator control panel buttons and indicators	Name	Description
	LOCK	Use the key to lock or unlock the operator control panel buttons. If the key is in the locked position all scanner buttons are disabled. If the key is in the unlocked position all scanner buttons are enabled.
	OFF	Press to power down the entire scanner. During shutdown, the light blinks until the shutdown task is complete.
	TRANSPORT	Press to activate Transport mode. Use the Rocker-Switch-Lift Up and Down buttons to put the scanner in Transport mode and the drive bar to transport the scanner.
	ON	Press to power up the scanner. During power-up, the light blinks until power-on task is complete.
	BATTERY BARS	Indicates the scanner's battery charge level. Each bar represents 10% of charge. If plugged into an AC outlet, the last bar blinks indicating the system batteries are charging.










Operator control panel buttons and indicators	Name	Description
	Rocker-Switch-Lift Up and Down buttons	<p>Press and hold the UP or Down Rocker-Switch-Lift button to raise or lower the scanner.</p> <p>Lowering the scanner to floor level, makes the scanner ready to scan. Raising the scanner, makes the scanner ready for transport.</p> <p>When the button is inactive, it is dim; when the button is active it is illuminated.</p>

Overview of the pendant

The pendant lets you move the scanner, turn on and off the lasers, zero reference the scanner, and program scan and rest positions for the scanner. See Table 21 for a list of what each button activates.

Table 21: Pendant buttons

Pendant	Button	Description	Action
		POWER	Illuminates when power is supplied to pendant.
		LASER	Turns on all three positional lasers. While the lasers are on, the scanner spins for the internal laser to be seen within the scanner opening.
		GO TO SCAN PLANE	Moves the scanner forward approximately 30cm. This is the distance between the internal and external lasers.

Pendant	Button	Description	Action
		ZERO REFERENCE	Sets the scanner to zero before starting a scout or a scan.
		MOVE BACKWARD (slow)	Pressing and holding moves the scanner backward 10mm per second.
		MOVE FORWARD (slow)	Pressing and holding moves the scanner forward 10mm per second.
		MOVE BACKWARD (fast)	Pressing and holding moves the scanner backward 60mm per second.
		MOVE FORWARD (fast)	Pressing and holding moves the scanner forward 60mm per second.
		SET MEMORY	Allows the user to program Scan and Rest positions for the scanner.
		SCAN POSITION	Moves the scanner to the Scan Position saved using the Set Memory feature.
		REST POSITION	Moves the scanner to the Rest Position saved using the Set Memory feature.

**DANGER**

Store the pendant in its holder when not in use to prevent accidental and/or unintentional contact by patient and/or users.



Figure 16: BodyTom Elite remote-control pendant

The scanner's position appears on the **positional display** on the front of the scanner. You can use the pendant to zero reference the scanner. The display shows a positive or negative positional number.



Figure 17: Scanner's positional display

Note The two light panels on either side of the scanner's position indicate x-ray is active, when lit. The scanner produces an audible alert during scanning.

Overview of the workstation




Figure 18: Workstation with leaded-glass shield (optional installation)

The **workstation** is an accompanying part of the scanner; it is the computer and control unit that operates most functions of the system. All basic information related to the workstation (for example, operating distance, warnings and cautions, connectivity, functionality, etc.) appear in Chapter 4 Basic Workstation Operations. The workstation includes the computer, monitor, ethernet connections, and the remote controls. The workstation also includes the **Uninterruptable Power Supply (UPS)**. The workstation can be installed with an optional leaded-glass shield (shown in Figure 18) for additional protection.

The workstation enables you to easily move it wherever you need it to go. The workstation is designed to let you navigate in and out of elevators, over doorway thresholds, or on any type of floor including carpet, with ease.

The workstation can be set up either wirelessly or hardwired to the scanner. The administrator makes sure wireless is enabled before you create a wireless connection between the workstation and the scanner (with **System Configuration > Scanner Setup**). See “Remote Support Setup” on page 134 to learn how the administrator enables the wireless connection.

Note Wireless connections can be slower than a hardwired connection. If an unexpected delay or disconnection, due to environmental/bandwidth interference, occurs when using wireless, it is recommended to hardwire an ethernet cable to the scanner for continued communication.

After respective peripherals are properly plugged in, make sure all applicable power switches are in the on position before using the system.

For any devices connected to an AC outlet, make sure the outlet is providing required power.



CAUTION When not in use, the scanner and workstation should always be plugged into power outlets to ensure maximum efficiency.

The system should be stored in an area with limited access to prevent inadvertent damage.

See “Powering on and off the BodyTom Elite system” on page 90.

Note Be sure to keep the workstation plugged in when it is not in use to charge the battery. Charge time while off is ~ (approximately) 2 hours; while in use, charge time is ~ 8 hours.



CAUTION Verify that the ethernet cable is unplugged from the workstation to avoid damage to the cable and outlet during transport.

Note Before using the workstation, be sure to read and understand how to clean and maintain it. See “Cleaning the scanner and workstation” on page 357.

Workstation types

There are two types of workstations available: Granite and Dell. For more information on powering on and off the workstation, see “Using the workstation” on page 101.

Note The Granite workstation does not allow the user to wait for the computer to shut down before communicating to the workstation to shut down; this is because the workstation sends a toggling signal to the computer. If the computer is off, it will turn back on. The workstation must already be in the process of shutting down by the time the computer is fully turned off.

The product safety coverage of the specified workstation (Safety Certified to IEC 60950 standards) was evaluated and deemed acceptable for use with the BodyTom Elite to appropriate IEC 60601-1 and/or IEC 60601-1-1 harmonized national standards.

The workstation will be kept outside the patient environment as defined by IEC 60601-1-1. BodyTom Elite is suitable for use inside patient environment.



CAUTION *Do not* connect or use equipment beyond what is specified by NeuroLogica Corp.; this practice may lead to a reduced level of system safety.

The recommended distances, provided in Figure 19, relate *only* to distances specified by IEC 60601-1 and *do not* relate to specific distances required for ionizing radiation and/or stray radiation protection for operators and bystanders.




Figure 19: Workstation safe distance location (two views)

Workstation remote power controls

The following are found on the workstation:

- Remote power controls
- Microphone and controls
- Speaker and controls



Figure 20: The remote power display
See “Using the workstation” on page 101.

Note The microphone and speaker are plugged in to the USB port on the back of the monitor. If the port is changed, you will need to reboot the workstation.

Workstation UPS

The workstation uses an **Uninterruptable Power Supply (UPS)** to supply power to the workstation for approximately (~) 6-8 hours when the workstation is unplugged. The **UPS** is located inside the workstation. This feature allows the new workstation to run when it is not plugged in to a wall outlet.

Workstation considerations before use

Before using the workstation, consider the following:

- If a problem is detected with the workstation, make sure repairs or adjustments are made to it **before** using it.
- Make sure the workstation operates easily and freely, and all parts work smoothly.
- Check for excess noise, vibration, or a change in the ease-of-use.
 - Noise, vibration, or change in ease-of-use can be signs of a problem and a need for servicing.
- Be sure to read its warnings carefully and completely **before** using the workstation. Do not attempt to service the workstation. Only skilled service personnel are permitted to service the workstation. See “Hazard information” on page 56.



CAUTION Users are not to perform service or maintenance on the system at any time. This includes battery maintenance.



CAUTION Failure to heed these warnings may cause injury to the user, to others, or damage to the equipment.

Note NeuroLogica recommends that a quarterly preventive maintenance be conducted by NeuroLogica's service personnel/trained facility bioengineer.

Instructions for replacing serviceable parts are identified in the Service Manual (1-NL4000-062).

Keyboard and mouse

The workstation comes with a keyboard and a mouse (**1** for left-mouse button; **2** right-mouse button).




Figure 21: Workstation keyboard and mouse

Workstation power cord

The power cord lets you power up the workstation.

Consider the following:

- Plug style will vary depending on factory installed elements based on geographic location and voltage requirements.



CAUTION The power cord selection must not be less than 110v/12A (USA) and 220v/7.5A (EU and Asia), made of 2.08mm (diameter) copper wire in accordance with local power supply cable standards.

Note Medical grade power cords should be used at all times.

Parts that potentially come into contact with the patient

While you use the system, be mindful that the patient may come into contact with the following parts:

- The BodyTom Elite system, especially the painted, external surfaces of the system's covers.
- Universal transfer board, if purchased.

Chapter 3 Basic Scanner Operations

Basic scanner skills include powering on and off the scanner, learning how to use and navigate the operator and scanner control panels, how to use **E-STOP**, and how to use the Rocker-Switch-Lift **UP** and **DOWN** buttons to lift and lower the scanner.

Scanning basics you should know before scanning a patient include how your system should be setup, how to position the scanner and the patient before the scan, and how to start a scan from the operator control panel.

Note It is recommended that the scanner is on for at least 60-90 minutes prior to performing the daily air calibration or scanning patients.

It is recommended practice that the scanner is plugged in and turned on even when it is not in use.

Powering on and off the BodyTom Elite system

The BodyTom Elite is not intended to be turned on and off; however, if the system should lose power, **it is advised** to power on the scanner first, to allow time for the scanner to power up.

Consider the following:

- Make sure the scanner is properly plugged in whenever possible; be sure the outlet(s) provide the required power. Plugging the electrical cord into the wall charges the batteries; the batteries are the only power source that allows the scanner to perform scans.
- When plugging in the scanner, make sure the cable lays flat on the floor to ensure the safety of hospital personnel. In addition, make sure that the floor behind the scanner is free of any obstructions or debris that could interfere with the centipedes during scanning.




Figure 22: AC cord and storage on scanner (120V left plug in or 240V right plug out)



Figure 23: Scanner's power cord receptacle for 120VAC and 240VAC

Be sure to keep the workstation plugged in when it is not in use to charge the battery. Charge time while off is ~2 hours; while in use, charge time is ~8 hours.

Checking a connection between the workstation and the scanner

To check if a wireless connection exists between the scanner and the workstation, look for the **Wireless connection** icon on the scanner's display screen: .

Note You must be logged into the workstation before this icon is visible on the scanner's display screen.

If the workstation is connecting to the scanner by hardwire, check if the supplied ethernet cable is connected between the workstation and scanner.




Figure 24: Scanner hardwired to the workstation with an ethernet cable to data-access ports

The administrator makes sure wireless is enabled when creating a wireless connection between the workstation and the scanner (with **System Configuration > Scanner Set Up**, through the workstation). See “Remote Support Setup” on page 134.

Note Make sure that the wireless setup in System Configuration has been done before making a wireless connection between the workstation and the scanner. See “Scanner Setup” on page 133 and “Remote Support Setup” on page 134.

Wireless connections can add lag time when compared with a hardwired setup. If an unexpected delay or disconnection occurs when using wireless, it is recommended to hardwire an ethernet cable to the scanner for continued communication.

Moving and transporting the scanner

To move the scanner, the scanner must be in **Transport** mode. The Rocker-Switch-Lift **UP** and **DOWN** buttons are located on the side of the scanner, on the operator control panel. These **UP** and **DOWN** buttons prepare the scanner to move up for transporting or down for positioning the scanner before scanning a patient. To set the scanner in **Transport** mode, go to the operator control panel and follow the procedure below.

Note If the scanner has been calibrated for multiple rooms, you must select the room prior to scanning to ensure the correct floor-calibration file is loaded. See “Selecting a room for the BodyTom Elite” on page 188.



CAUTION

Before transporting the scanner, verify that the ethernet cable is unplugged from the workstation to avoid damage to cable and receptacle. Verify that the power cable is unplugged from the wall to avoid damage to the cord and outlet.

1. Press the **Transport** button.




Figure 25: Transport button on the operator control panel

2. Press and hold the Rocker-Switch-Lift **UP** button to raise the scanner off the centipedes for transport.



Figure 26: Rocker-Switch-Lift UP (top) and DOWN (bottom) button

Note Hold the top Rocker-Switch-Lift **UP** button until the scanner is completely raised and on its castor wheels. The top Rocker-Switch-Lift **UP** button illuminates when it is in the up position as shown in Figure 26.

3. Grip the **drive bar** with both hands.



Figure 27: Drive bar front



Figure 28: Drive bar

Note The **drive bar** is equipped with an **enable bar**. The **enable bar** is a darker soft rubber material located on front of the **drive bar**. While holding the **drive bar**, squeeze the **rubberized enable bar** to activate system motion. If you let go of the **enable bar**, the scanner stops.

4. Move the scanner.

The following points explain how to hold the **drive bar** to move it forward, backward, left, and right:

- Push forward with both hands and use equal pressure to move the scanner forward.
- Push with the right hand to turn the scanner left; push with left hand to turn the scanner right.
- Pull back on the **drive bar** to move the scanner in reverse.

Note A three-point driving technique is required to turn in smaller spaces.

**CAUTION**

NeuroLogica recommends two people move the scanner within the facility: one to steer and another, in front of the scanner, to insure there are no obstacles. Two people ensures there are no collisions while maneuvering through tight hallways and around corners.

**WARNING**

If a loss of control is encountered while moving the system, release the **enable bar** to stop **ALL** movement.

**WARNING**

When transporting, use the video camera and scanner's display screen as guides to avoid hitting objects.

**CAUTION**

To prevent potential for injury from overbalancing and/or tipping, **do not** attempt to turn this system on an incline during transportation.

**WARNING**

Do not move the system right or left if transporting on an incline becomes necessary. **Always** keep the system in a straight motion. Contact **Technical Support** for assistance when movement is required on an incline.

**CAUTION**

The BodyTom Elite scanner is larger than most medical equipment. Therefore, NeuroLogica recommends proper training and practice.

Drive direction of scanner




Figure 29: Scanner drive direction (right side view)

Safety bumper system

The scanner's transport system is equipped with an active, **safety-bumper** system. Each safety bumper is electronically controlled to terminate motion in the direction in which the system is moving when the bumper is activated. If a bumper is pressed due to a collision, the drive system is disabled in that direction. All other directions are still enabled to allow you to reposition the scanner away from any impact area. The activation force needed to trigger the bumper system is approximately 7lbs.



Figure 30: Bumper system

Positioning the scanner before a scan



WARNING **Never** raise or lower the scanner when the patient is positioned in the system's bore. **Always** move the patient and table away from the bore before raising or lowering the system itself.

1. Ensure the scanner is in **Scan** mode.

To move the scanner, see "Moving and transporting the scanner" on page 92. See also "Performing a scan" on page 264.

2. To lower the scanner, press and hold the Rocker-Switch-Lift **DOWN** button until the scanner is completely lowered on its centipedes. The Rocker-Switch-Lift **DOWN** button illuminates when the scanner is completely lowered on the centipedes.
3. Align the patient and table with the scanner, ensuring the patient is in the center of the **Field Of View (FOV)**, also ensure that the scanner will not collide with the patient, patient support and/or any life supporting devices prior to scanning.

Note Be sure the floor is clear of debris or anything that can cause interference with the scanner's centipedes.

4. Make sure the patient scan table is locked.



CAUTION Make sure to lock the patient bed or scan table to prevent it from moving during the scan.

5. Adjust the scan table height so that the patient is centered within the bore.




Figure 31: Patient centered in bore (height positioning)

Positioning the patient using the laser lights



WARNING Before scanning, position the patient in such a way that extremities, hair, life-support equipment, etc. have sufficient clearance to prevent contact with scanner and or when used with accessories and options, such as head frames, scan boards, etc.



WARNING Make sure the patient is supported properly when positioned (both height and alignment) to prevent injury during scanning.



WARNING Make sure the foot pedal brake on the scan table/bed is engaged to prevent it from moving during the scan.



WARNING **Never** raise or lower the scanner buttons on the operator control screen when a patient is positioned in the system's bore. **Always** slide the scan table away (by disengaging its brake) from the system before raising or lowering the system itself.



CAUTION The following-instructions for patient positioning should be performed in accordance with NeuroLogica Corp.'s clinical training.



CAUTION If the patient becomes nauseated or is unable to be still (motionless), stop the scanner immediately using the **E-STOP** button.

Note To determine where personnel should stand during a scan, consult with the hospital physicist. NeuroLogica recommends a distance of 8-10 feet.

The scanner provides a rotating laser-light to guide you to properly position the patient. The laser indicates the center of the scan plane. The laser light is always centered on the actual x-ray beam. For multi-slice protocols, this means that the laser light will indicate the middle position of all simultaneous acquisitions being acquired. There is no offset between the laser-light plane and the actual scan-plane. The accuracy of the position of the laser plane, with respect to the scan plane, should be +/- 2mm.

There are two sets of laser lights: one set of external and one set of internal. The external set contains the following: 1 **Coronal** positioning, 2 **Transverse** laser, 1 mid-**Sagittal** laser. The internal set contains lasers affixed to the x-ray tube; the lasers rotate around the patient to designate internal zero-point reference. Keep the following in mind:

- To adjust vertical or horizontal positions, use table and bed controls, *only*.
- To adjust Z axis, walk the scanner to position using the pendant.

See “Overview of the pendant” on page 81 to see how buttons act.

1. Position the patient on the bed.
2. Align the bed to the scanner and make sure there is sufficient clearance around the patient by positioning the patient in the center of the **FOV**.



Figure 32: Phantom positioned in center of FOV

3. Use the pendant positioning buttons to center the patient within bore.



Figure 33: Pendant use for positioning lasers upon patient

4. On the pendant, press the **LASER** button to turn on all positional lasers. The following three lasers are available to help with positioning.
 - 1 set of external, lateral-side-positioning lasers (**Coronal** positioning for table height).
 - 1 external, crosshair-positioning laser (top, to designate the **mid-Sagittal** plane).

- 1 external crosshair positioning laser (side-to-side to designate the **Transverse** plane).




Figure 34: Positioning lasers upon patient
See the laser precautions in “Laser safety” on page 58

5. Use the scanner’s positional display to identify position of the scanner.



Figure 35: Positional display

6. Consider the following:
 - If the patient is conscious, request the patient remain still with eyes closed throughout the entire scan.
 - If the patient is unconscious, secure the patient.
 - Follow the appropriate facility guidelines when scanning unconscious patients if the patient’s eyes remain open.

Operating the E-STOP button




Figure 36: BodyTom Elite E-STOP locations (right and left)



Figure 37: E-STOP button on the scanner control panel on both the left and right sides of the scanner

1. Press the **E-STOP** button to perform the following:
 - Stop system (if it loses control).
 - Stop all system motion and x-ray.
 - Remove power to the gantry drives and x-ray system.
 - If the scanner starts to move unexpectedly.
2. Make sure to resolve the situation.

Note When **E-STOP** is activated, the moving gantry may overrun by less than 10mm.

Restoring the system from E-STOP

If you have pressed the **E-STOP** button to stop the system, follow these steps to restore the system, using **E-STOP**.

1. Make sure any hazard is removed.
2. Twist the **E-STOP** button clockwise until the button pops out to restore the system after pressing the **E STOP** button.

Chapter 4 Basic Workstation Operations

Basic workstation information includes understanding the different kinds of user access, learning how to power on and off the workstation, and getting familiar with the workstation screen, its menus, and tabs.

Note Whether you turn on the scanner or the workstation first (see Chapter 4) does not matter; however, it is advised to power up the BodyTom Elite system hardware first, to allow time for the scanner to warm up.

Understanding the types of users

There are three **User Levels** available on the workstation: administrator, limited operator, and restricted operator. User ID's and passwords can be created for individual users, and specific **User Levels** can be assigned to each user. The following define the access levels for each User Level.

Administrator	Full access to the system and it's configurations. Can create protocols, User ID's, and passwords, as well as access all functions of the system.
Limited	Modified access to the system. Users with Limited access can modify protocols during scanning but cannot create and save protocols; has no access to system configurations.
Restricted	Users with Restricted access can scan with the system but are unable to make any changes to protocol parameters while scanning, they also have no access to system configurations.

Using the workstation

The workstation uses an **Uninterruptable Power Supply (UPS)** to supply power to the workstation for ~6-8 hours when the workstation is unplugged. The workstation includes a computer, the remote power display, and a microphone and speaker.

Identifying the workstation's remote power display






The **Power On** button on the remote power display allows you to power on the Granite workstation.



Figure 38: Workstation remote power display

The following describes the action of each workstation power-control button.

Table 22: Workstation power-control buttons





Workstation power-control button	Button or indicator name	Description
	UPS battery level	Shows the battery usage; each LED represents 20 percent of battery power.
	Power-On	When the system has power, the LED light illuminates.
	Alarm	When the system is at or less than 20 percent battery power, an alarm sounds, and the LED light illuminates to warn you.
	Power-On and Power-Off	Press and hold the Power-On / Power-Off button for 3-5 seconds to turn on and off the workstation.
	Mute	Press this button to silence the alarm.

Identifying the microphone, speaker, and controls



Figure 39: Microphone, speaker, and controls

Table 23: Speaker control buttons

Microphone button	Name	Description
	Microphone mute	Press the Mic Mute button to mute the microphone. NeuroLogica recommends using the Mute button located on the bottom right of the workstation screen.
	Mute	Click the Mute button, on the bottom right of the workstation screen to mute the microphone and speaker.
	Speaker	The Speaker button appears on the bottom right of the workstation screen when the speaker is activated to hear the patient.
	Volume	Press the “-” button on the left side of the speaker to decrease the volume, press the “+” button on the right side of the speaker to increase the volume. The speaker has illuminated volume lights to indicate volume level.

Powering the workstation




Figure 40: Remote power display on workstation

Note Depending on the workstation, the **Power-On** and/or **Power-Off** button may not be in this location.

To power up a Granite workstation

1. Press and hold the Granite workstation's **Power-On/Power-Off** button on the workstation's remote power display until the green light blinks.

To power down a Granite workstation

1. Press and hold the Granite workstations **Power-On/Power-Off** button on the workstations remote power display until the green light blinks.

Note The Granite workstation does not allow you to wait for the computer to shut down before communicating to the workstation to shut down; this is because the workstation sends a toggling signal to the computer. If the computer is off, it will turn back on. The workstation must already be in the process of shutting down by the time the computer is fully turned off.

To power up a Dell workstation

1. Press the computer **Power-On** button on the side of the workstation.



Figure 41: Dell workstation power on button

To power down a Dell workstation

1. Click **File**
2. Click **Shutdown Computer**.
3. On the **Shutdown Computer?** popup box, click **Yes**.




Figure 42: Shutdown computer popup

Logging in to the workstation

To gain access to the workstation application, you must provide the system with credentials. These credentials consist of a **User ID** and **password**. Make sure you have a valid User ID and password before you log into the system.

If you do not have a valid User ID and password, ask your supervisor or administrator for one.

1. Click the User ID dropdown and select your User ID.



Figure 43: User ID dropdown box




Figure 44: User ID dropdown list

2. Click in the **Password** field and type your password. Passwords are case sensitive.



Figure 45: Password text box

3. Click the **Login** button.

If the User ID and/or password are invalid, a prompt appears to re-enter the information.

Note You have a limited number of login attempts before the system locks the account. An administrator can unlock the account. See “System and User Configuration and Setup” on page 121 for more details about unlocking an account.

When the User ID and password are verified, you are logged into the system.

The main screen to the workstation software appears with the **Patient Registration** tab open.



Figure 46: Patient Registration tab

4. Verify that the correct User ID appears at the top center of the screen.



Figure 47: User ID, current date, and time

If the correct User ID **does not** appear, contact your supervisor or administrative user to verify the account.

Locking and unlocking the workstation

The **Lock** button indicates whether your system is in **Lock** or **Unlock** mode. This button is a toggle, if the workstation is locked, the **Unlock** button shows, if the workstation is unlocked, the **Lock** button shows.

To prevent unwanted personnel from accessing the system, you should lock the workstation if you intend to leave the area for any period. When you lock the workstation, it remains **on**, but no one can access it without supplying a User ID and password.

Using Lock if you need to step away from the workstation

1. Click the **Lock** button located at the top-and center of the workstation screen.




Figure 48: System Lock button

The **Lock/Unlock System** popup appears.



Figure 49: Lock/Unlock System popup to lock the workstation

2. Enter the User ID and password by selecting the option from the dropdown and entering information in the **Password** field.
3. Do one of the following:

- Click the Lock button to lock your system.
When you select this option the **Lock** button changes to the **Unlock** button.
- Click the **Cancel** button to return to your work.

Using Unlock to view your work

1. Click the **Unlock** button.




Figure 50: Unlock button

The **Lock/Unlock System** popup appears.



Figure 51: Lock/Unlock System popup to unlock the workstation

2. Enter the User ID and password by selecting the options from the dropdown and entering information in the **Password** field.
3. Do one of the following:
 - Click the Unlock button to unlock your system.
When you select this option the button changes to the **Lock** button.
 - Click the **Cancel** button to return to your work.

Navigating around the workstation's main screen

The workstation screen shows similar information, even if the screen elements change when you perform different actions. These constant screen elements are as follows:

Main menu options	Appears in the top left corner of the screen and shows the File , Tools , Customize , and Help commands.
User, date, and time	Appears in the center of the screen and displays a greeting for the user, the signed-in user's name, the current date, and time.
Scanner and workstation status	Appears in the top right corner of the screen and displays status information for both the scanner and workstation.

The main screen to the workstation software always opens with the **Patient Registration** tab.

Brief overview of the main menu

The main menu provides you access to the basic functions from the commands: **File**, **Tools**, **Customize**, and **Help**.

The main menu appears on every screen and is always located in the same location, regardless of what you are doing.



Figure 52: Main menu

Note To select commands, click the first command and the subsequent commands. For example, if you see **Customize > System** that means, click **Customize** from the main menu and then click **System**.

Brief overview of the File menu

When you log off, restart, or shutdown the application and/or the workstation, you must use your User ID and password to log back in.

Logging off the system

1. Click **File** from the main menu.



Figure 53: File menu

The **File** dropdown appears.



Figure 54: File > Log Off

2. Click **Log Off** from the dropdown to shut down the software without shutting down the workstation.

The **Login** popup appears.



Figure 55: Login popup

This is also the login and shutdown portal. You can login or shutdown the workstation from this popup.

Restarting the application and/or computer

You can restart the computer in two ways: restarting the application only or restarting both the application and the computer.

1. Click **File** from the main menu.



Figure 56: File dropdown menu

2. Perform one of the following:

- Click Restart Application from the dropdown to restart only the application software.
- Click **Restart Computer** to restart both the application and the computer.

The following **Restart Application** or **Restart Computer** popup appears.



Figure 57: Restart Application or Restart Computer popup

3. Perform one of the following:

- Click the Yes button to restart the workstation.
- Click the No button to return to the screen.

Shutting down the computer

1. Click **File** from the main menu.



Figure 58: File dropdown menu

2. Click **Shutdown Computer** from the dropdown.
The **Shutdown Computer** popup appears.

3. Perform one of the following:

- Click the Yes button to shut down the workstation.
- Click the No button to return to the screen.

Brief overview of the Tools menu

This menu provides you with tools to store and print, set up protocols, and test your system to ensure it is operating as specified.




Figure 59: Tools dropdown menu

1. Click **Tools** from the main menu.
2. Click one of the following from the dropdown:

Store/Print Queue	The Store/Print Queue displays the status of studies being archived. You will learn more about how to store to various media later in this user manual; see page 121.
Protocol Manager	Allows users with Administrative privileges to create, modify, delete, and/or upload protocols to the scanner. You will learn more about how to use Protocol Manager later in this user manual; see page 189.
Quality Assurance (QA)	The tool that verifies the system is at its optimum performance. You will learn more about QA later in this user manual see page 222.

Brief overview of the Customize menu

This menu provides you with tools to set up the system as well as define user profiles.




Figure 60: Customize dropdown menu

1. Click **Customize**, then one of the following sub commands from the drop-down list:

System	Allows users with Administrative privileges to customize site-related settings; see "Chapter 5 System and User Configuration and Setup" on page 121.
User	Allows you to customize layouts in the system as well as set the password. See "Chapter 5 System and User Configuration and Setup" on page 121.
Select Room	Allows you to identify and select the room the scanner will be used in. See "Selecting a room for the BodyTom Elite" on page 188.

Getting Help from the Help menu

NeuroLogica Help includes an online user manual and information about the system. It also provides remote support from NeuroLogica **Technical Support** for file transfer, remote upgrades, or system review and support. When you enter a six-digit number, **Technical Support** will take control of the system to retrieve files or review the issue in question.




Figure 61: Help dropdown menu

Getting an online user manual

To open a .pdf version of this user manual:

1. Click **Help** from the main menu.
2. Click **User Manual** from the dropdown list.
A PDF version of this manual will be opened.

Getting remote support

1. Click **Help** from the main menu.
2. Contact NeuroLogica **Technical Support**.
See “Contact information” on page 25.
3. Click **Remote Support** from the dropdown list.
The **Support Connection** window appears.




Figure 62: Support Connection browser window

When connected to **Remote Support**, a NeuroLogica **Technical Support** representative will supply a six-digit code to start a remote support session, which allows the support representative to review your system and troubleshoot the issue.

Getting information about the product and NeuroLogica

To get additional information about the product and NeuroLogica:

1. Click **Help** from the main menu.
2. Click **About** from the dropdown list.
The **About Us** popup appears.




Figure 63: About Us popup

The following information is found:

Version(s)	Identifies the current software versions for the system.
Licensed To	Identifies who the product is licensed to.
Station AE Title	Identifies the title for your workstation (for PACS purposes).
Licensed Packages	Identifies any licensed packages available on the system.






Getting to know the status bar



The status bar appears in the top-right portion of the screen. The status bar provides a quick view of the systems current state. Details for the icons on the status bar are in the tables below.





Figure 64: Scanner and workstation status bar

Table 24: Status bar icons



Status bar icon	Status bar icon name	Status description
	Radiation status	Identifies x-ray as on or off. The icon changes from a gray/black icon when x-ray is off, to a rotating yellow/black icon when x-ray is on.
	System state	Identifies the system's current state. The orb changes color depending on the state the system is in. See Table 25 on page 118 for a list of the different orb colors and system states they identify.
	Scanner position	Identifies the system's current position relative to its zero reference.
	System E-STOP status	Identifies when E-STOP is engaged. The icon will flash when E-STOP is pressed.
	System tube heat status	Indicates the current X-Ray tube heat status. The values are color coded as follows: <div> <div>Blue</div> <div>0% - 25%</div> </div> <div> <div>Yellow</div> <div>26% - 50%</div> </div> <div> <div>Orange</div> <div>51% - 75%</div> </div> <div> <div>Red</div> <div>76% - 100%</div> </div>












Status bar icon	Status bar icon name	Status description						
	Scanner battery capacity status	<p>Indicates the remaining scanner battery percentage available. The capacity values are color coded as follows:</p> <table><tr><td>Green</td><td>100% - 51%</td></tr><tr><td>Yellow</td><td>50% - 25%</td></tr><tr><td>Red</td><td>24% - 0%</td></tr></table>	Green	100% - 51%	Yellow	50% - 25%	Red	24% - 0%
Green	100% - 51%							
Yellow	50% - 25%							
Red	24% - 0%							
	System air freshness status	<p>Indicates the air freshness status; it is recommended that an air calibration be performed:</p> <ul style="list-style-type: none">• Every eight (8) hours.• When the air freshness status falls below 50%.• If the scanner is moved to an area with a dramatic change in humidity and/or temperature. <p>The calibration status values are color coded as follows:</p> <table><tr><td>Green</td><td>100% - 51%</td></tr><tr><td>Yellow</td><td>50% - 25%</td></tr><tr><td>Orange</td><td>24% - 0%</td></tr></table> <p>After calibration it returns to 100%.</p>	Green	100% - 51%	Yellow	50% - 25%	Orange	24% - 0%
Green	100% - 51%							
Yellow	50% - 25%							
Orange	24% - 0%							

Status bar icon	Status bar icon name	Status description						
	Workstation battery capacity status	<p>Indicates the remaining workstation battery capacity available. The capacity values are color coded as follows:</p> <table><tr><td>Green</td><td>100% - 21%</td></tr><tr><td>Yellow</td><td>20% - 11%</td></tr><tr><td>Red</td><td>10% - 0%</td></tr></table> <p>You will be prompted to plug the workstation into an outlet to charge if the battery capacity is low; a scan cannot complete when the battery capacity is 10% or lower.</p> <p>When the workstation reaches the red capacity range, the system will shut down. A message informs you that the system will shut down due to low battery.</p> <p>The lightning bolt icon signifies that the workstation is currently charging and goes away when unplugged.</p>	Green	100% - 21%	Yellow	20% - 11%	Red	10% - 0%
Green	100% - 21%							
Yellow	20% - 11%							
Red	10% - 0%							
	Image storage space status	<p>Indicates the available disk space for image storage. The available space values are color coded as follows:</p> <table><tr><td>Green</td><td>100% - 51%</td></tr><tr><td>Yellow</td><td>50% - 20%</td></tr><tr><td>Red</td><td>19% - 0%</td></tr></table>	Green	100% - 51%	Yellow	50% - 20%	Red	19% - 0%
Green	100% - 51%							
Yellow	50% - 20%							
Red	19% - 0%							

The system changes states as it performs different actions. The following table indicates what state the system is in and the colored orb that correlates to that state.

Table 25: System state orbs

Orb	Color	State
	Dark gray	The system is in an unknown state.
	Light gray	The system is powering up or down.

Orb	Color	State
	Dark purple	The system is busy.
	Purple	The system is completing air calibration.
	Light purple	The system is archiving.
	Blue	The system is idle.
	Green	The system is ready to perform a scan.
	Light yellow	The system is planning.
	Dark yellow	The system is preparing.
	Light orange	The system is reconstructing.
	Dark orange	The system is scanning.
	Pink	The system is not ready.
	Red	The system is in fault.

The workstation tabs

To perform a patient examination, you will use the following five tabs on the workstation:



Figure 65: Workstation tabs to perform a patient examination

The tabs include active tabs which will be green, and inactive tabs which will appear gray. The active tabs are **Patient Registration**, **Post Reconstruction**, and **Patient Browser**. The **Acquisition** and **Viewing** tabs require additional steps to be performed before they become active. The following actions are available in each tab:

Patient Registration	Allows you to register a patient either manually or from the hospital's database sites.
-----------------------------	---

Acquisition	Allows you to select a protocol and perform the examination. This tab is inactive until a patient is registered.
Post Reconstruction	Allows you to manipulate raw data in different parameters and settings after your scan is completed.
Patient Browser	Allows you to view, manipulate, and archive scans already performed.
Viewing	Allows you to view patient images. . This tab is inactive until a study is loaded from Patient Browser.

In the following chapters, you will learn how to perform necessary steps to conduct a patient examination and learn how to manipulate and store the data you acquire.

Chapter 5 System and User Configuration and Setup

A user with administrative privileges must set up the BodyTom Elite system configurations for other users. System configuration is used to set up the scanner to meet site-specific needs. Most windows contain self-explanatory instructions and refer to elements that are known to the administrative user with radiological education and training.

Configuration includes setting up user permissions, or access privileges to manage other users, as well as servers, what is available through settings, presets, and other preferences for the use of the system at a site. Many system configurations are permissible to the administrator, *only*. While other configurations are permissible to users without administrator access.

Note You must have administrative access privileges and be logged in as an administrator to set configurations for the site.

Incorrect changes to the system configuration may make the system inoperative.

The following table shows the **System Configuration** options that appear when you click **Customize > System** and provides a brief description of each.

Table 26: System configuration tabs

Tab name	Description
General Settings	Allows the administrator to set system configurations such as hospital name, wireless settings, and dose report settings.
User Accounts	Allows the administrator to create and edit user accounts and permissions.
DICOM Servers	Allows the administrator to set up DICOM servers for archiving, such as PACS or HIS/RIS.
DICOM Setting	Allows the administrator to view and configure DICOM tags for HIS/RIS, MPPS, and Patient, Study, Series, and Image modules.
Audio Configuration	Allows the administrator to upload default audio files with protocols; also lets the administrator record, play, and remove audio files.
Dose Configuration	Allows the administrator to set up dose notifications, dose alerts, and configure dose limits for specific scans.
Windowing Presets	Allows the administrator to view and modify kernel and window width and window level presets.

Tab name	Description
Audit Trail Viewer	Allows the administrator to view and log all changes as well as actions in the system, including user logins, patient registrations, and series updates.
Image Orientation	Allows the administrator to view and modify how images are oriented in the system.
Filter Kernels	Allows the user to activate custom kernel options for both Axial and Helical scans.

Setting user accounts

Only a user with administrative access can update another user's account, add a new user, delete a user, and lock or unlock a user's access.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **User Accounts** tab.




Figure 66: User Accounts tab

3. For the **User ID** field, enter the User ID name.
4. For the **User Level** field, enter one of the following user levels:

Administrator	Full access to the system and its configurations. Can create protocols, User ID's, and passwords, as well as access all functions of the system.
----------------------	--

Limited operator	Modified access to the system. Users with Limited access can modify protocols during scanning but cannot create and save protocols; has no access to system configuration.
Restricted operator	Users with Restricted access can scan with the system but are unable to make any changes to protocol parameters while scanning, they also have no access to system configuration.

5. For the **Last Name** field, enter the user's last name.
6. For the **Enter First Name** field, enter the user's first name.
7. For the **Password** field, enter the user's password.

Note The password must contain 8 to 12 characters, and must include one number, one symbol and one letter.

8. For the **Verify Password** field, re-enter the user's password to confirm the password.

The screenshot shows the 'System Configuration' dialog box with the 'User Accounts' tab selected. On the left, there is a list of user levels: Administrator, Restricted Operator, and Limited Operator. The right side of the dialog contains several input fields: 'User ID' (filled with 'Administrator'), 'User Level' (a dropdown menu showing 'Administrator'), 'Last Name' (filled with 'System'), 'First Name' (filled with 'Administrator'), 'Password' (masked with dots), and 'Verify Password' (masked with dots). At the bottom of the dialog, there are buttons for 'Save', 'Update', 'Unlock', 'Delete', and 'Close'.

Figure 67: User account fields filled in

9. Click the **Save** button.
The user is added to the list.

The **Save Aborted** popup appears if your password does not meet the rule for passwords. If this is the case, return to the step above, and fulfill the password rule.




Figure 68: Save aborted popup message – Password requirements

10. Click the **Close** button to exit.

Setting or updating the user's information

1. Click **Customize > System**, from the main menu.
The **System Configuration** dialog box appears.
2. Click the **User Accounts** tab.
3. Select a user from the list of users.




Figure 69: List of users

4. Modify the user's information; for example, password, or user's first and last name.

Note The password must be 8-12 characters and include at least a letter, number, and symbol.

The user's changes take effect after clicking the **Update** button. It is recommended that you log off and log back on and check that the password is working.

5. Click the **Update** button to keep the change(s).

The **Update Aborted** popup appears if your password does not meet the rule for passwords. If this is the case, return to the step above, and fulfill the password rule.




Figure 70: Update Aborted popup message – Password requirements

6. Click the **Close** button to exit.

Unlocking the user

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **User Accounts** tab.
3. Select the user to unlock from the list of users in the panel.




Figure 71: List of users not selected

4. Click the **Unlock** button.

The user's changes take effect after clicking the **Update** button. It is recommended that you log off and log back on and check that the password is working.

Deleting the user

Note The administrator user cannot be deleted.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **User Accounts** tab.
3. Select the user to delete from the list of users.




Figure 72: List of all available users

4. Click the **Delete** button.
The **Action Succeeded** popup box appears.
5. Click the **Ok** button in the **Action Succeeded** popup.
6. Click the **Close** button to exit.

Modifying the order of the user in the accounts list

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **User Accounts** tab.
3. Select the user order to modify from the list of users in the panel.
4. Click the **Down** arrow to move the user down the list.




Figure 73: Down arrow

5. Click the **Up** arrow to move the user up the list.




Figure 74: Up arrow

The **Up** arrow will not activate until you move down the list of users.

6. Click the **Save** button under the user list to keep the new user list order.




Figure 75: Save button for list order

7. Click the **Save** button next to the **Update** button to keep other changes.
8. Click the **Close** button to exit.

Assigning general settings

Perform the following to set how the hospital, workstation's application, scanner, and remote support are configured.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **General Settings** tab.




Figure 76: General Settings tab

3. The following tabs are available to the administrator, *only*:

Hospital Setup	Sets up hospital information specific to the site.
Application Setup	Sets up application information the user will experience.
Scanner Setup	Sets up scanner IP address information.
Remote Support Setup	Sets up IP address information to allow NeuroLogica remote support.

The following sections provide detailed information related to the options available in the **General Settings** tabs.

Hospital Setup sub tab

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **General Settings** tab.
3. Click the **Hospital Setup** sub tab.




Figure 77: General Settings > Hospital Setup subtab

4. Enter the institution name in the **Institution Name** field.

Note The name appears on all images.

5. Enter the department name in the **Department Name** field.
6. Enter the institution address in the **Institution Address** field.
7. Click the following options that are applicable:
 - When **Auto Store** is selected and you Finalize a scan, the system will automatically send the images to the **Default PACS** server defined on the **DICOM Servers** tab. (see **System Configuration > DICOM Servers > Servers**). The status of your export will appear in the bottom right of the **Patient Browser** screen.
 - Click the **Use Device AE in Server Setting** option to apply the BodyTom Elite AE title tag to DICOM when the operator archives the image data to PACS.
 - If enabled when sending images to **PACS**, this option will include the BodyTom Elite AE setting as a **DICOM** tag.
8. Enter the system name (for example BodyTom Elite) in the **Station AE Title** field.
9. Enter the name you want to call the archive option in the **“Push To” Text** box.

10. Click the **Save** button to keep your changes.
The **Save Successful** popup appears.
11. Click the **Ok** button in the **Save Successful** popup.
12. Click the **Close** button to exit.

Application Setup subtab

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **General Settings** tab.
3. Click the **Application Setup** sub tab.




Figure 78: General Settings > Application Setup subtab

4. Click and/or enter information for the following options that apply:

Auto Lock	The Auto Lock feature can be used to automatically lock the workstation screen after a user defined Idle Time is reached.
Automatically Show Store/Print Queue	Allows the user to automatically display the Store/Print Queue status when images are set to archive to a network device.
Show Dose Report	Displays the Dose Report on the screen when the Finalize button is clicked.

	A dose report will not be generated until the operator clicks the Finalize button on the Acquisition tab.
Automatically Generate Dose SR	Generates a Dose SR (Structured Report) along with the dose report when the Finalize button is clicked.
Show Stored Query Results	Displays the Stored Results at the bottom of Patient Registration .
Prompt for Auto Reconstruction	Allows the user to automatically start any additional reconstructions added to a scan series after the exam is completed.
Allow Manual Patient Registration	Allows the user to manually register a patient.
Default Pediatric Patient Type	If selected, the Protocol Manager will default to the pediatric protocols.
Change Password Interval (Days)	Sets the number of days before a password change is required.

5. Click the **Save** button to keep your changes.
6. Click the **Close** button to exit.

Scanner Setup subtab

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **General Settings** tab.
3. Click the **Scanner Setup** sub tab.



Figure 79: General Settings > Scanner Setup subtab

4. Click and/or enter information for the following options that apply:

Scanner IP address	Sets the scanner's IP address in the field(s).
Point to Point Wireless	<p>Sets up wireless information regarding the connection from the workstation to the scanner.</p> <p>For Scanner Network Adapter, enter the adaptor, for example, Wireless to Scanner.</p> <p>For Scanner Wireless Profile, enter the wireless identifier in the field.</p>

5. Click the **Save** button to keep your changes.
6. Click the **Close** button to exit.

Remote Support Setup subtab

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **General Settings** tab.
3. Click the **Remote Support Setup** subtab.




Figure 80: General Settings > Remote Support Setup subtab

4. Click the **Remote Support Network Connection** dropdown to select one of the following network connections:
 - Wireless to Scanner
 - Wired to Network
5. Click the **Use Static IP** option to enter the **Static IP Setup** data for the following:
 - System IP Address
 - System Subnet
 - System Gateway
6. Enter the **DHCP Setup** to enter IPs for the following:
 - DNS Server (Preferred)
 - DNS Server (Alternate)
7. Click the **Save** button to keep your changes.
8. Click the **Close** button to exit.

Managing DICOM servers

Digital Imaging Communication in Medicine is the definition of the acronym **DICOM**. **DICOM** servers are used to export images from the scanner. The **System Configuration > DICOM Servers** tab allows a user with administrative rights to access all the **DICOM** devices connected to the scanner.

See **DICOM** standards on the **NEMA.org** website for a full description of settings and actions that are available.

Note You must have administrative privileges and be logged in as an administrator to access and modify DICOM servers.

Incorrect changes to the DICOM servers may make the system inoperative.

DICOM servers are set up by the **field-service engineer** and the appropriate IT person at the hospital.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **DICOM Servers** tab.




Figure 81: DICOM Servers tab

The following tabs appear:

Servers	Lists existing servers based on type: Store and Worklist Store: Identifies a storage server. Worklist: Identifies servers in a database you can query from. Also displays server details and options, with controls for saving, updating, deleting, and echoing servers.
PACS List	Displays a list of PACS by Server Name , Type , and In List – to send to by default.
Options	Displays controls for PACS Options and HIS/RIS Options .



Figure 82: DICOM Servers tabs

- Go to the following sections to assign specific actions to the **DICOM** server.

Assigning a server as a store or worklist server

- Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
- Click the **DICOM Servers** tab.
The **Servers** tab is active.



Figure 83: DICOM Servers > Servers tabs

3. Click one of the following options:

Store	A storage server, typically a PACS server that archives images and patient information. The images and data can later be imported from the server to another system or the same system. It is partly a backup and partly a waypoint for transferring data from the system the scans were acquired on to another system for viewing.
Worklist	A database of patient information that can be queried to generate a list of patients based on name, modality, procedure date, and other variables. These patients are usually hosted on a server within the facility and the information can be imported from the server when trying to acquire all a patient's information before a procedure is run. This eliminates the need to manually enter patient information at the time of the procedure.

4. Under **Details**, enter the server name in the **Server Name** text box.
5. For the **HIS/RIS** and **PACS Network Connection**, click the dropdown to identify the following:
 - Wireless to Scanner
 - Wired to Network
6. Click the **Point to Point Wireless** check box if it applies.

7. Click the **Wireless Profile** dropdown to select the appropriate profile.
8. In the **Server AE Title** text box, enter the server AE title.
9. In the **Device AE Title** text box, enter the BodyTom Elite AE title.
10. In the **Server IP Address** text boxes, enter the server IP address.
11. In the **Port** text box, enter the port identifier.
12. In the **Connection Timeout** text box, enter the number of seconds before the connection timeout is activated.
13. Click the **Use Static IP** option to identify the following:
 - System IP Address
 - System Subnet
 - System Gateway
14. To set the server as a default server, click the **Default Server** check box.
15. To set the server as the default surgical navigation server, click the **Default Navigation** check box.
16. To enable the storage commitment, click the **Enable Storage Commit** check box to send a message back to system that confirms the storage was successful; it is an extra confirmation from **PACS** that the images were received.
17. To gather responses, click the **Listen for Responses** check box.
For every image that is sent, the system will wait for acknowledgement before sending the next image.
18. To gather 12-bit images, click the **12-bit Images** check box.
19. Click the **Save** button, to keep your work.



Figure 84: Action Succeeded popup message – Server saved
The new server should appear in the list box to the left.

20. Click the **Ok** button.
21. Click the **Close** button to exit.

Modifying a server

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **DICOM Servers** tab.
3. Click server type: **Store** or **Worklist**.
4. Select a server from the list.
5. Modify the parameters.
See "Assigning a server as a store or worklist server" on page 137.
6. When all your changes are made, click the **Update** button.
A message appears that explains the update was successful and includes the update(s).



Figure 85: Action Succeeded popup message – Server updated

7. Click the **Ok** button.
8. Click the **Close** button to exit.

Echoing a server

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **DICOM Servers** tab.
3. Click a server type: **Store** or **Worklist**.
4. Select the server to echo from the list.
5. Click the **Echo** button.
The status of the server appears.



Figure 86: Echo Successful and Echo Failed popups

6. Click the **Ok** button.
If the echo was unsuccessful, determine why and repeat step 5 until you are successful.
7. Click the **Close** button to exit.

Deleting a server

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **DICOM Servers** tab.
3. Click a server type: **Store** or **Worklist**.
4. When the server you want to remove is available, click the server from this list.

- Click the **Delete** button.
The server is removed from the list; the **Action Succeeded** popup appears.
- Click the **Ok** button.
- Click the **Close** button to exit.

Moving a server up and down the server list

- Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
- Click the **DICOM Servers** tab.
- Click a server type: **Store** or **Worklist**.
- Select the server to move up or down the list.
- Click the **Up** arrow to move the server up the list; click the **Down** arrow to move the server down the list.




Figure 87: Up and Down arrows to move up and down server list

- Click the **Save** button to save the server order.




Figure 88: Save button

The **Action Succeeded** popup appears.

7. Click the **Ok** button.
8. Click the **Close** button to exit.

Saving DICOM servers to a PACS list

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **DICOM Servers** tab.
3. Click the **PACS List** tab to view available servers.



Figure 89: DICOM Servers > PACS List tab

4. Double-click the light-gray checkmark under **In List**. Each checkmark adds the server to the **PACS** listing. The checkmark turns green when active.
5. Click the **Save** button. The **PACS List Saved** popup appears.



Figure 90: PACS List Saved popup message – PACS saved

6. Click the **Ok** button.
7. Click the **Close** button to exit.

Selecting PACS options

1. Click **Customize > System** from the main menu. The **System Configuration** dialog box appears.

2. Click the **DICOM Servers** tab.
3. Click the **Options** tab.





Figure 91: DICOM Servers > Options tab

4. Under **PACS Options**, click the **Enable Storage Commit** to verify that patient information and data archived to the **PACS** server was received.
 - If **Enable Storage Commit** is not selected, it is assumed and accepted that the data was received.
 - If **Enable Storage Commit** is selected, the workstation sends a request to the PACS server to verify that the data was received.
5. If the **Enable Storage Commit** check box is selected, perform the following:
 - Enter the appropriate title in the System Called AE Title text box.
 - Enter the port identifier in the Port text box.
 - Enter the number of seconds before a connection timeout in the Connection Timeout (secs) text box.
 - Enter 30 seconds or 60 seconds.
6. Under **Purge Options**, make selections based on the understanding that any studies archived to the **PACS** server are deleted from the workstation on a regular basis, depending on the selected interval; identify the following:

Purge Interval	Select one of the following from the dropdown: Everyday , Weekly , or Monthly . When you select the first Purge Interval dropdown and select Weekly or Monthly , the inactive dropdown is active to let you select the day of the week or the first of the month.
Time	Changes the hour, minute, and second interval; use the buttons to increase those time elements.  Figure 92: Time (increase and decrease time) arrows

7. Under the **HIS/RIS Options**, enter the maximum number of results sent back from a query worklist in the **Query Max** text box.
There is no maximum limit.
8. Click the **Send MPPS** check box to apply a service that allows a modality to better coordinate with image storage servers by giving the server a list of objects to send before or while sending such objects.
9. Under **Import Options**, click the **Use C-GET** check box to pull information from a **PACS** server when importing *from* the server (as opposed to archiving to it).

The administrator sets this to pull from **PACS** from anywhere, so the machine does not have to be set up as a reliable destination on the **PACS** machine. **PACS**, typically needs to equate a computer's IP address with an AE title; however, **C-GET** accepts that the calling IP is a legitimate device.

The NeuroLogica BodyTom Elite scanner automatically uses **C-Move** when importing from **PACS**. If the operator wants to use **C-GET** instead, the user can select **C-Get**.

10. Click the **Save** button.
The **PACS List Saved** popup appears.



Figure 93: PACS List Saved popup

11. Click the **Ok** button.
12. Click the **Close** button to exit.

Assigning DICOM settings

DICOM settings include many kinds of settings. The administrator can add or remove optional information to be displayed using actions described in this section.

See **DICOM** standards on the **NEMA.org** website for a full list and description of **DICOM** tags.

Note You must have administrative privileges and be logged in as an administrator to access and modify DICOM settings.

Incorrect changes to the DICOM settings may make the system inoperative.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **DICOM Settings** tab.




Figure 94: DICOM Settings tabs (six)

- Click the **HIS/RIS Query** tab to select the types of **HIS/RIS** query results the user will see when performing a **HIS/RIS** query.
See "Selecting PACS options" and go to step 6, for more information.

Note Even numbered **DICOM** tags are public **DICOM** tags as per the **DICOM** standard. Odd numbered **DICOM** tags are vendor specific.

Green checkmarks are optional **DICOM** tags and orange checkmarks are required per the **DICOM** standard and cannot be modified.




Figure 95: DICOM Settings > HIS/RIS Query

- Click the **Modality Performed Procedure Step (MPPS)** tab to select the types of **MPPS** information the user will see.




Figure 96: DICOM Settings > MPPS

- Click the **Patient Module** tab to select the types of **Patient Module** information the user will see.



Figure 97: DICOM Settings > Patient Module

- Click the **Study Module** tab to select the types of **Study Module** information the user will see.




Figure 98: DICOM Settings > Study Module

- Click the **Series Module** tab to select the types of **Series Module** information the user will see.



Figure 99: DICOM Settings > Series Module

- Click the **Image Module** tab to select the types of **Image Module** information the user will see.




Figure 100: DICOM Settings > Image Module

- Click the **Save** button to save your changes.
- Click the **Close** button to exit.

Assigning audio configuration

Default audio files are installed on the workstation. Audio files can be attached to protocols and sent to the scanner. Each audio file has an indication if it has been sent to the scanner.

Note You must have administrative privileges and be logged in as an administrator to access and modify audio configurations.

Incorrect changes to the audio configurations may make the system inoperative.

Finding and listening to audio files

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Audio Configuration** tab.




Figure 101: Audio Configuration tab

3. Review the audio files that exist on the **Workstation**.




Figure 102: Audio files list

4. To listen to an audio file, select the name from the workstation list and click the Play button.
5. To exit the Audio Configuration, click the **Close** button.

Recording and saving an audio file

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Audio Configuration** tab.
The available audio files appear under **Workstation**.




Figure 103: Audio files list

3. Enter the name of your new audio file in the File name text box.
4. Click the **Record** button.
5. Record your audio file.
6. Press the **Stop** button.
7. To review, highlight your new recording and press the **Play** button.




Figure 104: New audio file

8. When you like your recording, highlight the file and press the **Transmit** button to copy the file into the audio folder for your scanner protocols.
9. Click the **Close** button to exit.

Transmitting an audio file

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the Audio Configuration tab.
The available audio files appear under **Workstation**.
3. Select the audio files under **Workstation** audio files to transfer to the scanner.
4. Click the **Transmit** button.




Figure 105: Audio files transmitted to save to the scanner

5. Click the **Close** button to exit.

Deleting an audio file

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the Audio Configuration tab.

The available audio files appear under **Workstation**.

3. Select the audio file you want to delete from the list.
4. Click the **Delete** button.
5. Click the **Close** button to exit.

Assigning dose report

The **dose report** is created at the end of the scan and can be customized to include **DICOM** specific tags.

Note You must have administrative privileges and be logged in as an administrator to access and modify dose report settings.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Dose Configuration** tab.
3. Click the **Dose Report** tab.



Figure 106: Dose Report tab

4. Select the **DICOM** tags you want to see in the **dose report**.
A **dose report** is generated after the exam is finalized; the black area includes dose report information like the following.



Figure 107: Generated dose report

5. Click the **Save** button.
6. Click the **Close** button to exit.

Applying dose configuration

Note You must have administrative privileges and be logged in as an administrator to access and modify dose configurations.

Incorrect changes to dose configuration settings may make the system inoperative.

Note You can check the **Audit Trail** to review the audit log that details what dose limit was removed, by whom, and the date and time it took place.

Dose configuration consists of both **Dose Notifications** and **Dose Alerts**.

Dose Notification	Notifies the user when the planned CTDI _{vol} and/or DLP value of a single series will exceed the defined value.
Dose Alert	Notifies the user when the planned CTDI _{vol} and/or DLP value from the combination of all planned series will exceed the defined value set in System Configuration . Dose Alerts represent a value which would be well above an institutions established CTDI/DLP range to the given examination and warrant a more stringent review and consideration before proceeding.

Setting Dose Check

See Appendix A on page 389 for information on protocols, CTDI_{vol}, and DLP.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Dose Configuration** tab.
3. Click the **Dose Check** tab.




Figure 108: Dose Configuration > Dose Check

4. Click one of the following **Dose Check Type** options:

Dose Notification	Notifies the user when a pre-defined CTDI _{vol} or DLP value will be exceeded on a series-by-series basis.
Dose Alert	Notifies the user when a pre-defined CTDI _{vol} or DLP value will be exceeded from a combination of all planned series or scans.

Note The default **Dose Alerts** which are set at 1000mGy CTDI and 2000mGy*cm DLP are designed to prevent the patient from receiving any possible deterministic effects due to excess dose. However, the system allows these values to be modified by the user. Any modifications to the **Dose Alerts** should be done by qualified medical personnel.

5. Click one from the following **Scan Type** options.

All	Identifies all scan types.
Axial	Identifies only Axial scan types.
Helical	Identifies only Helical scan types.
Dynamic	Identifies only Dynamic scan types.

6. Define the **Dose Limit** by entering the following:
- Enter the CTDI_{vol} (mGy) value in the text box.
 - Enter the DLP (mGy.cm) value in the text box.

7. Click the **Save** button.
The **Save Successful** popup appears.




Figure 109: Save Successful popup – Dose Check successfully saved

8. Click the **Ok** button.
9. Click the **Close** button to exit.

Assigning Dose Configuration to a patient protocol

Dose Configuration limits are used to prevent users from selecting kV or mA values that are not appropriate for the given patient types, such as pediatrics etc.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Dose Configuration** tab.
3. Click the **Dose Configuration** subtab.





Figure 110: Dose Configuration > Dose Configuration for adult and pediatric

4. Click one of the following:

Adult	Selecting Adult shows the pre-defined adult protocols, stored by anatomical area.
Pediatric	Selecting Pediatric shows the pre-defined pediatric protocols, stored by anatomical area.
Trauma 	The Trauma orb can be used to store protocols commonly used for emergency scans.

5. Click the colored orb that marks the anatomical region to apply the dose to.




Figure 111: Anatomical orbs

6. Click a scan type from the following list:

Axial	Identifies only Axial scan types.
Helical	Identifies only Helical scan types.
Dynamic	Identifies only Dynamic scan types.

7. Enter a description for the **Dose Configuration** in the **Description** text box.
8. For **Pediatrics** enter the **Minimum** and **Maximum Weight** and **Length** information.
9. Under the **Max Power** settings, click the **kV** dropdown and select the maximum allowed kV.
10. Click the **mA** dropdown and select the maximum allowed mA.




Figure 112: Pediatric Dose Configuration Parameters

11. Click the **Save** button to save your work.
If the level overlaps an existing level, you are prompted to adjust.



Figure 113: Invalid Parameter popup message – Dose setting kV already exists

If the save is successful, the **Save Successful** popup appears.




Figure 114: Save Successful popup message – Maximum dose saved

12. Click the **Ok** button.
13. Under the **Saved Doses** list box, check that your dose configuration appears, if so, go to the next step.
14. Click the **Close** button to exit.

Updating saved dose

To modify a saved **Dose Configuration**.

See “Setting Dose Check” page 158 and/or “Assigning Dose Configuration to a patient protocol” on page 159.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Dose Configuration** tab.
3. Select the desired **Saved Dose** from the **Saved Doses** list.




Figure 115: Saved Doses List

4. Modifying values causes the **Update** button to become active.
5. Click the **Update** button.




Figure 116: Save Successful popup message – Maximum dose saved

6. Click the **Ok** button.
7. Click the **Save** button to save your work.

Note If the level overlaps an existing level, you are prompted to adjust.

8. Click the **Close** button to exit.

Deleting a saved dose limit

To remove a saved **Dose Configuration**.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Dose Configuration** tab.




Figure 117: Dose Configuration > Dose Check tab

3. Click an already-saved dose, under the **Saved Doses** list.
4. Perform one of the following:
 - To delete a saved dose from the Saved Doses list and clear the restriction, select the dose, and click the Delete button.
 - To delete all the saved doses in the Saved Doses list and clear all restrictions saved, click the Delete All button, which returns all settings for that selection to the maximum scanner default.

Note If there are no saved doses or limits, the operator will be able to scan using the maximum 140kV and 300mA available on the scanner.

5. Click the **Save** button.
The save success message appears and, because the **CTDIvol (mGy)** and **DLP (mGy.cm)** are empty, there is no longer a limit applied.
The **Save Successful** popup appears.




Figure 118: Save popup message – Maximum dose saved

6. Click the **Ok** button.
7. Click the **Close** button to exit.

Applying Windowing Presets

Windowing presets allow you to define window width and center presets for specific anatomical locations as well as specific reconstruction kernel presets. An Administrative User can delete or update the default Window Presets as well as create new Windowing Presets.

Note You must have administrative privileges and be logged in as an administrator to access and modify the windowing presets.

Editing kernel presets

Note Kernel presets are pre-installed in the system; kernel presets can be set and modified.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Windowing Presets** tab.




Figure 119: Windowing Preset tab

3. Click the **Kernel Presets** tab.




Figure 120: Windowing Presets > Kernel Presets tab

4. Click the **Scan Type** dropdown to select one of the following scan types:

- Axial
- Helical

- Click the **Sharpness** dropdown to select a sharpness from the list. **Sharpness** is the reconstruction algorithm, and the available **Sharpness** values are based on the scan type.




Figure 121: Sharpness dropdown

- Enter the Window Width in the **Width** text box. **Window Width** describes the range of Hounsfield units, or shades of gray, displayed across the image. The **Window Width** controls the contrast of the image. Low Hounsfield numbers below the range are displayed as black, while High Hounsfield numbers above the range are displayed as white.
- Enter the Window Center in the **Center** text box. **Window Center** describes the Hounsfield number in the center of the Window Width. **Window Center** controls the brightness or density of the image.
- Click the **Save** button to save your work. The **Action Succeeded** popup appears.



Figure 122: Action Succeeded popup message – Preset saved

9. Click the **Ok** button.
10. Click the **Close** button to exit.

Setting Window Presets

Window presets allow you to define **Window Width** and **Window Center** presets for specific anatomical locations, such as bone, brain, lung, and soft tissue.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Windowing Presets** tab.
3. Click the **Window Presets** subtab.




Figure 123: Window Presets tab

4. Enter the name of the window preset in the **Name** text box.

The screenshot shows the 'System Configuration' application with the 'Windowing Presets' tab selected. The 'Window Presets' sub-tab is active, displaying a list of kernel presets on the left: Abdomen, Angio, Bone, Brain, Chest, and Lungs. On the right, there are three text input fields labeled 'Name', 'Width', and 'Center'. A yellow arrow points to the 'Name' field. At the bottom, there are buttons for 'Save', 'Update', 'Delete', and 'Close'. A 'Save' button is also located at the bottom left of the kernel list.

Figure 124: Window Presets > Name

5. Enter the width of the window preset in the **Width** text box.

This screenshot is identical to Figure 124, showing the 'Window Presets' dialog box. In this view, a yellow arrow points to the 'Width' text input field.

Figure 125: Window Presets > Width

6. Enter the center for the window preset in the **Center** text box.




Figure 126: Window Presets > Center

7. Click the **Save** button to save your work.
The **Action Succeeded** popup appears.




Figure 127: Action Succeeded popup message – Preset saved

8. Click the **Ok** button.
9. Click the **Close** button to exit.

Editing a window preset

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Windowing Presets** tab.

3. Click the **Window Presets** subtab.
4. Click a preset that exists in the **Window Presets** listing.




Figure 128: Listing update

5. To edit the preset, make your changes in the **Name**, **Width**, and/or **Center** text boxes.
6. Click the **Save** button to save your changes. The **Action Succeeded** popup appears.



Figure 129: Action Succeeded popup message – Preset saved

7. Click the **Ok** button.
8. Click the **Close** button to exit.

Deleting a preset

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Windowing Presets** tab.
3. Select the preset.
4. Click the **Delete** button.
The **Action Succeeded** popup appears.



Figure 130: Action Succeeded popup message – Preset deleted

5. Click the **Ok** button.
6. Click the **Save** button to exit.

Setting up the Audit Trail Viewer

The **Audit Trail Viewer** gives a user with administrative access the ability to view all activities performed by anyone logged into the system. This includes changes to protocols, deletion of images, as well as acknowledgement of alerts etc.

Note You must have administrative privileges and be logged in as an administrator to access the Audit Trail Viewer.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Audit Trail Viewer** tab.



Figure 131: Audit Trail Viewer tab

3. From the **Audit Trail Viewer** tab, use the calendar to select a date range to view.
 - To select a single date, position the mouse pointer in the top box and click the date on the calendar to find audits for that date.
 - To select a date range, click the desired start date on the calendar. This will automatically populate the top box of the **From** range. Click the desired end date on the calendar, which automatically fills the bottom box of the **From** range.



Figure 132: Adding a date or a date span

4. Click the **Audit Type** dropdown to select the type of audit you are searching.




Figure 133: Audit Trail Viewer > Audit Type dropdown

5. From the **User ID** dropdown, click the type of user to track.




Figure 134: Audit Trail Viewer > User ID dropdown

6. Click the **View** button to see the result of audits that met your criteria.




Figure 135: Audit results

7. Perform one of the following:
 - Click the **Clear** button to remove the audit results.
 - Click the **Export** button to export the audit result that you selected to the audit backup file on the system.
 - Click the **Export All** button to export the audit results to the audit backup file on the system.
8. Click the **Close** button to exit.

Setting image orientation

NeuroLogica describes patient orientation as if the viewer were looking towards the front of the gantry. In other words, if the patient is lying face up with their head in the gantry, the image orientation displays the patient's Right side on the Left side of the Viewer. If the patient's feet are going into the gantry, the image orientation displays the patient's Left side on the Left side of the Viewer.

Note You must have administrative privileges and be logged in as an administrator to modify image orientation settings.

Changes to image orientation settings will modify the displayed orientation markers on the images.

1. Click **Customize > System** from the main menu.

The **System Configuration** dialog box appears.

2. Click the **Image Orientation** tab.




Figure 136: Image Orientation tab

The top half of the **Image Orientation** screen shows the available patient orientations. The black image orientation square represents the viewing area and shows four different orientation markers: A = anterior, L = left, P = posterior and R = right. If you do not see letters in the image orientation box, select an image orientation from the list.

3. Select the appropriate orientation from the list.
For example, select, **Head First Supine**. In the figure below, the highlighted selection shows the current orientation in the **Current Flip Orientation** field, which is not changeable; however, the **New Flip Orientation** lets you change the orientation.
4. Click the **New Flip Orientation** dropdown to select one of the following new-flip orientations:
 - NONE
 - LEFT_AND_RIGHT
 - TOP_AND_BOTTOM
 - BOTH




Figure 137: Image Orientation > New Flip Orientation dropdown

5. Click the **Save** button to save changes.
6. The **Settings Saved** popup appears.



Figure 138: Settings Saved popup message – Image orientation settings saved

7. Click the **Ok** button.
8. Click the **Close** button to exit.

Setting Filter Kernels

Filter kernels allow you to activate custom kernel options for both **Axial** and **Helical** scans to control the sharpness and smoothness of the images.

1. Click **Customize > System** from the main menu.
The **System Configuration** dialog box appears.
2. Click the **Filter Kernels** tab.




Figure 139: Filter Kernels tab

3. Perform one of the following:
 - To add new Axial kernels, select the **Axial** radio button, then double-click the **Value** cell next to the desired Axial kernel.




Figure 140: Selected Axial kernel

To add new Helical kernels, select the **Helical** radio button, then double-click the **Value** cell next to the desired Helical kernel.




Figure 141: Selected Helical kernel

4. Click the **Save** button.
The **Success** popup appears.
5. Click the **Ok** button.
6. Click the **Close** button to exit.

User configuration

User configuration allows users with either **Administrator** or **Limited Operator** access to change the password for their own account.

Updating your user account

1. Click **Customize > User** from the main menu.
The **User Configuration** dialog box appears.

Note You must have administrative or limited operator privileges to access and modify user accounts.

Incorrect changes to user accounts may cause users to be unable to access the system.

The **Update Accounts** tab is the default selection.

2. Modify the following fields associated with your user account by entering relevant information:
 - Last Name
 - First Name
 - Password
 - Verify Password




Figure 142: Last Name, First Name, Password, and Verify Password fields

3. Click the **Update** button.

4. The **Update Succeeded** popup appears.



Figure 143: Update Succeeded popup message – Account updated

5. Click the **Ok** button.
6. Click the **Close** button to exit.

Applying column settings to HIS/RIS Query

Allows you to customize the columns of information that appear when viewing the **Hospital Information System (HIS)** and/or **Radiology Information Systems (RIS)** information that is queried.

All users can access user configuration and make changes to the column settings; however limited and restricted operators cannot make their changes to the column headings the default. Only users with administrative access can make column settings a default, using the **Make Default** option.

1. Click **Customize > User** from the main menu.
The **User Configuration** dialog box appears.
2. Click the **Column Settings** tab.
3. There may be no entries that appear, initially.

Note When an option is selected (for example **HIS/RIS Query** or **Patient Browser**), a table is created that lists the columns, along with a check box to indicate whether it will be displayed within the table. Required columns **cannot** be unchecked and are colored orange instead of the default green.

4. Click the **HIS/RIS Query** option.




Figure 144: Column Settings dialog box with HIS/RIS Query option

- Double click the Checkmark under **Displayed** column for the row you want to display or remove if it is already selected.




Figure 145: Column Settings with a selected query in HIS/RIS

- Notice that the **Up** and **Down** buttons are active when a row is selected.
- Click the **Up** or **Down** buttons to move the active selection up or down the list.




Figure 146: Column Settings with HIS/RIS Query option using Up and Dwn buttons

8. If you have administrative privileges, click the **Make Default** option to make the selected column display the default for all users.




Figure 147: Make Default option

9. Click the **Save** button to keep changes.
10. Click the **Close** button to exit.

Applying column settings to Patient Browser

Allows you to configure the columns of information seen in the **Patient Browser**.

1. Click **Customize > User** from the main menu.
The **User Configuration** dialog box appears.
2. Click the **Column Settings** tab.
3. Click the **Patient Browser** option.



Figure 148: Column Settings with Patient Browser option

4. Click one of the following options:

Patient/Study	Information that appears on the top portion of the Patient Browser that defines patient specific information.
Series	Information that appears on the lower portion of the Patient Browser that defines series specific information.

5. Double click the Checkmark under **Displayed** column for the row you want to display or remove if it is already selected.
6. Notice that the **Up** and **Dwn** buttons are active when a row is selected.
7. Click the **Up** or **DWN** buttons to move the active selection up or down the list.




Figure 149: Column Settings with Patient Browser Series option – using Up and Dwn buttons

8. If you have administrative privileges, click the **Make Default** option to make the selected column display the default for all users.




Figure 150: Make Default option

9. Click the **Save** button to keep changes.
10. Click the **Close** button to exit.

Viewing Scan Dosage Report

1. Click **Customize > User** from the main menu.
The **User Configuration** dialog box appears.
2. Click the **Scan Dosage Report** tab.



Figure 151: Scan Dosage Report tab

3. To view **Dosage Reports**, click the date or date range on the calendar.

Note If the date is left blank – all doses for all dates are retrieved.

4. Click the **Protocol** dropdown to select a protocol.

Note If the protocol is left blank – all doses for all protocols are retrieved.

5. Click the **mA Range** dropdown to select the mA range.
The **mA Range** default is 20-30; it can be changed after data is retrieved.




Figure 152: Date, Protocol, and mA Range filled

- Click the **View** button to display a graph showing dosages performed by the scanner using the selected filters.




Figure 153: Scan Dosage Report results

Note If you adjust the mA range, the graph displays those ranges within the retrieved data.

7. To clear the filters selected, click the **Clear** button.

Selecting a room for the BodyTom Elite

Selecting a room ensures that the correct calibration is loaded or used when the scanner has been calibrated in more than one location.

1. Click **Customize > Select Room**.
2. A list of the rooms available appear in the cascading menu.




Figure 154: Available rooms before moving the scanner

3. Click the room in which the scanner will be used.
4. Move scanner to the selected room.

Chapter 6 Protocol Manager

Note You must have administrative privileges and be logged in as an administrator, to access the Protocol Manager.












Protocol Manager allows a user with administrative privileges to create new protocols, modify existing protocols and delete protocols from the system. **Protocol Manager** provides three patient options: **Adult**, **Pediatric**, or **Trauma**  patient.

Table 27: Protocol Manager command buttons

Button	Action
	Allows you to create a new protocol.
	Allows you to create a new protocol from an existing protocol.
	Allows you to modify protocols.
	Deletes a saved protocol.
	Imports previously exported protocols into the workstation.
	Exports protocols to a media device.
	Closes the Protocol Manager dialog box.
 	Moves a protocol up or down the ordered list.
	Saves the order of the protocol list.

Note Different patient sizes and attenuations require specific protocols where technical parameters have been adjusted according to the physical characteristics of the patient. It is recommended that you select pediatric protocols when scanning children, rather than selecting adult protocols and modifying technique factors.

Creating a new protocol

1. Click **Tools >Protocol Manager** from the main menu.
The **Protocol Manager** dialog appears.



Figure 155: Protocol Manager for adult and pediatric

2. Click one of the following:

Adult	To create and/or scan with adult scan protocols, which are stored by anatomical location.
Pediatric	To create and/or scan with pediatric scan protocols, which are stored by anatomical location.
Trauma	The Trauma orb can be used to store protocols commonly used for emergency scans.
	By selecting either an Adult or Pediatric patient, the corresponding list of saved protocols becomes available.

Note Adult and pediatric protocol parameters are customized to meet your requirements in conjunction with local and nationally recognized published guidelines. These protocols **must be** approved by your facility physicist **before** the system's acceptance.

This scanner comes with reference protocols for pediatric patients that vary according to age or weight. These protocols are meant to be a guide for routine scanning. An individual patient may require higher or lower doses than in the reference protocol to achieve the diagnostic goals for the patient's medical condition. These reference pediatric protocols were established as a reduction from the adult protocols to achieve the same image signal to noise. Currently there is no agreement among scientific and medical organizations as to what the proper protocol is for pediatric patients. It is the responsibility of the user to decide whether these protocols are adequate for achieving their diagnostic goals on a patient-by-patient basis.

The user can modify any of the existing protocols on the machine and/or create new protocols as deemed necessary. To create these protocols, Administrator privileges are required.

3. Click the colored orb corresponding to the appropriate body part.




Figure 156: Adult and pediatric anatomical orbs, with Chest orb selected

Existing protocols in the selected Orb will appear in the **Protocol** list box as seen below. The **New** button will become active.




Figure 157: Adult and pediatric protocol lists

4. Click the **New** button to create a new protocol.
The **New Protocol** dialog box appears.




Figure 158: New Protocol dialog box

Enter **Protocol Information** in the text boxes:

Protocol Description	The name of the protocol as it will be displayed in protocol manager or when selecting a protocol to use for a scan.
Anatomical Reference	References the anatomy that will be scanned.
Default Study Description	This will appear as the DICOM image tag, this description also appears in PACS as a Study Description DICOM tag (0008,1010).

5. Under **Patient Position**, select one of the following scanning positions from the dropdown:

HFS	Head First-Supine
HFP	Head First-Prone
HFDL	Head First-Decubitus Left
HFDR	Head First-Decubitus Right
FFS	Feet First-Supine
FFP	Feet First-Prone
FFDL	Feet First-Decubitus Left
FFDR	Feet First-Decubitus Right

You can also move the rotating positional handles by hovering the mouse pointer over the handle and clicking to select a position.

- The arrows at the feet rotate the patient orientation from **Head First** to **Feet First**.
- The arrow above the patient rotates the patient orientation from **Supine** to **Prone** to **Decubitus**.

The **Patient Position** settings are identical whether it is for an adult, pediatric, or emergency patient.




Figure 159: Patient position handles

6. Under **Protocol's Series**, click the **New** button.
The **New Series** dialog box appears.

Figure 160: New Series dialog box

7. For **Scan Type**, select one of the following:

- Axial
- Helical
- Dynamic
- Reference
- Scout

Note For **Helical** scanning, the exposed area is extended by at least $\frac{1}{2}$ rotation to 1 full rotation at the start and end of the planned scan, based on the sharpness selected.

8. For **Scout Type**, select one of the following:

- AP
- PA
- Lateral

Scout Type is not available for **Axial**, **Helical**, **Dynamic**, and **Reference** scan modes.

9. For **kV** (scan voltage), select one of the following:

- 80 To set the scan kV to 80.
- 100 To set the scan kV to 100.
- 120 To set the scan kV to 120.
- 140 To set the scan kV to 140.

kV is not selectable when using the **Dynamic** and **Reference** scan modes.

See “Identifying load factors” on page 235.

10. For **mA** (scan current), select the desired value from the dropdown.
mA is not selectable when using the **Reference** scan mode.
The available **mA** range is 30 to 300 with an increment of 5.
11. For **Slice Thickness/Spacing**, select the desired value from the dropdown.
Slice Thickness/Spacing is not available for **Scout** and **Reference** scan modes.
12. For **Sharpness**, select the image reconstruction kernel from the following list of kernels:
 - Low Noise QA – Not for Clinical Use
 - Soft Tissue – Available for Axial Scan Type only
 - Soft Tissue – Head – Available for Helical Scan Type only
 - Soft Tissue – Abdomen – Available for Helical Scan Type only
 - Pos. Fossa/Vessel – Available for Axial Scan Type only
 - Sharp – Available for Axial Scan Type only
 - Bone (No AEC) – Available for Axial Scan Type only
 - Bone – Head (No AEC) – Available for Helical Scan Type only
 - Bone – Abdomen (No AEC) – Available for Helical Scan Type only
 - Sharp Lung (No AEC) – Available for Axial Scan Type only
 - High-Res QA (No AEC) – Not for Clinical UseSharpness is not selectable when using the **Reference** and **Scout** scan modes.

Note The **Low Noise QA** and **High-Res QA (No AEC)** options **should not** be used for clinical scanning.

13. For **Resolution**, which also refers to scan time, select one of the following options:
 - 1 Second(s)
 - 2 Second(s)**Resolution** is only available for **Axial Scan Types**.
14. For **Pitch**, which describes how fast the scanner is moving during one rotation of the x-ray tube, select one of the following options:
 - **0.4** where the scanner will move 16mm per second.
 - **0.8** where the scanner will move 32mm per second.**Pitch** is only available for **Helical Scan Types**.

15. For **Body Part Examined**, select the appropriate Body Part from the drop-down menu.
16. For **Window Width**, enter the range of CT numbers that are distributed over the viewable gray scale of the display device or film.
17. For **Window Center**, enter the CT number in the center of the viewable gray scale.
18. For **Description**, enter the desired study description.
19. For **Start Position**, enter the start scan position.
20. For **End Position**, enter the end scan position.
21. **Coverage** is a calculated value that automatically fills based on the **Start** and **End** position values.
22. For **Contrast**, enter the type of contrast given for example
23. For **Contrast Volume**, enter the amount of the contrast given.
Contrast is not available for **Reference** and **Scout** scan modes.
24. For **Delay**, enter the delay time that will occur after clicking the **START** button on the scanner control panel.
25. **Number of images** is a calculated value based on the Slice Thickness/Spacing and length of the scan.
26. **Scan Time** is a calculated value based on the protocol parameters selected. Scan time is affected by **Resolution**, **Pitch** and **Scan Length**.
27. For **CTDI_{vol} (mGy)**, if applicable, the calculated number appears here, depending on other selections.
CT Dose Index Volume (CTDI_{vol}) represents the dose for a specific scan protocol, which considers gaps and overlaps between the radiation-dose profiles from consecutive rotations of the x-ray source. The CTDI_{vol} is calculated differently for both the **Axial** and the **Helical** scan modes:
 - For Axial scan mode: $CTDI_{vol} = [(N \times T)/I] \times CTDI_w$
 - For Helical scan mode: $CTDI_{vol} = 1/pitch \times CTDI_w$
28. **Dose Length Product (DLP (mGy.cm))**, is the measure of ionizing radiation exposure during the entire acquisition of images.
Therefore, $DLP (mGy.cm) = CTDI_{vol} (mGy) \times irradiated\ length (cm)$.

29. Select the following options, if applicable. See Scanning with special features for more details.

Step & Shoot	Allows you to manually start the Axial scan acquisition from the workstation when scanning a patient who is unable to remain still.
Bolus Tracking	A CT angiography technique that allows you to monitor the administration of contrast to initiate the scan at peak contrast enhancement.
Enable AEC	Allows you to automatically adapt the tube current according to the patient's body habitus to achieve the specified image quality at the lowest possible dose.

30. To add a secondary reconstruction for the protocol, click the **New** button in the **Recons** section.

The **New Reconstruction** popup appears.

The screenshot shows a 'New Reconstruction' dialog box. It contains the following fields and controls from top to bottom:

- Description:** A text input field.
- Slice Thickness/Spacing:** A dropdown menu.
- Sharpness:** A dropdown menu.
- Window Width:** A text input field.
- Window Center:** A text input field.
- FOV Width:** A text input field.
- FOV Top Left X:** A text input field.
- FOV Top Left Y:** A text input field.
- Buttons:** 'Save', 'Reset', and 'Cancel' buttons at the bottom.

Figure 161: New Reconstruction popup

31. Complete the following in the **New Reconstruction** popup:
- Enter a description in the **Description** text box to identify the new reconstruction.
 - Click the **Slice Thickness/Spacing** to select a slice thickness and spacing.
 - Click the **Sharpness** dropdown to select a sharpness from the list.
 - Enter the window width in the **Window Width** text box.

- Enter the window center in the **Window Center** text box.
- If needed, enter the FOV width in the **FOV Width** text box.
- If needed, enter the FOV top left x location in the **FOV Top Left X** text box.
- If needed, enter the FOV top left y location in the **FOV Top Left Y** text box.

32. Perform one of the following:

- Click the **Save** button to save the reconstruction protocol to the list.
 - The dialog box closes, and your changes are added to the **Recons** area.
- Click the **Reset** button to reset the fields to their original data.
- Click the **Cancel** button to remove your changes and return to the previous dialog box.

33. Click the **Save** button on the **New Series** dialog box.

The 'New Series' dialog box contains the following parameters and controls:

- Series Parameters:**
 - Scan Type: Helical
 - Description: CHEST
 - FOV Width: [text box]
 - Scout Type: [text box]
 - Start Position: 0
 - FOV Top Left X: [text box]
 - kV: 120
 - End Position: 300
 - FOV Top Left Y: [text box]
 - mA: 175
 - Coverage: 300
 - Bolus Tracking: [checkbox]
 - Slice Thickness/Spacing: 1.25 x 0.625
 - Contrast: [text box]
 - Auto-Start: [radio]
 - Auto-Stop: [radio]
 - Test Bolus: [radio]
 - Sharpness: Soft Tissue - Abdomen
 - Contrast Volume: [text box]
 - Bolus Scan Time: [text box]
 - Threshold: [text box]
 - Resolution: 1 Second(s)
 - Delay: [text box]
 - AEC: [checkbox]
 - Enable AEC: [checkbox]
 - Pitch: 0.8
 - Number of Images: 480
 - Minimum mA: [text box]
 - Noise Level: [text box]
 - Maximum mA: [text box]
 - Body Part Examined: CHEST
 - Scan Time: 9.375
 - Recons: [list area]
 - Window Width: 400
 - CTDIvol (mGy): 16.91
 - Window Center: 40
 - DLP (mGy.cm): 507.3
- Buttons:**
 - Step & Shoot
 - Use Breathe Indicator Audio
 - Choose...
 - Save (highlighted)
 - Reset
 - Cancel
 - New
 - Edit
 - Remove

Figure 162: Edit Series dialog box

34. Repeat the steps 6 thru 33 to add additional scans to the protocol.

35. When all required series have been created click the **Save** button on the **New Protocol** dialog box.



Figure 163: Save New Protocol

36. Click the **Close** button to exit.






Figure 164: Close Button

Using Build From to create a new protocol

The **Build From** button is used in **Protocol Manager** when you want to create a new protocol from an existing protocol.

1. Click **Tools >Protocol Manager** from the main menu.
The **Protocol Manager** dialog appears.
2. Click one of the following:

Adult	To create and/or scan with adult scan protocols, which are stored by anatomical location.
--------------	---

Pediatric	To create and/or scan with pediatric scan protocols, which are stored by anatomical location.
Trauma 	The Trauma orb can be used to store protocols commonly used for emergency scans.
	By selecting either an Adult or Pediatric patient, the corresponding list of saved protocols becomes available.

Note Adult and pediatric protocol parameters are customized to meet your requirements in conjunction with local and nationally recognized published guidelines. These protocols **must be** approved by your facility physicist **before** the system's acceptance.

- Click the colored orb corresponding to the appropriate body part.




Figure 165: Anatomical orbs

- Click the protocol you will **Build From** in the **Protocol** list. The **Build From** button will become active.




Figure 166: Build from protocol selected

5. Click the **Build From** button.




Figure 167: Build From button

6. The **New Protocol** dialog box appears.




Figure 168: New Protocol dialog box

7. Modify the protocol parameters to meet your needs. Click the **Update** button on the **Edit Series** dialog box to save your changes.

Note Different patient sizes and attenuations require specific protocols where technical parameters have been adjusted according to the physical characteristics of the patient. It is recommended that you select pediatric protocols when scanning children, rather than selecting adult protocols and modifying technique factors.

See “Creating a new protocol” on page 190 to learn how the fields and options perform to make informed choices on what to change.

Note Be sure to assign the **Build From** protocol a new **Protocol Description** before you make your additional changes.

8. When all required series have been modified click the **Save** button on the **New Protocol** dialog box.




Figure 169: Build from save

9. Click the **Close** button to exit.






Figure 170: Build from close

Editing an Existing Protocol

The **Edit** button is used in **Protocol Manager** when you want to modify the parameters of an existing protocol.

1. Click **Tools >Protocol Manager** from the main menu.
The **Protocol Manager** dialog appears.
2. Click one of the following:

Adult	To create and/or scan with adult scan protocols, which are stored by anatomical location.
Pediatric	To create and/or scan with pediatric scan protocols, which are stored by anatomical location.
Trauma 	The Trauma orb can be used to store protocols commonly used for emergency scans.
	By selecting either an Adult or Pediatric patient, the corresponding list of saved protocols becomes available.

Note Adult and pediatric protocol parameters are customized to meet your requirements in conjunction with local and nationally recognized published guidelines. These protocols **must be** approved by your facility physicist **before** the system's acceptance.

- Click the colored orb corresponding to the appropriate body part.




Figure 171: Edit protocol orbs

- Click the protocol you wish to **Edit** in the **Protocol** list. The **Edit** button will become active.




Figure 172: Edit protocol selected

5. Click the **Edit** button.




Figure 173: Edit button

6. The **Edit Protocol** dialog box appears.




Figure 174: Edit Protocol dialog box

7. Modify the protocol parameters to meet your needs. Click the **Update** button on the **Edit Series** dialog box to save your changes.




Figure 175: Edit series update button

Note Different patient sizes and attenuations require specific protocols where technical parameters have been adjusted according to the physical characteristics of the patient. It is recommended that you select pediatric protocols when scanning children, rather than selecting adult protocols and modifying technique factors.

See “Creating a new protocol” on page 190 to learn how the fields and options perform to make informed choices on what to change.

Note Be sure to assign the **Build From** protocol a new **Protocol Description** before you make your additional changes.

8. When all required series have been modified click the **Update** button on the **Edit Protocol** dialog box.




Figure 176: Edit protocol update button

9. Click the **Close** button to exit.






Figure 177: Edit protocol close button

Copying and pasting protocols

To copy and paste protocols from one body part orb to another, including the Trauma orb.

1. Click **Tools >Protocol Manager** from the main menu.
The **Protocol Manager** dialog appears.

2. Click one of the following:

Adult	To create and/or scan with adult scan protocols, which are stored by anatomical location.
Pediatric	To create and/or scan with pediatric scan protocols, which are stored by anatomical location.
Trauma 	The Trauma orb can be used to store protocols commonly used for emergency scans.
	By selecting either an Adult or Pediatric patient, the corresponding list of saved protocols becomes available.

Note Adult and pediatric protocol parameters are customized to meet your requirements in conjunction with local and nationally recognized published guidelines. These protocols **must be** approved by your facility physicist **before** the system's acceptance.




Figure 178: Protocol Manager for Adult and Pediatric

3. Click the colored orb corresponding to the appropriate body part.




Figure 179: Anatomical orbs, in this case the chest orb

4. Review the protocol you would like to copy.
5. Highlight the protocol, right-click to see the floating menu, and click **Copy**.




Figure 180: Copy right-click floating menu

6. Go to body part orb you want to paste the protocol to, which can include the Trauma orb.

- Right-click to see the floating menu and click **Paste**.




Figure 181: Paste right click floating menu

- Click **Close** button to exit.

Deleting a protocol

- Click **Tools >Protocol Manager** from the main menu.
The **Protocol Manager** dialog box appears.
- Click one of the following:

Adult	To create and/or scan with adult scan protocols, which are stored by anatomical location.
Pediatric	To create and/or scan with pediatric scan protocols, which are stored by anatomical location.
Trauma	The Trauma orb can be used to store protocols commonly used for emergency scans.
	By selecting either an Adult or Pediatric patient, the corresponding list of saved protocols becomes available.

- Click the colored orb corresponding to the appropriate body part.
Select the protocol from list to be deleted.




Figure 182: Protocol Manager with a protocol selected

4. Click the **Delete** button.
The **Delete Confirmation** popup appears.




Figure 183: Delete Confirmation popup message – Yes or No to delete selection



5. Perform one of the following in the **Delete Confirmation** popup:
 - Click the **Yes** button to delete the selected protocol.
 - Click the **No** button to return to the Protocol Manager dialog box.The **Delete Confirmation** dialog box disappears, and the **Protocol Manager** dialog box appears.
6. Click the **Close** button to exit.

Adding breathing instructions to your protocol

Default audio files are installed on the workstation. Audio files can be attached to protocols and sent to the scanner. Each audio file has an indicator whether it was sent to the scanner.

1. Click **Tools >Protocol Manager** from the main menu.
The **Protocol Manager** dialog box appears.

2. Click one of the following:

Adult	To create and/or scan with adult scan protocols, which are stored by anatomical location.
Pediatric	To create and/or scan with pediatric scan protocols, which are stored by anatomical location.
Trauma 	The Trauma orb can be used to store protocols commonly used for emergency scans.
	By selecting either an Adult or Pediatric patient, the corresponding list of saved protocols becomes available.

3. Click the colored orb corresponding to the appropriate body part.
4. Click the protocol you would like to add **Breathing Instructions** to.
5. Click the **Edit** button.




Figure 184: Edit button

The **Edit Protocol** dialog box appears.




Figure 185: Edit Protocol dialog box

6. Select the **Protocol's Series** you want to add breathing instructions to.

The **Edit** button is enabled.




Figure 186: Add breathing edit button

7. Click the **Edit** button.
The **Edit Series** dialog box appears.

The 'Edit Series' dialog box contains the following parameters and controls:

Series Parameters		
Scan Type	Description	FOV Width
Helical	CHEST	
Scout Type	Start Position	FOV Top Left X
	0	
kV	End Position	FOV Top Left Y
120	300	
mA	Coverage	Bolus Tracking
175	300	Bolus Tracking
Slice Thickness/Spacing	Contrast	Auto-Start Auto-Stop Test Bolus
1.25 x 0.625		
Sharpness	Contrast Volume	Bolus Scan Time Threshold
Soft Tissue - Abdomen	0	
Resolution	Delay	AEC
1 Second(s)	10	Enable AEC
Pitch	Number of Images	Minimum mA Noise Level
0.8	480	
Body Part Examined	Scan Time	Maximum mA
CHEST	9.375	
Window Width	CTDIvol (mGy)	Recons
400	16.91	
Window Center	DLP (mGy.cm)	
40	507.3	

Buttons: Update, Reset, Cancel

Figure 187: Edit Series dialog box

- Click the **Use Breathe Indicator Audio** option and click the **Choose** button.

The 'Edit Series' dialog box is shown with the 'Use Breathe Indicator Audio' checkbox and the 'Choose...' button highlighted with an orange box.

Figure 188: Use Breathe Indicator Audio option

The **Breathe Indicator Audio Files** popup appears.




Figure 189: Breathe Indicator Audio Files popup

9. Select an audio file to use for the **Breathe** instruction.
The **Breathe** and **Hold** buttons are activated.
10. Click the **Breathe** button to place the file with **Breathe** files.
11. Select an audio file to use for the **Hold** instruction.
12. Click the **Hold** button to place the file with **Hold** files.
13. Click the **Apply** button to keep the files you selected to use with protocols.
14. Click the **Close** button to exit.

Importing protocols from a storage device

1. Click **Tools >Protocol Manager** from the main menu.
The **Protocol Manager** dialog box appears.
2. Click the **Import** button




Figure 190: Import button

3. The **Select File** popup appears.




Figure 191: Select File popup

4. Double click the Drive Letter that contains the protocols you want to import.




Figure 192: Select file

5. Click the file in the **Select File** popup.




Figure 193: Import button active in Select File when file(s) selected

6. Click the **Import** button.



Figure 194: Protocols Imported popup message – Protocols imported

7. The **Protocols Imported** popup appears.
8. Click the **OK** button.
9. Check that the required files have been imported.
10. Click the **Close** button to exit.

Exporting protocols to a storage device

1. Click **Tools >Protocol Manager** from the main menu.
The **Protocol Manager** dialog box appears.
2. Click the **Export** button.




Figure 195: Export button

3. The **Select Directory** dialog box appears.




Figure 196: Select Directory popup

4. If more than one **Media Device** is available, select the device to use.
5. Click the **Select** button.
6. The **Protocols Exported** popup appears.





Figure 197: Protocols Exported popup message – Protocols exported

7. Click the **Ok** button.
8. Check that the exported files are exported.
9. Click the **Close** button to exit.

Changing the order of protocols in the list

1. Click **Tools >Protocol Manager** from the main menu.
The **Protocol Manager** dialog box appears.

2. Click one of the following:

Adult	To create and/or scan with adult scan protocols, which are stored by anatomical location.
Pediatric	To create and/or scan with pediatric scan protocols, which are stored by anatomical location.
Trauma 	The Trauma orb can be used to store protocols commonly used for emergency scans.
	By selecting either an Adult or Pediatric patient, the corresponding list of saved protocols becomes available.

3. Click the colored orb corresponding to the appropriate body part.
4. Click the protocol to move up or down the list.
5. Click the **Up** arrow to move the protocol up the list; click the **Down** arrow to move the protocol down the list.




Figure 198: Changing protocol order with Up and Down (arrow) buttons

6. When you are finished ordering your protocols, click the **Save** button to save the new ordered list.




Figure 199: Protocol Save button

7. Click the **Close** button to exit.

Chapter 7 Daily Calibration and Quality Assurance

In this chapter, you will learn how to perform a daily air calibration and use the **Quality Assurance (QA)** tool that verifies the system is working as specified.

Keep in mind that **before** using the BodyTom Elite system, you **must** conduct a **Quality Assurance (QA)** test to verify the system is working as specified. Performing a daily (air) calibration

Note NeuroLogica recommends that an air calibration is performed every 6-8 hours. If the air freshness falls below 50%, or the scanner is moved to an area with a dramatic change in humidity or temperature, perform another air calibration to ensure optimum image quality during patient scanning.

If room-temperature fluctuations have occurred, you may need to perform more than one air calibration. In addition, scanners can drift out of alignment; make sure you perform a **QA Test** with the test phantom **before** scanning a patient.

Note It is recommended that the scanner is on for at least 60-90 minutes prior to performing the air calibration.

It is recommended practice that the scanner is plugged in and turned on even when it is not in use.

1. Make sure that nothing is in the bore before starting the daily air calibration.
2. Click **Tools > Perform Daily Cal** from the main menu.
The **Perform Daily Cal** popup appears.




Figure 200: Perform Daily Cal popup

Colors identify previous air calibrations outcomes:

Green	Indicates the calibration was successful
Yellow	Indicates the calibration is soon to expire
Orange	Indicates the calibration has expired
Red	Indicates the calibration failed.

3. Select one of the following options:
 - Click **Select All** to perform all calibration steps.
 - Click **Clear All** and individually select the calibration step(s) to perform.
4. Click the **Start** button.
The **Perform Daily Cal** popup appears, and the timer counts down.




Figure 201: Perform Daily Cal popup with count down

A warmup period begins, and the countdown begins; when completed the daily calibration will perform the calibration(s) you selected.

Note To stop the calibration, click the **Cancel** button to end the daily (air) calibration(s).

When the calibrations are completed the **Daily Calibration Summary** will display showing the status of the steps performed.



Figure 202: Perform Daily Cal summary popup

The following are the status indicators:

Green	Identifies the calibration completed successfully.
Yellow	Identifies the calibration is in progress.
Red	Identifies the air calibration failed.

- Click the **Close** button to exit the **Perform Daily Cal** popup.
The **Daily Cal** icon will change to green when it reaches a 100% air freshness.



Figure 203: Air freshness icon changes as the air quality drops from green to yellow to red

The QA phantom overview

The **QA phantom** is a device that measures parameters that characterize image quality; these parameters are as follows:

- Uniformity
- Noise
- High-contrast resolution
- Slice width
- Low-contrast resolution
- Sensitometry (contrast scale)




Figure 204: QA phantom

The **QA phantom** is a 20cm diameter disk consisting of a substrate made of **poly methyl methacrylate (PMMA)** containing specific inserts. The uniform area of the disk is used to measure uniformity and noise. Four other parameters are measured by the inserts in the substrate.

The QA phantom goes onto the phantom holder when performing a QA. See “Storing the QA phantom” on page 360.



Figure 205: Phantom holder

The **Axial** resolution wire, also called the **Modulation Transfer Function (MTF)** wire, is intended for measuring resolution in the **Axial** plane. Resolution is defined as the ability to distinguish small objects. It is expressed in line pairs per millimeter.

The **slice width wires** are the two inclined wires. They are intended to determine scanner resolution along the Z axis, that is, in the direction that is perpendicular to the **Axial** plane. Resolution along the Z axis is expressed in terms of slice width in millimeters. Although one wire is sufficient to measure the Z axis resolution if its position is accurately known, a second wire is included to confirm the alignment. If the alignment was incorrect, the results of the slice width test would not be accurate.

The low-contrast insert is a compound insert. It is made of two half cylinders of different materials with a known contrast difference between them. The low-contrast insert is intended to measure the contrast resolution of the scanner. The contrast resolution is the ability to measure slight differences in x-ray attenuation.

The sensitometry inserts are an air bore and cylinder made of different materials. They are intended to measure the contrast scaling of the scanner.

The QA scan protocols appear in the following table.

Table 28: Scan protocols used by the QA

Scan voltage	120 kV
Scan current	200 mA
mAs	400 mAs
Scan time	2 second
Kernel	Pos. Fossa/Vessel
Slice thickness	10mm

Starting Quality Assurance

To ensure the system is at its optimum, factory-specification level, the workstation provides **QA** tools to verify the system's state and to perform image-quality verification. To maintain consistent image quality over the system's lifetime, you should establish and maintain a regular **Quality Assurance (QA)** program. **QA** results are stored in the **Patient Browser**. Contact your local service representative to delete **QA** results.

The **QA protocol** is shipped with the system and appears when you click **Quality Assurance** from the main menu. You cannot customize or modify the **QA protocol**.

Note The QA test should be conducted per hospital requirements.

Before you begin this section, be sure to run a fresh **Daily Calibration** on the system. See "Performing a daily (air) calibration" on page 222.

Also, before you start the **QA protocol**, make sure the **QA phantom** is available and ready to install on the phantom holder. The phantom serial number label should be facing the front of the scanner and be positioned at the top – as shown in the figure below. The red insert should be on the bottom right when facing the scanner. The position of the phantom will greatly affect the **QA** results.




Figure 206: Phantom on the phantom holder

1. Move the QA stand to the front of the gantry, so the prongs point into the bore as shown.



Figure 207: Place QA phantom

2. Ensure the QA stand is centered in the bore using the sagittal laser as shown. If needed, adjust the prongs side to side.




Figure 208: Proper QA stand positioning

3. Place the QA phantom on the prongs as shown.




Figure 209: QA phantom positioning

1. The serial number label should be facing the front of the gantry and at the top.
2. The red insert should be on the bottom right.
3. The phantom should be in the middle of the carbon fiber post/prongs.
4. The two wires in the phantom

4. On the pendant, press the **Laser** button and align the internal laser to the etched line in the center of the phantom.



Figure 210: Laser button

See the laser precautions in “Laser safety” on page 58.

5. Click **Tools > Quality Assurance** from the workstation main menu.
6. The following **Quality Assurance** popup appears.




Figure 211: Quality Assurance popup

7. Click the **Prepare** button to begin the QA procedure.
8. The **System Ready to Scan** popup appears.




Figure 212: System Ready to Scan popup message – System is ready to begin

9. Go to the scanner and press the **START** button.
The system will scan the phantom and display the **QA Results** image.

QA Results				
Name	Value			
Radial Resolution At 10%	PASSED: 6.98,	HIGH LIMIT: 7.75,	LOW LIMIT: 6.25	
Radial Resolution At 50%	PASSED: 4.30,	HIGH LIMIT: 4.75,	LOW LIMIT: 3.25	
Tangential Resolution At 10%	PASSED: 6.98,	HIGH LIMIT: 7.75,	LOW LIMIT: 6.50	
Tangential Resolution At 50%	PASSED: 4.30,	HIGH LIMIT: 4.75,	LOW LIMIT: 3.25	
Slice Width	PASSED: 10.11,	HIGH LIMIT: 11.00,	LOW LIMIT: 9.00	
Noise	PASSED: 2.86,	HIGH LIMIT: 3.50,	LOW LIMIT: 2.50	
Low Contrast Resolution	PASSED: 6.00,	HIGH LIMIT: 6.00,	LOW LIMIT: 4.00	
Uniformity	PASSED: 4.44,	HIGH LIMIT: 5.00,	LOW LIMIT: 0.00	
CT of Air	PASSED: -991.73,	HIGH LIMIT: -950.00,	LOW LIMIT: -1030.00	
CT of Teflon	PASSED: 981.66,	HIGH LIMIT: 1004.00,	LOW LIMIT: 924.00	
CT of Acrylic	PASSED: 114.05,	HIGH LIMIT: 155.00,	LOW LIMIT: 75.00	

Figure 213: QA results of QA image

Note Items in green are passed results. Items in red are failed results. Often positional issues cause the failure; reposition your phantom and perform another scan. If you try multiple times and failures persist, call your service representative or **Technical Support**.

- Click the **Close** button on the **QA Results** popup when finished reviewing. The image of the phantom appears.



Figure 214: Phantom image

- Click the **Finalize** button on the workstation to exit the protocol.
- The **QA** appears in the **Patient Browser**; however, it is locked.

Locked	Study Date	Study Time	Acquisition	Referring Physician's Name	Study Description	Patient Name	ID	Patient Date of Birth	Patient
Locked	24/04/2017	8:39:50 AM			QUALITY ASSURANCE		21364401134376400001390949	20/01/1960	
	24/03/2017	8:06:05 PM			Age 15, No. Inside Head 7.3		21364401134376400001390949	20/01/1960	
	24/03/2017	3:55:57 PM			Age 15, No. Inside Head 7.3		21364401134376400001390949	20/01/1960	
	24/03/2017	3:25:58 PM			Age 15, No. Inside Head 6.5		21364401134376400001390949	20/01/1960	

Figure 215: Locked QA results shown in Patient Browser

Note See your service representative to remove locked QA results.

Ensuring good image quality

To produce consistent image quality over the system's lifetime, you should establish and maintain a regular **Quality Assurance (QA)** program. **QA** results are stored in the Patient Browser. Contact your local service representative to delete **QA** results.

- Compare the results to previous or optimum values and repeat these tests on a regular basis to detect changes in image quality values **before** any problem becomes visible.

Note If you notice degradation in image quality or a change in QA values, schedule a site visit and let your service representative or imaging physicist run more detailed tests.

Early intervention could prevent a major breakdown.

QA begins with baseline performance data that is acquired during system installation or after the repair or replacement of an x-ray generator-assembly, collimator, detector, Data Acquisition System (DAS) or main power circuitry.

2. Compare subsequent **QA** results against the baseline.

Baseline images can be saved for a visual comparison with QA checks, but measurement values provide a more objective way to monitor quality.



Name	Value
Radial Resolution At 10%	PASSED: 6.98, HIGH LIMIT: 7.75, LOW LIMIT: 6.25
Radial Resolution At 50%	PASSED: 4.30, HIGH LIMIT: 4.75, LOW LIMIT: 3.25
Tangential Resolution At 10%	PASSED: 6.98, HIGH LIMIT: 7.75, LOW LIMIT: 6.50
Tangential Resolution At 50%	PASSED: 4.30, HIGH LIMIT: 4.75, LOW LIMIT: 3.25
Slice Width	PASSED: 10.11, HIGH LIMIT: 11.00, LOW LIMIT: 9.00
Noise	PASSED: 2.86, HIGH LIMIT: 3.50, LOW LIMIT: 2.50
Low Contrast Resolution	PASSED: 6.00, HIGH LIMIT: 6.00, LOW LIMIT: 4.00
Uniformity	PASSED: 4.44, HIGH LIMIT: 5.00, LOW LIMIT: 0.00
CT of Air	PASSED: -991.73, HIGH LIMIT: -950.00, LOW LIMIT: -1030.00
CT of Teflon	PASSED: 981.66, HIGH LIMIT: 1004.00, LOW LIMIT: 924.00
CT of Acrylic	PASSED: 114.05, HIGH LIMIT: 155.00, LOW LIMIT: 75.00

Figure 216: Results of QA image after the QA test

Identifying filtration accuracy

Values of attenuation equivalent, half-value layer, and quality-equivalent filtration are expressed as a thickness of aluminum at the minimum of 99.9% purity. Attenuation of items in the x-ray beam should not be higher than 2mm of equivalent Aluminum (Al).



CAUTION Do not put anything in the x-ray beam that exceeds 2mm of equivalent AL as it may produce adverse effects to the image.

Using Axial plane to determine image resolution

The method to determine resolution in the **Axial** plane is to measure the modulation-transfer function of the scanner. A planar section of the **MTF** wire provides a point, called an **impulse**. The impulse is blurred by the imaging system, and the measurement of the blur quantifies the resolution. The blur is quantified by the **MTF** amplitude, which provides a measure of gain for a given object size in line pairs per centimeter (lp/cm).

The **MTF** is measured in two directions, called **radial** and **tangential** directions. The **radial** direction is along the line that joins the wire to the scanner isocenter. The **tangential** direction is perpendicular to the radial direction. The **MTF** along each direction produces a curve. The points at which each curve's amplitude is 50% and 10% of its amplitude at zero lp/cm are reported.

The expected results are given below.

Table 29: Modulation Transfer Function (MTF) direction

Direction	50%	10%
Radial (lp/cm)	4.7	7.2
Tangential (lp/cm)	4.7	7.2

Table 30: QA results

	Low limit	High limit
Slice width (mm)	9	11
Noise (HU)	2.5	3.5
Low-contrast resolution (mm)	4	6
Uniformity (HU)	0	5
CT of air (HU)	-1030	-950
CT of Teflon (HU)	924	1004
CT of acrylic (HU)	75	155

Using slice width

The method for determining the slice width for the **Axial** mode QA is to take an image of the inclined wire. The scanned section of the inclined wire is a line segment. The scanner blurs a scanned object in the **Axial** plane as well as in the direction perpendicular to it. The image of the inclined wire includes both the **Axial** plane blurring (**MTF**) of the scanner as well as the blurring in the z-direction. The slice width is determined by removing the component of in-plane blurring, by measuring the length of the wire segment and by using the known angle of wire inclination. The range is noted in QA results; see Table 29 on page 233.

Measuring noise

Noise is measured as the standard deviation of pixel values in a 1cm **Region of Interest (ROI)** at the center of the phantom. The **ROI** selection is automatic. The CTDI₁₀₀ center dose in a standard CTDI head phantom is 145 mGy for this scanning technique. The BodyTom Elite noise measurement is performed on a 10mm slice.

Measuring low contrast

Low-contrast resolution is measured as the difference between the mean CT values in each half of the low-contrast insert. An **ROI** is automatically selected around the low-contrast phantom and is automatically segmented into halves. Within each **ROI**, the mean pixel value is computed. The two mean values are subtracted.

The expected difference in the mean values is given in the electronic report. The low contrast should be: 5.0 ± 1.0 HU.

Finding uniformity

A **ROI** is automatically selected in each of five locations in the phantom. One **ROI** is at the center. Four outer **ROIs** are 60 to 70mm from the center of the phantom and spaced 90 degrees apart. A mean value is calculated in each **ROI**. The maximum difference between the means is calculated. The maximum allowable difference between the means is 3 HU.

Separate, independent measurements of the CT numbers of air, Teflon, and acrylic (see “Identifying CT contrast scale” on page 235 below) ensure that the overall CT number-scaling of the system is in order. If the uniformity portion of the QA protocol passes (maximum difference less than or equal to 3 HU), but the mean values of the five **ROIs** are all uniformly high or low in CT number, it is expected that the tests for air, Teflon, and/or acrylic will fail. See Table 39 on page 240.

Identifying CT contrast scale

Contrast scale represents the attenuation scaling of the scanner. The mean CT numbers of each of the sensitometry objects is calculated and reported.

Using performance phantoms

The phantoms (in the following sections) were used for measurement of dose and imaging performance.

Measuring dose

The dose is measured using the standard CTDI body phantom.

Identifying load factors

Table 31: Load factors

Protocol description	kV	mA	Time (seconds)
Axial	80-140	30-300	1
Helical	80-140	30-300	1 per rotation



CAUTION When conducting multiple or repeat scans, ensure that the total exposure does not exceed 1Gy CTDI.

Note The highest x-ray tube current is 300mA and the highest x-ray tube voltage selection at this current is 140kV.

The nominal x-ray output power is 42kW when operating at an x-ray tube voltage of 140kV and x-ray current of 300mA for 4 seconds.

The x-ray tube voltage/current tolerance is $\pm 10\%$.

The nominal x-ray tube voltage is 140kV with 100cfm minimum cooling flow.

The BodyTom Elite dose information (21 CFR 1020.33 c)

Dose is measured using standard CTDI head and body phantoms. Surface and center CTDIs were both measured. Weighted CTDI is computed using surface and center CTDIs:

$$CTDI_w = \left(\frac{2}{3} CTDI_{surf} + \frac{1}{3} CTDI_{cen} \right)$$

Measured values are normalized to scan current, for example, CTDI values are in mGy/100 mAs. For any given scan protocol $CTDI_w$ can be estimated using following equation and data from Table 32 and Table 33 on page 238:

$$CTDI_w(kV, m, S) = \left(\frac{m}{100.0} \cdot S \right) CTDI_w(kV, 100mAs) mGy$$

$CTDI_w$ can also be computed using data from Table 34 and Table 35 on page 238, and the following equation:

$$CTDI_w(kV, m, S) = \left(W(kV) \cdot \frac{m}{100.0} \cdot S \right) CTDI_w(120_{kV}) mGy$$

Where **W** is the **kV** relative dose ratio with respect to 120 kV. **m** is the x-ray tube current in mA and **S** is the scanning time in seconds. If scan kV matches measured scan voltage, then **W** is equal to **1.0**. For **Helical** scans, $CTDI_{vol}$ is calculated as follows:

$$CTDI_{vol} = \frac{CTDI_w}{Pitch}$$

For **Axial** scans, $CTDI_{vol}$ is calculated as follows:

$$CTDI_{vol} = \frac{CTDI_w}{Scan Increment}$$

For example, $CTDI_w$ of 2 seconds, 100 kV **Axial** head scan, with 200 mA can be evaluated as follows using data from Table 32 and Table 33 on page 238:

$$\begin{aligned} CTDI_w(100kV, 200mA, 2 Sec) &= \left(\frac{200.0}{100.0} \times 2 \right) CTDI_w(100_{kV}, 100mAs) \\ &= 4 \times (100_{kV}, 100mAs) mGy = \\ &4 \times 14.4 = 57.6 mGy \end{aligned}$$

Using data from Table 34 and Table 35 on page 238, $CTDI_w$ can be computed as follows:

$$\begin{aligned} CTDI_w(100kV, 200mA, 2 Sec) &= \left(0.626 \times \frac{200.0}{100.0} \times 2 \right) CTDI_w(120_{kV}, 100mAs) \\ &= 2.504 \times (120_{kV}, 100mAs) mGy = \\ &2.504 \times 23.0 = 57.59 mGy \end{aligned}$$

$CTDI_{vol}$ is evaluated as follows:

$$CTDI_{vol} = \frac{CTDI_w}{Scan Increment} = \frac{57.6}{1.0} = 57.6 mGy$$

In another example, we can compute $CTDI_w$ of 1 second, 120 kV **Helical** abdomen scan, with 200 mA and a pitch of 0.8. Dose can be evaluated as follows, using data from Table 32 and Table 33 on page 238:

$$\begin{aligned} CTDI_w(120kV, 200mA, 1 Sec) &= \left(\frac{200.0}{100.0} \times 1 \right) CTDI_w(120_{kV}, 100mAs) \\ &= 2 \times (120_{kV}, 100mAs) mGy \\ &= 2 \times 7.56 = 15.12 mGy \end{aligned}$$

Using data from Table 34 and Table 35 on page 238, $CTDI_w$ can be computed as follows:

$$\begin{aligned} CTDI_w(120kV, 200mA, 2 Sec) &= \left(1.0 \times \frac{200.0}{100.0} \times 1 \right) CTDI_w(120_{kV}, 100mAs) \\ &= 2.0 \times (120_{kV}, 100mAs) mGy = \\ &2.0 \times 7.56 = 15.12 mGy \end{aligned}$$

$CTDI_{vol}$ is evaluated as follows:

$$CTDI_{vol} = \frac{CTDI_w}{Scan Pitch} = \frac{15.12}{0.8} = 18.9 mGy$$

Body CTDI_w phantom

CTDI_w using CTDI body phantom is listed in the following table. Data was measured using the 32 rows collimation with the phantom placed on the phantom holder. Dose measurements were taken using raw data acquisition in **Service** mode.

Note Performing scans in different **Acquisition** modes can cause slight variations in measured dose.

Table 32: Body CTDI_w (mGy/100mAs)

	140 kV	120 kV	100 kV	80 kV
CTDI ₁₀₀ Center (C)	6.81	4.54	2.65	1.20
CTDI ₁₀₀ Surface (S)	14.1	9.84	6.49	3.48
CTDI _w	11.7	8.07	5.21	2.72

Head CTDI_w phantom

Weighted average Computed Tomography Dose Index (CTDI_w) using the CTDI head phantom is listed in the following table. Data was measured using the 8 rows collimation. Dose measurements were taken using raw data acquisition in **Service** mode using phantom holder.

Table 33: Head CTDI_w (mGy/100mAs)

	140 kV	120 kV	100 kV	80 kV
CTDI ₁₀₀ Center (C)	34.3	23.6	14.6	7.55
CTDI ₁₀₀ Surface (S)	38.1	26.7	17.0	9.33
CTDI _w	36.8	25.7	16.2	8.74

Normalized CTDI tables are listed below. CTDI is normalized with respect to a typical 120kV scan protocol:

Table 34: Normalized CTDI of body phantom

	140 kV	120 kV	100 kV	80 kV
CTDI ₁₀₀ Center (C)	1.50	1.00	.584	.264
CTDI ₁₀₀ Surface (S)	1.43	1.00	.660	.354
CTDI _w	1.45	1.00	.646	.337

Table 35: Normalized head CTDI

	140 kV	120 kV	100 kV	80 kV
CTDI ₁₀₀ Center (C)	1.45	1.00	.619	.320
CTDI ₁₀₀ Surface (S)	1.43	1.00	.637	.349
CTDI _w	1.43	1.00	.630	.340

The BodyTom Elite dose in air

Dose measurements were taken using raw data acquisition in **Service** mode.

Table 36: CTDI air (mGy/100mAs)

	140 kV	120 kV	100 kV	80 kV
32 rows	31.2	22.5	14.9	8.56
8 rows	47.3	33.4	21.9	12.6

Dose is measured using a typical head protocol and a typical abdomen protocol. Dose in air was also measured for repeatability over 10 scans. Average value and standard deviation are noted below:

Table 37: Mean and standard deviation of CTDI air

	8 rows	32 rows
Mean mGy	19.3596	13.18493
Standard deviation mGy	0.058716	0.017831

QA measurements

The QA phantom is typically used to monitor the scanner on site; however, the following phantoms can be used for measuring the imaging performance of the scanner.

Note Actual results on installed units can vary 20% due to machine and test tolerances.

Tube accuracy

Table 38: Tube accuracy

kV	Tolerance
140	± 5%
120	± 5%
100	± 5%
80	± 5%

ACR testing procedure

Most sites use the ACR phantom for evaluating the QA parameters of the scanner. Furthermore, each scanner is evaluated using the ACR phantom prior to shipping. Due to the special tube filtration some of the limits for CT values may be different from those set by the ACR committee. Table 39 lists the NL limits for the CT number and linearity of each insert in the ACR phantom. It also lists the limits set by ACR. The difference in CT

number is mainly due to the difference of the x-ray beam quality due to the tube filtration and the ACR committee acknowledge this in published papers³.

Table 39: The CT number and linearity of the different inserts in the ACR phantoms

Insert Material	NL Limits	ACR Limits
Air	-1005 to -970	-1005 to 970
Polyethylene	-110 to -85	-110 to -85
Water	-7 to 7	-7 to 7
Acrylic	110 to 135	110 to 130
Bone	1010 to 1110	850 to 970

The scan protocols are typically selected by the site physicist in the CT manager. However, ACR recommends the use of standard head and abdomen protocols. NeuroLogica uses the protocols listed in Table 40.

Table 40: The NeuroLogica head and abdomen ACR scan protocols

Protocol	Head/Abdomen	Abdomen
kVp	120	120
mA	200	250
Time per rotation (seconds)	1	1
Dose (Weighted)	45 to 50 mGy	22 to 25 mGy
Scan FOV	59.5 cm	59.5 cm
Display FOV (minimum)	25.0 cm	25.0 cm
Reconstruction sharpness	Soft tissue	Soft Abdomen
Scan type	Axial	Helical
Z-axis collimation	10 mm	40 mm
Table increment (mm) or Table speed (mm/rot)	10 mm	32 mm/rot
Slice thickness	5 mm	5 mm
Scan time (seconds)	1	1
Slice separation	5 mm	5x5 mm
Number of images per scan	2	N/A

³ The ACR committee will often change their limits, as such, the limits listed in Table 39 may have been changed.

Identifying high-contrast resolution

The high-contrast resolution phantom is a wire placed at the center of a uniform disk. The wire provides an impulse function in the **Axial** plane when it is placed parallel to the scanner-gantry axis-of-rotation. The high-contrast resolution is measured from the **Modulation Transfer Function (MTF)**. Typical **MTF** curves are shown in the following figures. Variations of 10% may occur in measurements due to phantom placement error and measurement inaccuracies.




Figure 217: MTF

Table 41 lists the 50%, 20% and 10% cutoffs of the most commonly used kernels on the scanner. The cutoffs were measured using the MTF curve for each kernel like the one displayed in Figure 217.

Table 41: The cutoffs of some of the common reconstruction kernels

Scan Type	Kernel	MTF50%	MTF20%	MTF10%
Axial	Soft Tissue	3.4038	4.9538	5.7187
Axial	Pos. Fossa/Vessel	3.9930	5.7816	6.7455
Axial	Sharp	6.6819	8.2549	9.2587
Axial	Bone	7.7757	9.6942	10.9681
Axial	Sharp Lung	6.1620	8.9318	11.8802
Axial	High-Res QA	7.9286	10.5069	12.4050
Helical	Bone Head	5.2370	6.5734	7.1577
Helical	Soft Tissue-Abd	6.2422	7.2252	7.7846
Helical	Soft Tissue-Head	3.0327	4.4322	5.1328
Helical	Bone-Abdomen	5.5166	6.8081	7.4292

Noise, uniformity, and mean CT number of water

One of two phantoms may be used in these tests. These are Catphan® 412 or a cylindrical 20cm diameter water cylinder.

The variation in standard deviation may be $\pm 10\%$ due to variations between systems.

Noise is measured as the standard deviation at isocenter. The value is 2.1 ± 2 HU when the imaging protocol is 140 kV, 42 mAs and standard kernel. This protocol gives a CTDI₁₀₀ center dose of 160 mGy.

Noise is measured as the standard deviation of pixel values in a 1cm **ROI** at phantom's center. Range is as noted in QA results (see Table 30 on page 233).

Table 42: Uniformity and Mean CT Number using Water Phantom

Description	Noise (HU)
Body (cp300 mm Water phantom)	11.9HU
Head (cp200mm Water Phantom)	2.89HU

Uniformity and mean CT number

The mean CT number of air is -1000 HU and that of water is 0 HU. The tolerance of the mean CT number will be ± 3 HU. For mean CT numbers measured at different points of the water phantom, the maximum difference in the means will be less than 4 HU⁴.

An **ROI** is automatically selected in each of five locations in the phantom. One **ROI** is at center. Four outer **ROIs** are 60 to 70mm from phantom center and spaced 90 degrees apart. A mean value is calculated in each **ROI**. The maximum difference between means is calculated. The maximum allowable difference is as noted in QA results (see Table 30 on page 233).

Maximum difference between periphery **ROIs** and center **ROI** mean CT values in an image is less or equal to 4 HU. The maximum error in CT number of water is ± 3 HU.

Table 43: Uniformity and Mean CT Numbers using Water Phantom

Description	Uniformity (HU)	Mean CT Number (HU)
Body (cp300 mm Water phantom)	1.91	0.32
Head (cp200mm Water Phantom)	1.21	1.21

⁴ For States that required the sites to perform water phantom testing, please contact customer service for assistance.

Low-contrast resolution

The phantom used for low-contrast resolution measurement is CTP 515 section of the Catphan 600.

The **low-contrast resolution** is 4mm rod at 0.3% contrast when the center CTDI_w dose is 71 mGy. The imaging protocol is 120 kV, 300 mAs with 1 rotation, 5mm slice thickness, and using the soft tissue filter.




Figure 218: Catphan 515 using 120kV, 300mA, 1 rotation, and 5mm slice

Geometric efficiency in the Z axis direction notification

This is the ratio of the integral of the dose profile integrated over the detector width in Z divided by the total CTDI₁₀₀. The geometric efficiencies for the two available collimations are listed in Table 47 on page 251. The scanner has two scan modes: Axial and Helical. The Helical scan mode uses the wide collimation of 32 rows therefore its accuracy is around $83 \pm 5\%$. The Axial scan uses the narrow collimation of 8 rows only with its geometric efficiency of $60 \pm 5\%$. When this occurs, the operator will be given an opportunity to continue/affirm the prescribed scan. See Figure 219.




Figure 219: Geometric efficiency confirmation

Half-value layer

Table 44: Half-value layer

Scan voltage (kV)	100	120	140
Half value	6mm	7mm	8mm

Allowable variations

The following are allowable variations:

Dose	A ±5% variation in dose may occur due to variations between systems and measurement differences. The maximum variation is ±10%.
High Contrast Resolution	The variation in values on the MTF curve may be ±10%. These will occur mainly due to phantom placement errors, measurement inaccuracies and system variations.
Noise	The variation in standard deviation may be ±10% due to variations between systems.
Uniformity	The maximum difference between ROI means in an image is 4 HU. The maximum error in the CT number of water is ±3 HU.

Dose: Maximum variation is ±10%.

Variation in values on the **MTF** curve may be ±10%. These occur mainly due to phantom placement errors, measurement inaccuracies and system variations.

Scatter radiation

Reference the following radiation scatter plots identifying proper distances to protect from radiation exposure. The scatter plots provide scattered radiation dose in air-kerma, per current-time product in both standard and SI units for nominal technique of 120kV (μRad/100mAs and μGy/100mAs respectively). This information is given so the facility physicist and/or **Radiation Safety Officer (RSO)** can use these charts to calculate exposure with the following formula:

$$\text{Stray radiation (scan current, scan voltage)} = \text{stray radiation (100, 120)} \times \left(\frac{\text{scan current}}{100} \right) \times \left(\frac{\text{scan voltage}}{120} \right)^{2.3}$$

In addition, per IEC 60601-2-44, 3rd Edition, “Particular requirements for the basic safety and essential performance of x-ray equipment for computed tomography”, scatter plots are provided as shown in Figure 222 on page 247 and Figure 223 on page 248 for the maximum techniques settings of 140kV and 300mA (standard and SI units respectively).

This information is specifically intended for the facility Physicist and/or an **RSO** to perform a safety and shielding analysis such as described in NCRP 147, “Structural Shielding Design for Medical X-Ray Imaging Facilities.”



WARNING Exposure to secondary radiation can be harmful, and scanner usage should only be done under the direct supervision of the facility’s qualified **Radiation Safety Officer (RSO)** in compliance with site, local, state, provincial, and national regulations. Only this **RSO** can perform the calculations necessary to determine what additional safety precautions are necessary, such as shielding, personal protections, and so on.

Note The BodyTom Elite scanners are compatible with IRR1999 and EU Directive 96/29/EURATOM.

Typical application environment and radiation safety

The BodyTom Elite is an advanced radiation protection mobile CT. There is an effective x-ray shielding that is equivalent of 0.75mm of lead within the gantry. Under normal circumstances, no additional radiation shielding is needed per recommendations of the facility physicist. The scanner can be used in a mobile environment and/or within an enclosed environment.

The scatter plot below shows the dose map during a normal scan:







Figure 220: Scatter plot (120kV, 100mA in μRad)

Note In compliance with IEC 60601-2-44:2009, section 203.11, the above figure shows the scatter radiation measured at the edge of the gantry in the tomographic plane is $\approx 20\%$ of the scatter radiation measure at the same distance along the axis of rotation in the horizontal plane.

The black box (located at 0 on the X axis and -2 on the Y axis) represents an approximate (24 x 24 x 79in. or 60 x 60 x 200cm) zone of occupancy. The system in **Scan** mode stands at 78.5in., which meets the 79in. (200cm) requirement.

Figure 221: Scatter plot (120kV, 100mA in μGy)Figure 222: Scatter plot (140kV, 300mA in μRad)

Figure 223: Scatter plot (140kV, 300mA in μGy)

Note The anatomical body phantom was placed on a scan table inside the gantry to consider scatter through patient. Measurements were made using the following scan protocol: 140kV, 200 mA, and 5 sec. The following figures show measurement points in vertical X—Y and perpendicular Y—Z planes, followed by corresponding tables detailing resulting data.




Figure 224: Scatter measurements (X—Y plane)

Table 45: Scatter measurements (X—Y plane) ($\mu\text{Rad}/100\text{ mAs}$)

	A'	B'	C'	D'	E'	F'
A	12.2	15.3	17.7	18.3	15.7	12.7
B	11.9	14.8	16.4	16.7	16.0	12.6
C	11.4	13.6	11.1	12.9	15.2	12.5
D	9.40	10.4	8.87	9.57	12.4	10.4
E	6.09	6.26	4.09	4.26	7.66	7.66




Figure 225: Scatter measurements (Y—Z plane)

Table 46: Scatter measurements (Y—Z plane) ($\mu\text{Rad}/100\text{ mAs}$)

	A'	B'	C'	D'	E'	F'
A	300	676	852	591	870	437
B	218	539	2320	2940	1090	465
C	137	328	931	1120	461	233
D	27.8	1080	2790	2240	765	282
E	8.96	844	441	85.3	538	345

Dose profile/Geometric Efficiency

A graphical presentation of the **dose profile** along a line – Z perpendicular to the **tomographic plane** and centered at the **isocenter**, determined in free air for one **Axial** scan, in the center location of the head-dosimetry phantom, and the center location of the body-dosimetry phantom – is given in the accompanying documents for each selectable value of $N \times T$. When more than three different values of $N \times T$ are available, the information is provided for at least the minimum, maximum and one mid-range value. The **dose profile** is presented on the same graph and to the same scale as the corresponding **sensitivity profile** required by 203.111.

Dose profile was measured for 32 rows by taking a stationary scan with the radio-chromic film, centered on top of the detector array. The scan protocol was as follows: 120 kV, 200 mA, 5 sec. After the scan was taken, the radio-chromic film was scanned, and the profile extracted. The following figure shows the dose profile and detector array.

Geometric efficiency was calculated as the ratio of the detector array to FWHM of dose profile using the following formula:

$$Geom_{Efficiency} = \frac{N_{rows} \cdot w_z \cdot M_f}{FWHM}$$

where N_{rows} is the number of detector rows; w_z is the width of the detector in the z-direction; M_f is the magnification factor; $FWHM$ is the full width at half maximum of the profile. Table 47 on page 251 lists the measured geometric efficiencies for the two existing collimations of the BodyTom Elite.




Figure 226: Dose profile for 8 rows




Figure 227: Dose profile for 32 rows

Table 47: The geometric efficiency of the two different collimations of the BodyTom Elite

Collimation	Geometric Efficiency
32 rows collimation	$83 \pm 5\%$
8 rows collimation	$60 \pm 5\%$

Chapter 8 Patient Registration

Patient Registration is the first step in the patient scan process. You can register a patient in the following ways:

- Manually register a patient from the **Patient Registration** tab.
- Perform a query to acquire already-entered patient data from the **Hospital Information System (HIS)** or **Radiology Information System (RIS)**.

It is assumed that the workstation is connected to the site's **HIS/RIS** system. If you are not connected, you can always manually register a patient.






Figure 228: Activated Patient Registration tab






Navigating the Patient Registration screen

Make sure the **Patient Registration** tab is selected.

Notice the buttons at the bottom of the **Patient Registration** dialog box. Many of these buttons are active *only* if you are already connected to the site's **HIS/RIS** or if you clicked the **Query** button, and the list of patients populated in the **Query Results** list. When a patient is selected, the buttons are active.

Table 48: Patient Registration buttons

Patient Registration button	Action
	Searches the HIS/RIS server for scheduled patients. The population of patient information could take several minutes to appear, depending on the number of patient entries the query retrieves after clicking the Query button.
	Cancels the current query. Entries retrieved prior to cancellation appear in the Query Results list and Stored Results if they are moved there.
	Registers the selected patient and takes you to the Acquisition tab to select a protocol to be used for scanning.

Patient Registration button	Action
	Shows selected patient details.
	Searches queried patient entries for specific information.
	Selects patient(s) from query results and moves them into the Stored Results list.
	Removes patient(s) from the Stored Results list.
	Manually enters a new patient and, when completed, takes you to the Acquisition tab to select a protocol to be used for scanning.

Registering the patient

The following procedures show you how to register or enter patient information into the system before scanning a patient. Patients are registered manually or queried from the **Hospital Information System/Radiology Information System (HIS/RIS)**. The system can be configured to add or create specific patient information when the patient is registered.

Querying patient information

1. If necessary, click the **Patient Registration** tab on the main screen.



Figure 229: Patient Registration tab

- Click the **Query** button at the bottom of the screen.
The **Query Information** dialog box appears.

The dialog box titled "Query Information" contains a dropdown menu for "HIS/RIS Server" set to "NeuroLogica HIS/RIS". Below this is a table with two columns: "Name" and "Value". The table lists several query fields with their corresponding values. At the bottom of the dialog are three buttons: "Query", "Reset", and "Cancel".

Name	Value
(0040,0002) Scheduled Start Date	20170427
(0008,0060) Modality	CT
(0010,0020) ID	
(0010,0010) Patient Name	
(0008,0050) Accession	
(0040,0001) Scheduled Station AE Title	
(0040,1001) Requested Procedure ID	

Figure 230: Query Information dialog box

- Click the **HIS/RIS Server** dropdown and select the worksite to pull data from.
The default worksite appears at the top. If there is no list, see your site administrator to set it up.
- Double-click any of the named **Query Fields** you would like to use to query for patients by entering the value in the **Value** column.
A popup associated with the **Query Field** you selected appears. For example, if you double-click the **Scheduled Start Date** row, the **Edit Value** popup appears. Enter the desired start date. Another example would be to click the **Patient Name** value row. Again, the **Edit Value** popup appears; however, this time **Patient Name** text boxes are provided so you can type the patient's name to query. You can click any of the **Value** rows to fill in data to help query for the patient you are looking for. You can enter as much or as little information as needed. If no information is available, leave the value blank.




Figure 231 shows a 'Edit Value' dialog box. It contains the text 'Please enter new value for the following:' followed by '(0010,0010) Patient Name'. Below this are five text input fields labeled 'Prefix', 'First Name', 'Middle Name', 'Last Name', and 'Suffix'. At the bottom of the dialog are two buttons: 'Update' and 'Close'.

Figure 231: Edit Value popup for name

5. When you are finished filling in query selections, perform one of the following:
 - Click the **Update** button to query based on the newly entered data to help narrow down your search.
 - Click the **Close** button to remove any changes and return to the previous **Query Information** popup.
6. A list of patients matching your selected criteria variables populates in the **Query Results** list on the **Patient Registration** tab.




Figure 232 shows a screenshot of the 'Patient Registration' tab in the software. At the top, there is a 'Query Results' table with several columns. Below this table is a large, empty area for displaying patient data. At the bottom of the window, there is a row of buttons including 'Query', 'Filter', 'Sort', 'Print', 'Export', 'Import', 'Refresh', and 'Close'.

Figure 232: Patient Registration Query Results table

7. Select a patient and click the **Register** button to register the patient for the exam.
The system enables and opens the **Acquisition** tab. To perform the acquisition steps, see "Performing a scan" on page 264.

Storing patients in the Stored Results list

This list is helpful when multiple patients need to be scanned and a connection to a worksite like **HIS/RIS** is unavailable at the exam location.

1. If necessary, go to the **Patient Registration** tab to query the patients(s).
2. Perform steps 2 through 5 in “Querying patient information” on page 253.
3. Click the **Query** button.
Let the criteria you selected populate into the **Query Results** list area.
4. Select one or more patient entries from the **Query Results** list.
Select patients in the following ways:
 - To select one patient, click anywhere in the patient’s row.
 - To select more than one patient at a time, press and hold the **Ctrl** key and click patient entries until finished and release the **Ctrl** key.
 - To select all the patients, press and hold the **Shift** key, click the first patient in the list, then click the last patient to highlight all patients between the first patient selected and the last.
5. Click the **Store** button.
The patient information you selected appear in the **Stored Results** list at the bottom of **Patient Registration**.




Figure 233: Patient Registration Stored Results table

6. Click the patient you want to select from the **Stored Results** table.
7. Click the **Register** button to register the patient for the exam.
The system enables and opens the **Acquisition** tab. To perform the acquisition steps, see “Performing a scan” on page 264.

Manually registering a patient


You manually register a patient for examination when the **HIS/RIS** server is unavailable, the patient cannot be found, and/or was never entered into the system.

1. If necessary, go to the **Patient Registration** tab.



Figure 234: Patient Registration tab

2. Click the **Manual** button at the bottom of **Patient Registration**. The **Exam Information** dialog box appears with the **Patient** tab open.



Notice that the value appears with the last name first, separated by ^, the first name next, separated by ^, the middle name (if you supplied that information), and any other information you entered.

Note The **Patient ID** is automatically generated by the system; you can replace this automated identifier with your patient's hospital ID number.

Name	Value
00000000 Patient Name	
00000000 Patient ID	000000000000000000000000
00000000 Patient Date of Birth	
00000000 Patient's Sex	
00000000 Patient's Age	

Register Cancel

Figure 236: Patient ID field


5. For the **Patient Date of Birth**, perform the following:
 - Double-click the **Patient Date of Birth** field.
 - Enter the patient's birth date in the **Patient Date of Birth** field. Be sure to move the cursor to the far left to ensure two digits are included for the month and the day; four digits are required for the year.
 - Perform one of the following:
 - Click the **Update** button to save your work and close the **Edit Value** dialog box.
 - Click the **Close** button to close the **Edit Value** dialog box without saving your work.
6. For the **Patient's Sex**, perform the following:
 - Double-click the **Patient's Sex** field.
 - Enter the patient sex in the field by entering the appropriate letter:
 - **F** for Female
 - **M** for Male
 - **O** for Other
 - Perform one of the following:
 - Click the **Update** button to save your work and close the **Edit Value** dialog box.

- Click the **Close** button to close the **Edit Value** dialog box without saving your work.

7. Perform one of the following:

- Click the **Register** button to register your patient data.
- Click the **Cancel** button to exit without entering your data.

When you click the **Register** button, the system enables and opens the **Acquisition** tab.



	Name	Value
(0010.0010) Patient Name	Stevens, John	
(0010.0020) ID	216840.114379.4000.1.20170901.133202.4700	
(0010.0030) Patient Date of Birth	19971126	
(0010.0040) Patient's Sex	M	
(0013.0020) Patient ID		

Figure 237: Patient data filled in

After your patient is registered, you can view the **Patient Exam Details** to ensure your data is correct.

If it is not correct, go to the next step to make the necessary changes.

8. Click the **Expand** link.

The **Exam Information** popup appears.




Figure 238: Expand link in context and close up

9. Make your changes in the **Exam Information** popup.




Figure 239: Exam Information popup

10. Click the **Update** button to save changes.

Viewing patient information

This procedure lets you view, but not change, the patient information.

1. If necessary, click the **Patient Registration** tab on the main screen.



Figure 240: Patient Registration tab

2. Select a patient from the **Query Results** list or the **Stored Results** list.
3. Click the **View** button.
4. Review the patient's information.
This popup presents static information that you cannot change.
5. Click the **Close** button to exit the **View Entry Information** popup.

Deleting patients from the Stored Result list

Patient information can be manually deleted from the **Stored Results** list, you cannot delete patients from the **Query Results** list.

1. If necessary, go to the **Patient Registration** tab.



Figure 241: Patient Registration tab

2. Select one or more patients from the **Stored Results** list to delete.
Select patients in the following ways:
 - To select one patient, click anywhere in the patient's row.
 - To select more than one patient at a time, press and hold the **Ctrl** key and click patient entries until finished and release the **Ctrl** key.
 - To select all the patients, press and hold the **Shift** key, click the first patient in the list, then click the last patient to highlight all patients between the first patient selected and the last.
3. Click the **Delete** button.
The patients you selected are removed from **Stored Results** list.

Chapter 9 Patient Scanning

After you register the patient, the **Acquisition** tab automatically opens. The **Acquisition** tab lets you check that the selected patient information is accurate before you perform the scan. The **Acquisition** tab is also where you select protocols for the scan before you scan the patient. A protocol determines the parameters used to acquire patient images.



Figure 242: Active Acquisition tab









After the protocol is selected, you can scan the patient. See “Performing a scan” on page 264.



The following table provides information on the buttons on the **Acquisition** tab and what they are used for. Later you will learn how to set protocols for the scan.



CAUTION When conducting multiple or repeat scans, make sure the total exposure does not exceed maximum limit of 1Gy.

Table 49: Acquisition buttons

Acquisition button	Action
	Allows you to modify the protocol selected or choose a new protocol.
	Authorizes the scanner to move to the next step if applicable.
	Begins any post-reconstructions that were defined during the protocol setup.
	Allows you to repeat a portion or all the scan.
	Allows you to extend the currently active scan.
	Allows you to pause the scan acquisition.
	Allows you to pause the entire multi-step protocol acquisition.
	Cancels the current scan within a protocol.

Acquisition button	Action
	Cancels the entire exam.
	Completes the examination.

The following shows what appears in **Acquisition**:



Figure 243: What appears on Acquisition

Identifying Scan Types

Scan types identify how images are acquired during a scan. The following Scan types are available.

Axial

The **Axial** scan type lets you scan in the **Transverse** plane. Data is acquired as the x-ray tube rotates around the patient.

Helical

The **Helical** scan type acquires data continuously as the x-ray tube rotates around the patient and the scanner translates over the patient in the Z axis.

Dynamic

The **Dynamic** scan type acquires data at multiple time points over the same anatomic location while the scanner remains stationary; x-ray exposure can be continuous or intermittent.

Reference

The **Reference** scan type acquires a single 10mm slice to review anatomical position or place the **Region of Interest (ROI)** for **Bolus Tracking** scans. **Reference** scanning can only be used in conjunction with **Helical** and **Dynamic** scanning during a **CT Angiography (CTA)** or **Perfusion** protocol.

Scout

The **Scout** scan type acquires data continuously as the x-ray tube remains stationary at a designated angle and the scanner translates over the patient in the Z axis. The resulting **2D** projection is used during scan planning.

Performing a scan

You cannot complete this procedure without a registered patient.

Note If the scan needs to be stopped, perform the following:

For an immediate or hard stop, press the **E-STOP** button. This stops x-ray, centipede movement, and gantry rotation immediately.

For a controlled stop, press the **Cancel Scan** button.

Note Be sure the scanner is calibrated for the room you will scan in. See "Selecting a room for the BodyTom Elite" on page 188 and "Performing a daily (air) calibration" on page 222.



1. From the workstation, go to the **Patient Registration** tab to assign the patient to the scan in one of the following ways:
 - Query an already existing patient from the HIS/RIS.
 - Manually register the patient.
See "Chapter 8 Patient Registration" on page 252.The **Acquisition** tab will be activated when the patient is registered.
2. From the **Acquisition** tab, click the **Protocol** button to open the **Exam Planner** dialog box.
3. Move scanner and align patient as needed.

See “Positioning the scanner before a scan” on page 96.

4. On the pendant, press the **Laser** button to turn on the laser and use it to align the patient to the scanner.

See “Positioning the patient using the laser lights” on page 97.

5. On the workstation screen, click the appropriate option:

Adult	To create and/or scan with adult scan protocols, which are stored by anatomical location.
Pediatric	To create and/or scan with pediatric scan protocols, which are stored by anatomical location.
Trauma 	The Trauma orb can be used to store protocols commonly used for emergency scans.
	By selecting either an Adult or Pediatric patient, the corresponding list of saved protocols becomes available.

Note Adult and pediatric protocol parameters are customized to meet your requirements in conjunction with local and nationally recognized published guidelines. These protocols **must be** approved by your facility physicist **before** the system’s acceptance.




Figure 244: Exam Planner for Adult and Pediatric

6. Click the colored orb corresponding to the appropriate body part you will scan.




Figure 245: Anatomical orbs, with the Chest orb selected

7. Click the appropriate protocol from the list.
8. Click the **Edit** button to review the selected protocol.




Figure 246: Protocol selected and Edit button active

The **Edit Protocol** dialog box appears.




Figure 247: Edit Protocol dialog box

The **Protocol Information** tab displayed on the left and the Protocol's **Series** boxes displayed on the right show the series that are already created. The **Patient Position** appears identical whether it is for an adult, pediatric, or trauma patient.

Note You can modify a protocol; however, changes you make from **Acquisition** will not be saved permanently. Permanent changes to protocols can only be made by in **Protocol Manager**.

Assuming you have the proper user privileges, you can modify protocol parameters such as, kV, mA, and coverage at the time of the scan, but the modifications will not be saved for future use.

9. To edit an existing protocol, perform the following:
 - In the **Edit Protocol** dialog box, go to the **Protocol's Series** list and select the series to modify.
 - Click the **Edit** button.
The **Edit Series** dialog box appears

Figure 248: Edit Series dialog box

- Make desired changes. Select the **Update** button in the **Edit Series** dialog box.
- Alternatively, click the **Reset** button to remove any changes and return to the previous settings or click the **Cancel** button to return to the previous dialog box.

10. Click the **Update** button on the **Edit Protocol** dialog box.

Figure 249: Update button

- Click the **Begin** button from the **Exam Planner** dialog box.
- When the **Is Scanner Properly Positioned?** popup appears, click the **Continue** button.




Figure 250: Is Scanner Properly Positioned? popup

Note If zero reference is not selected when starting your scan, the scanner considers the last known zero reference point to be the origin and start-point for the next scan. **Always** make sure to zero reference the scanner, when you set up a scan.

The system state orb will change color from yellow to green. The **System Ready to Scan** popup appears.




Figure 251: System Ready to Scan



WARNING **Do not** stand in either the forward or reverse paths of the scanner during the scan.

Note The scanner's side panels permit a low radiation exposure rate of $<0.01\text{mR/sec}/100\text{mAs}$ – when x-ray is emitted.

The **START** button on the scanner control panel turns green when it is enabled.



Figure 252: Scanner control panel – START button

13. Press the **START** button on the scanner control panel to acquire your scan.

The pre-set scan delay countdown begins. The green light turns off when the **START** button is pressed.



Figure 253: Countdown popup

You can press **CANCEL** on the scanner to end the current scan operation. If pressed when lit, the system cancels the current scanning operation. If pressed during scanning, 1 current scan rotation, or 1 second, completes and then the scan is terminated. Alternatively, you can press the **Cancel Scan** button on the screen to cancel the entire scan or **Cancel Exam** button to cancel the entire exam.



Figure 254: Scanner control panel – CANCEL button

Note During the scan, observe the following:

Yellow lights on top of the scanner, and an audible beep identify that radiation is being emitted.

The patient's scan results appear; approximately one image per second.

When scanning begins, the **Continue**, **Repeat**, **Extend**, **Pause Scan**, **Pause Exam**, and **Cancel** buttons are enabled.

When you click the **CANCEL** button, the message "Scan is terminated" appears on both scanner and workstation.

14. If applicable, set your parameters and **Field Of View (FOV)** on your scouts.

Note **FOV** can only be adjusted when two scouts are acquired.



Figure 255: Scouts and FOV button

Scan coverage can be modified by selecting the drag boxes and adjusting the lines and can be centered by clicking on the small green circle and dragging the plan box.

15. Click **Continue** to proceed with your planned scan.



Figure 256: Continue button

The **Pending Scanning Movement** popup appears.




Figure 257: Pending Scanner Movement popup message

16. Click the **Continue** button to scan.
Click the **Cancel** button to cancel the scan.

17. The **System Ready to Scan** popup appears.




Figure 258: System Ready to Scan popup message – System is ready to begin scan

18. The **START** button on the scanner control panel turns green when it is enabled.



Figure 259: Scanner control panel – START button

19. Press the **START** button on the scanner control panel to acquire your scan.

The pre-set scan delay countdown begins. The green light turns off when the **START** button is pressed.

20. If the **Perform Reconstructions** popup appears, do one of the following:

- Click the **Yes** button to perform post reconstructions now.
- Click the **No** button to pause the reconstructions until a later time. When ready, click the Start Recons button.




Figure 260: Perform Reconstructions popup message – To perform post reconstructions

21. Use the **Viewing** tools to review the scan.
See “Examining the scanned image with tools” on page 299.
22. Click the **Finalize** button when finished.
The dose report if **Show Dose Report** in **System Configuration** is enabled appears. In addition, the examination details are saved.

Note You must press the **Finalize** button before you can send the patient’s data to **PACS**.




Figure 261: Dose report

Repeating an image

The **Repeat** function can be used to repeat a scan if necessary. The entire scan can be repeated, or after reviewing the images, a new start position and coverage can be selected if only a portion of the scan needs to be repeated.

1. While the **Acquisition** tab remains active, click the **Repeat** button. The **Protocol Viewer** dialog box appears.




Figure 262: Protocol Viewer – Start Position and End Position

2. Review the protocol parameters.
3. Click the **Repeat** button from the **Protocol Viewer** dialog box.
4. The **Repeat Protocol** popup appears.




Figure 263: Repeat Protocol popup

Note You can change the start and end position or use what appears.

5. Click the **Repeat** button on the **Repeat Protocol** popup.
Scout lines appear in blue, which indicates the second scan and modifications of the start and end points – if made.
The scanner will move to the start position.
6. Press the **START** button on the scanner's control panel to begin the repeat scan.

Extending a scan

The **Extend** function can be used to add additional scan length to **Axial** scans only.

Note If scouts have been acquired, **Extend** can only be used within the boundaries of the Scout. If scouts are not acquired, there are no limits to extending the scan.

1. While the **Acquisition** tab remains active, click the **Extend** button.

The **Protocol Viewer** dialog box appears.

Protocol Viewer

Protocol	
Protocol Name	Brain-Axial
Patient Region	HEAD
Patient Type	ADULT
Patient Position	HFS
Weight (kg)	
Anatomical Reference	Head

Exam Series		
AXIAL BRAIN	Scan Type	Axial
	kV	120
	mA	275
	mAs	275
	CTDIw	65.065
	CTDIvol (mGy)	65.07
	DLP (mGy.cm)	455.49
	Phantom	Head 16cm
	Slice Thickness/Spacing	5 x 5
	Sharpness	Soft Tissue
	Resolution	1 Second(s)
	Start Position	-220
	End Position	-150
	Coverage	70
	Contrast	
Contrast Volume	0	
Delay	10	
Number of Images	14	
Scan Time	7	
Scan Region	Field of View	X0: 173.15 Y0: 173.15 W: 250.00
	Step & Shoot	No
	Body Part Examined	BRAIN

Series Recons	
ST THINS	
BONE	

Extend Close

Figure 264: Protocol Viewer dialog box

2. Review the protocol parameters

3. Click the **Extend** button.

The **Extend Protocol** popup appears.

Extend Protocol

Number of Millimeters

Extend Cancel

Figure 265: Extend Protocol popup

4. Enter the length of the extension in the **Number of Millimeters** text box. Length must be in 10mm increment.




Figure 266: Extend Protocol popup

5. Click the **Extend** button.
6. Press the **Start** button on the scanner's control panel to perform the **Extended** scan.

Scanning with special features

The following features are available for use in protocols.

Using the step-and-shoot option

The **Step & Shoot** option in the protocol lets the user control the start of the scan acquisition. This is helpful in the case of an uncooperative or ill patient where motion is an issue.

1. If necessary, change the **Scan Type** to **Axial**.
2. Click the **Step & Shoot** option in the **Edit Series** dialog box.




Figure 267: Step & Shoot option in the Edit Series dialog box

3. Click the **Update** button in the **Edit Series** dialog box.
4. Click the **Update** button in the **Edit Protocol** dialog box.
5. Click the **Begin** button in the **Exam Planner** dialog box to begin the scan.
The system state orb will change color from yellow to green. The **System Ready to Scan** popup appears.




Figure 268: System Ready to Scan popup

6. To continue the scan, go to the scanner and press the **Scan** button on the screen.
The first set of images are acquired at this position.
The **Step & Shoot** popup appears for you to control the next acquisition.




Figure 269: Step & Shoot popup

7. Click the **Shoot** button to start the scan.
To cancel the scan, click the **Cancel** button.
8. Continue for the length of the scan.
9. Click the **Finalize** button when finished.

Performing a scan with Automatic Exposure Control

Note Depending on system's configuration, not all functions may be available to perform this procedure.

Computed Tomography (CT) is responsible for the largest contribution to the collective effective dose to patients in radiology. The challenge to radiologists and medical physicists is to establish adequate image quality with the lowest radiation exposure to the patient.

Tube current (mA) is one of the key technical scanning parameters for adjusting radiation dose in CT. To optimize radiation dose in CT, users can adjust mA either with manually selected values or with the application of Automatic Exposure Control (AEC). **AEC** refers to the automatic adaption of mA based on user specified image quality and attenuation characteristics of the scanned body region.

Scout scans provide a graph of mA values based on object density and desired noise level. Axial or Helical scans in the protocol can utilize AEC, limiting the mA value of each slice to the minimum necessary to achieve the desired image quality. This ability to modulate the mA throughout the scan to achieve the desired noise level can reduce patient dose.

When using AEC, it is vitally important that the patient is well-centered in the gantry. AEC aims to deliver the specified image quality across a range of patient sizes. The use of AEC may change the planned $CTDI_{vol}$ and DLP values. It tends to increase $CTDI_{vol}$ for large patients and decrease it for small patients relative to a reference patient size.

Note Ensure patient is accurately centered in the gantry.

Do not use AEC when any type of metal is going to be scanned.

Do **not** use **AEC** with a small **FOV**, that is, tiny neonatal patients.

An automatic adjustment of the tube's current cannot occur when the tube potential is changed.

AEC can only be used with the following Sharpness options:

- Helical Soft Tissue – Head
- Helical Soft Tissue – Abdomen
- Axial Soft Tissue
- Axial Pos. Fossa/Vessel
- Axial Sharp

1. Under **AEC**, click the **Enable AEC** option.




Figure 270: Edit Series dialog box with AEC options selected

2. Select the **Minimum mA** dropdown to set the minimum allowed mA value to be used for scanning.
3. Select the **Maximum mA** dropdown to set the maximum allowed mA value to be used for scanning. The available range is 30 to 300.
4. Select the **Noise Level** to set the standard deviation of noise value for the acquired images. The noise range is 1-200.

Note Consult with the site's physicist for guidance specific to the department.

5. Click the **Update** button in the **Edit Series** dialog box.
6. Click the **Update** button in the **Edit Protocol** dialog box.
7. Make sure your patient and scanner are properly positioned.
8. Click the **Begin** button to begin the scan.
9. Press the **START** button on the scanner control panel.
10. After the scouts are acquired and the scan region is set, click the **AEC** tab.
11. Click the **Toggle Graph** button to view the graph on the scout.



Figure 271: Toggle Graph button

The graphs will now appear on the scout(s). Review the mA modulation to ensure it meets your clinical needs.




Figure 272: Graphs on the scout(s)




Figure 273: AEC modulation graph

12. To return to the scout parameter view, click the **Toggle Graph** button, again.
13. If changes to the mA or Noise levels are required, you can modify the Minimum, and Maximum mA and noise as needed on the image.




Figure 274: Minimum mA and maximum mA; noise level

14. When the desired level is achieved according to your department policy, click the **Continue** button to start the scan.
15. Press the **START** button from the scanner.

Note While reviewing the scan you will see mA modulation as per your graphs.





16. Click the **Finalize** button.

Performing a CT angiography scan with Bolus Tracking

CT angiography is a technique that uses contrast to visualize arterial and venous vessels throughout the body. This ranges from arteries serving the brain to those bringing blood to the lungs, kidneys, arms, and legs.

Table 50: Bolus tracking parameters and tools

Option	Description
Bolus Scan Time	The amount of time allowed to monitor the bolus.
Auto-Start	Begins the scan after the specified bolus scan time if no bolus is detected.
Auto-Stop	Stops the scan after the specified bolus scan time if no contrast is detected.

Option	Description
 Test Bolus	A small amount of contrast is injected, and a timing graph is displayed after specified bolus scan time.
 Threshold	Hounsfield Unit detection at the area being monitored – ROI.
 Use Calculated	When performing the test bolus, click the Use Calculated option to use the bolus timing calculated from the test bolus scan.
 Use Manual	Allows you to select a manual timing of bolus after the test bolus has calculated a timing.

1. Perform steps 1 through 8 in “Performing a scan” on page 2644.
2. Click the **Bolus Tracking** option and set parameters, such as **Auto-Start**, **Auto-Stop**, **Bolus Scan Time**, and **Threshold**.




Figure 275: Bolus Tracking option

3. Click the **Update** button in the **Edit Series** dialog box.
4. Click the **Update** button in the **Edit Protocol** dialog box.
5. Click the **Begin** button to acquire the scout(s).
6. Acquire the scout(s).
7. To move the scout or reference line go to the **Viewing** tools and select **Active Scan Region**.




Figure 276: Active Scan Region – Bolus Reference or Helical CTA

Note If the protocol contains one or more scouts, the system places the scan region for each series on the scout based on start and end positions that were entered for each series in the protocol. The **Active Scan Region** dropdown populates the data with the ID that corresponds to each region. You can adjust the scan region prior to continuing the exam, as described below.

8. On the scout, use the Blue lines to plan the helical scan.




Figure 277: Scout line (blue) and Reference line (green)

Note The distance between the **Reference** and **Bolus Tracking** scans cannot exceed 100mm.

9. Position the green **Reference** line at the desired anatomical location.
10. Click the **Continue** button.
The scanner will move to the reference line noted on the scout.
11. Go to the scanner and press the **START** button on the scanner control panel.
12. Click the **ROI** tool and draw the **ROI** on the **Reference** image.




Figure 278: Bolus ROI tool




Figure 279: ROI on the Reference scan

13. Click **Continue**.
14. Load the injector and set your desired flow and rate and arm the injector.
15. When the scanner is ready, press **START** on the scanner and manually start the injector at the same time.
16. The Helical scan will trigger when the threshold value set for the ROI is exceeded.




Figure 280: Scan triggers when bolus enters reference point/ROI

17. Review the completed scan.



Figure 281: Scan at peak enhancement

18. Press the **Finalize** button when complete.

Performing Test Bolus

Test Bolus tests the timing of the injected contrast.

A small amount of contrast is injected, and a timing graph is displayed after the specified bolus scan time. When the contrast is detected, the system stops scanning and a report on the recommended delay-time for CTA protocols appear.

1. Perform steps 1 through 8 in “Performing a scan” on page 2644.
2. Click the **Bolus Tracking** option, click the **Test Bolus** option and set parameters.




Figure 282: Test Bolus option

3. Click the **Update** button in the **Edit Series** dialog box.
4. Click the **Update** button in the **Edit Protocol** dialog box.
5. Click the **Begin** button to acquire the scout(s).
6. Acquire the scouts(s).
7. To move the scout or reference line go to the **Viewing** tools and select **Active Scan Region**.




Figure 283: Active scan region

8. On the scout, use the blue lines to plan the helical scan.




Figure 284: Scan planning lines

9. Position the green **Reference** line at the desired anatomical location.
10. Click the **Continue** button.
The scanner will move to the reference line noted on the scout.
11. Go to the scanner and press the **START** button on the scanner control panel.
The reference image will be scanned and displayed.
12. Click the **ROI** tool and draw the **ROI** on the **Reference** image.




Figure 285: Bolus ROI




Figure 286: ROI on the Reference scan

13. Click **Continue** button.

14. Load the injector and set your desired flow and rate and arm the injector.
15. When the scanner is ready, press **START** on the scanner and manually start the injector simultaneously.
16. The **Bolus Timing** graph appears and shows the calculated, bolus-tracking time.



Figure 287: Bolus timing graph

17. Select from one of the following:

Use Calculated	Uses the bolus timing calculated from the Test Bolus scan.
Use Manual	Allows you to manually set the delay time prior to the start of the helical scan.
Auto Start	Begins the scan after the specified bolus scan time if no bolus is detected.
Auto Stop	Stops the scan after the specified bolus scan time if no contrast is detected.

18. Click the **Continue** button.
Review completed scan.
19. Press the **Finalize** button when complete.

Performing a CT Perfusion Scan

CT perfusion (CTP) is a technique used to evaluate cerebral perfusion of the level of blood flow in the brain, by monitoring the initial phase of iodinated contrast media through the vasculature of the brain.

1. Perform steps 1 through 8 in “Performing a scan” on page 2644.

- After selecting CTP protocol, review the Dynamic CTP options including Scan Time and make your selections.

The screenshot shows the 'Edit Series' dialog box with various parameters for a Dynamic CTP scan. The 'Scan Time' parameter is highlighted with an orange box and is set to 60. The 'Update' button is located at the bottom left of the dialog box.

Figure 288: Edit Series CTP Scan Time

- Click the **Update** button in the **Edit Series** dialog box.
- Click the **Update** button in the **Edit Protocol** dialog box.
- Click the **Begin** button to acquire the scout(s).
- Perform scout and set Dynamic CTP scan location.
- To move the Dynamic CTP location use the green circle inside the reference line.
- Click the **Continue** button.
The scanner will move to the reference line noted on the scout.
- Load the injector and set your desired flow rate. When the scanner is ready, press **Start** on the scanner and manually start the injector at the same time.
- Review completed scan.




Figure 289: Brain Perfusion Image

11. Select **Finalize**.

Calculating and creating perfusion maps

1. Select a perfusion patient from **Patient Browser** and select the series to view.
2. To open the image, click the **View Images** button or double-click the series.

The **Viewing** tab is enabled, and the **CTP** tab automatically opens.










Figure 290: CTP tools

Table 51 CTP Tools

CTP Viewer tools		
	Perfusion Artery/Vein Selection	Select to place the arterial and venous ROIs on the images.
	Calculate CBF, CBV, MTT Map	Select to calculate the CT Perfusion maps.
	Clear Perfusion Map	Cancels the calculations and returns to Calculation mode.
	Show Artery/Vein Flow Graph	Displays the Arterial Venous Flow graph. 
	Peak Image	Displays the image that has the highest HU value based on the arterial ROI placement.

- Click the **Perfusion Artery/Vein Selection** tool.
- Place an **Arterial ROI**.
- Place a **Venous ROI**.
- Click the **Calculate CBF, CBV, MTT Map** tool.
- The Perfusion Maps are calculated and displayed:




Figure 291: Perfusion maps

The calculations produce three maps and the perfusion image:

Perfusion Images	
Cerebral Blood Flow (CBF)	Top, right corner
Cerebral Blood Volume (CBV)	Bottom, left corner
Mean Transit Time (MTT)	Bottom, right corner

8. Clicking the **Artery/Vein Flow Graph** displays the following:




Figure 292: Arterial Venous Flow

This graph displays arterial and venous flow rates with respect to time.

9. Click the **Peak Image** tool to display the **peak image**.
The peak image displays the image that has the highest HU value based on the arterial ROI placement.
10. Click the **Clear Perfusion Maps** tool to cancel calculations and return to **Calculation** mode.
11. Click the **Reset** button to reset images back to the original setting(s).
You cannot undo this action.

Using the Interventional Package

The Interventional Package is designed to make Interventional procedures quick and easy for the technologist and physician. When activated the gantry will spin continuously to make the transition to a scan as fast as possible.

When enabled the **Interventional Tab** will be displayed as seen below.




Figure 293: Interventional Tab

When activated the **Interventional Tab** provides a streamlined workflow specifically for interventional cases.




Figure 294: Interventional Tab - Patient exam details

The Interventional Tab includes the following options:

Option	Description
	Allows a typed in value, when Move is clicked the scanner goes to that location
	Moves the scanner to the Laser location
	Toggles the laser lights on or off
	Allows protocol parameters for the Interventional procedure to be modified
	Initiates X-Ray when the user is ready
	Repeats the last scanned protocol
	Allows a selected raw data set to be Post Reconstructed without leaving the Interventional workflow
	Stops the gantry rotation
	Cancels the currently planned scan
	Ends the current examination

The tab also includes a **Scan Tree**, **Current Scan**, and an updated **Dose Gauge**.




Figure 295: Scan Tree, Current Scan, and Updated Dose Gauge

The **Scan Tree** allows access to acquired scans for quick image loading or repeat scanning using those parameters.

- When a series in the **Scan Tree** is selected, the images associated with it will be displayed, allowing them to be reviewed to ensure the proper start and end locations as well as the thickness and spacing are selected for repeat scans.
- Clicking any of the scans in the **Scan Tree**, then clicking 'Repeat', moves the scanner to the Start location from the selected series and automatically initiates scanning.

The **Current Scan** box shows the protocol parameters used for the most recent active scan.

The **Dose Gauge** displays the accumulated CTDI and DLP in mGy for the current procedure and is updated each time a new scan is initiated.

Viewing Tools, Windowing, and Zoom options are available by clicking the appropriate line below the Dose Gauge.




Figure 296: Viewing Tools, Windowing, and Zoom options

A right click menu on the active image, also allows access to a full array of viewing and measuring tools.




Figure 297: Interventional Workflow

Interventional Workflow

Any protocol can be used as an Interventional protocol.

1. Register a patient
2. Click '**Protocol**'
3. Select the desired protocol
4. Ensure the '**Interventional Scan**' box is selected.




Figure 298: Interventional Workflow - Protocol Information dialogue box

5. Click '**Begin**' to accept the protocol and perform the scout scan(s).
6. Plan your scan over the region of the procedure and acquire that scan.
7. Once the scan has been acquired the system will automatically move to the Interventional tab and show the following:



Figure 299: Interventional Workflow - Scan acquired

8. Review the acquired images from the previous scan to determine the start and end locations of the Interventional scan.



Figure 300: Modify protocol parameters

9. Click '**Protocol**' to modify the protocol parameters for the Interventional procedure.
10. Once parameters are selected click '**Begin**'
The scanner will move to the selected Start position.
11. Click the '**Start**' button to initiate the scans set in the protocol.



Figure 301: Initiate Scans - Interventional protocol

12. When satisfied the anatomy of interest is included in the images, you can enter a location in the **'Go To Position'** box and click **'Move'**.
- This moves the scanner out of the way to allow the physician access to the patient.




Figure 302: Move the Scanner

13. When ready to perform another scan, click **'Repeat'**
- The scanner will move back to the previously defined start location and perform additional scans as defined in the protocol.




Figure 303: Repeat Scans - Interventional protocol

14. You can now use the **'Move'** and **'Repeat'** options to move the scanner back and forth to allow the physician access to the patient and perform repeat imaging to confirm device position in the patient.
15. If changes to the protocol are necessary, click **'Protocol'**.
16. Modify the parameters to the desired values and click **'Begin'**.
17. You must then click **'Start'** to initiate the new scan.

Examining the scanned image with tools

The image tools can only be used when the **Acquisition** tab is enabled and an image is displayed, or when images are loaded to the viewer from the Patient Browser tab.

From the **Acquisition tab**, you can zoom, pan, modify window width, and level, and change layout; see the table below to understand the basics of what each tool looks like and how it performs.

Using tools on the Acquisition tab

The following table describes some of the tools available to you when the **Acquisition** tab is active. For a comprehensive list, see Table 55: 2D, MPR, 3D, and CTP image tools.

Table 52: Image tools








Image tool	Tool name	Action
	Clear Tool	Resets the tool to the default pointer.
	Window Width/Center	Adjusts window width and center of image.
	Zoom	Magnifies the image.
	Pan	Adjusts image on X or Y axis
	Toggle Scouts	Display's or removes scouts from Acquisition.

Image tool	Tool name	Action
	Toggle Layout	Changes the layout to 2x2. Repeat process to return to 1x1.
	Scan Region Re-Draw	If scout lines and the scan region is deactivated, allows you to reactivate.

Chapter 10 Patient Browser










The **Patient Browser** lets you view patient information and images associated with the patient information after the patient's scan.



Figure 304: Active Patient Browser tab

The following table identifies the buttons found on the **Patient Browser**.

Table 53: Command buttons

Button	Action
	Selects the archive destination for selected information.
	Imports exam information from PACS or Media .
	Deletes the selected exam information from the Patient Browser tab.
	Reregisters a patient who is already in the system.
	Generates a Dose Report and Dose SR for the selected patient.
	Combines two different image sets.
	Shows patient, study, series, and image information.
	Displays selected patient images in Viewing .
	Allows you to compare two series.

Navigating the Patient Browser

The **Patient Browser** lets you perform tasks on existing series, for example archiving and previewing the series. This section will introduce you to **Patient Browser** and identify the symbols, areas, and buttons you can use.

The **Patient Browser** can be broken down into the following sections:

- Exam table
- Series table
- Selected protocol table
- Patient Browser's active buttons
- Preview window




Figure 305: Patient Browser sections

Identifying symbols on Patient Browser

Symbols may appear next to series in the exam and/or series tables. These symbols are more vivid when active; they are identified below as active symbols.

Locked 	Indicates the series is locked and cannot be deleted.
Read/Mark 	Indicates the series is marked to be read or read by the physician.
PACS 	Indicates the series has been sent to PACS .
Stored 	Indicates the series has been sent to PACS ; the archived series appear below the initial table
Media 	Indicates the series has been sent to a media device, for example USB.




Figure 306: Patient Browser locked, read, PACS and Stored (archived), and media symbols

Using the vertical and horizontal scroll bars on Patient Browser

Navigation scroll bars let you move the lists in the Patient Browser sections up, down, and horizontally to view all available exam information.




Figure 307: Patient Browser horizontal and vertical scroll bars

Locking a study

1. Click the **Patient Browser** tab.
2. Select the study to lock.

3. Right-click the mouse button.
4. Click **Lock** on the floating menu.



Figure 308: Floating menu - Lock

The **Lock** symbol appears for any selected series.



Figure 309: A locked series

A series cannot be deleted while in the **Lock** mode.

5. To unlock right-click and click **Unlock**.

Note All **QA** series are locked to prevent deletion. The **QA** series can only be unlocked by your field service representative.

Marking a series to read

1. Click the **Patient Browser** tab.
2. Select the study to mark.
3. Right-click the mouse button.

- Click **Mark** on the floating menu.



Figure 310: Floating menu - Mark
The **Mark** symbol appears for any selected series.

- To unmark any series, right-click **Unmark**.

Using the preview window

- Click a patient in the exam table.
- Click a series in the series table.
- Click the **Preview** button to the right of the **Series Protocol Table**.



Figure 311: Preview Button

- The selected series will appear in the **Preview** window.



Figure 312: The series appears in the preview window

Archiving patient series

You can archive patients and studies (or series) to **PACS**, media (USB or CD), or surgical navigation devices.

Archiving to PACS

1. Click the **Patient Browser** tab.
2. Select the patient study for **PACS** in the following way:
 - To select one patient and all associated series, click the patient, and click the **Archive** button.
 - To select specific series for a patient, press and hold the **Ctrl** key, then click each individual series from the **Series** table, and click the **Archive** button.

The **Archive Destination** popup appears.




Figure 313: Archive Destination popup

3. Click the **PACS** button.
The **Archive to Server** popup appears.
4. Click the **Select Archive Location** dropdown and select a site.




Figure 314: Archive To Server popup


5. Review the **Image Range** items to make sure all those items you selected in step 2 are captured.
6. Click the **Archive** button to begin the archive process.
- If enabled the **Store/Print Queue** dialog box will appear to show the status of your image transfer. You can also activate the **Store/Print Queue** dialog box by clicking **Tools > Store/Print Queue** from the main menu.








Figure 315: Store/Print Queue dialog box

7. Watch the status of each series:
- **Pending** informs you that the series is paused because you clicked the **Pause** button.
 - **Connecting** informs you that the series is in process of archiving to its targeted location.
 - Each series will move from the top portion of the popup to the bottom portion of the **Store/Print Queue** popup when it has been processed.
8. While the archiving is in process, you can perform one of the following from the buttons in the **Store Print Queue** dialog box.

Table 54: Store and Print Queue buttons

Store and Print Queue button	Action
	When you select one or more series, temporarily stops the series from being stored. This is a toggle button with the Resume button.

Store and Print Queue button	Action
	When you select one or more series, deletes either a series to be stored, or a series that failed to store.
	Stops the archive to USB or a drive.
	When you select one of more series, tries to archive the selections.
	When you select one of more series, displays an explanation of why a series failed to store.
	Closes the Store/Print Queue popup.

9. If the series is not successfully stored to its targeted destination, the “Store Failed” message appears in the **Failure** column. This means the series was not successfully archived.
10. If there are failed archived series; click the **Retry** button to attempt to archive the series you selected.

Note Any **Storing Failure** status appears in the bottom of the popup to inform you why the failure occurred. If an archive job fails, it will be sent to the **Failed Store/Print Jobs** list.

11. When the archiving is complete, click the **Close** button to exit the **Store/Print Queue**.

You can also click the **Close** button and the archiving process will continue as you do other tasks.

Archiving to Media

1. Click the **Patient Browser** tab.
2. Select the patient study to archive following way:
 - To select one patient and all associated series, click the patient and click the **Archive** button.

- To select specific series for a patient, press and hold the **Ctrl** key, then click each individual series from the **Series** table, and click the **Archive** button.

The **Archive Destination** popup appears.




Figure 316: Archive Destination popup

3. If you are archiving a USB device, insert the USB drive into the USB port.
4. Click the **Media** button.

The **Archive to Media** popup appears.




Figure 317: Archive to Media popup

5. Click the targeted drive and path destination.
6. The **Archive** button is active.




Figure 318: Archive Button active

7. Click the appropriate check boxes for your archive process:

Anonymous Patient	Makes the patient's information anonymous for HIPAA standards.
Anonymous Institution	Makes institutional information anonymous for HIPAA standards.
Include CD Viewer (requires CD viewer software installed)	Includes a CD viewer application to view images from the media.
Export as JPEGs	Exports image files in .JPG format.
Show Demographics	Includes the demographic information in archive if you clicked the Export as JPEGs check box.

8. Click the **Archive** button to begin the archive process.
- The **Cancel** button is active after clicking the Archive button; click the **Cancel** button to stop the archive.
9. The **Archive to Media** dialog will update the status when archiving is complete.




Figure 319: Archiving complete

- Click the **Refresh** button to remove any messages that appear in the **Message** box.
- Click the **Minimize** button to continue working in other areas while the archiving process runs.
 - A disk appears at the bottom; click it to maximize the **Archive to Media** popup.
- Click the **Close** button to exit the **Archive to Media** popup after the archive process is complete.

Archiving to Navigation

1. Click the **Patient Browser** tab.
2. Select the patient study or series.
3. Click the **Archive** button.
4. Click the **Navigation** button.
5. Click the **Select Archive Location** dropdown and select the location.




Figure 320: Archive to Server popup

6. Under **Image Range**, select the image(s) you want to send to navigation.
To return to the default selections, click the **Reset** button.
7. Perform one of the following:
 - Click the **Archive** button to send the image to **Navigation**.
 - Click the **Cancel** button to return to **Patient Browser**.

Import

Import allows you to add patient images to the Patient Browser.

Importing from PACS

1. Click the **Patient Browser** tab.
2. Click the **Import** button to import data.
The **Import Location** popup appears.




Figure 321: Import Location popup

3. Click the **PACS** button.
The **Import from PACS** dialog box appears.




Figure 322: Import from PACS dialog box

4. Click the **Query** button.
The **Query Information** dialog box appears.




Figure 323: PACS Query Information dialog box

5. Perform the following:
 - Select a **HIS/RIS** server from the dropdown.
 - Set the values to search in your query.
 - Click one of the following buttons:
 - Click **Query** to save the search results
 - Click **Reset** to clear the query information.
 - Click **Cancel** to exit the **Query Information** popup.
6. From the **Queried** results, select a patient and click the **Get Series** button.




Figure 324: Import PACS dialog box with active Get Series button

- Click the **Import** button.
- Click the **Close** button to exit the **Import from PACS** dialog box.

Importing from media

- Click the **Patient Browser** tab.
- Click the **Import** button on **Patient Browser**.
The **Import Location** popup appears.



Figure 325: Import Location popup

- Click the **Media** button.
The **Import from Media** popup appears.




Figure 326: Import from Media popup

4. Click the drive and path where images were previously stored.
The **Import** button is active.




Figure 327: Active Import button

- If necessary, click **Subfolders** to see the entire path.

5. Click the **Import** button.
The imported images appear in **Patient Browser**.

Delete

1. Click the **Patient Browser** tab.
2. Select the study or the series to delete.

- Click the **Delete** button.
The **Confirm Deletion** popup appears.




Figure 328: Confirm Deletion popup

- Click the **Delete** button on the **Confirm Deletion** popup.
The patient data will be deleted from the **Patient Browser**.

Registering a patient from Patient Browser

If additional scans must be performed on a patient that is listed in the **Patient Browser**, you can register them by performing the following:

- Click the **Patient Browser** tab.
- Select the patient to register.
- Click the **Register** button.



Figure 329: Patient browser register button

- The **Create New Study** popup appears.




Figure 330: Create New Study popup

- Perform one of the following:

- Click the **Patient/Study** button to use all exam information selected, including the accession number.
- Click the **Patient** button to use only patient information.
- Click the **Cancel** button to exit the **Create New Study** popup.

Building dose from Patient Browser

The **Build Dose** button in the Patient Browser, allows you to manually create a Dose Report and Dose SR image which will appear in the Series table when completed.

1. Click the **Patient Browser** tab.
2. Select the patient to use for the **Build Dose**.
3. Click the **Build Dose** button.




Figure 331: Build dose button

4. The **Build Dose Please Wait** popup appears.




Figure 332: Please Wait popup

5. The Dose Report and Dose SR image will be saved to the Series Table.
6. If the **Dose Build Failed** popup appears, click the OK button, make the necessary changes, and try again.



Figure 333: Dose Build Failed popup

Note **Dose Structured Reports (Dose SR)** cannot be viewed in the BodyTom Elite system; **Dose SR** can be viewed in **PACS** with the appropriate viewer.

Merge

The merge function is used to combine two different series for review.

Note Only series with same patient name, ID, slice thickness, slice spacing, kernel, and pixel spacing can be merged.

1. Click the **Patient Browser** tab.
2. Select the study to merge.
The study must include several scans.
3. Select two series to merge by pressing and holding the **Ctrl** key and highlighting the two series.
The **Merge** button appears and is active.



Figure 334: Two series selected to merge

4. Click the **Merge** button to merge the series.




Figure 335: Active Merge button

5. The **Please Wait** popup appears while the application loads the series to be merged.
6. When the **Merge** window appears, the selected series are displayed side by side in separate windows.
7. The bottom 'filmstrip' is empty when the window first displays.



Figure 336: Image viewer without image(s) selected

8. Images can be drag from the series displays on the sides of the **Merge** window to the Image Viewers.
 - These images will be outlined in green in the series list.
 - Multiple images in the series lists can be selected by clicking each image.
 - You can also right click on the series lists and choose **Select All**.
 - To Deselect a single image, click the image, the green outline will change to white.
 - You can also right click on the series lists and choose **Deselect All**.

9. Drag the images you wish to merge from the series lists to the bottom 'filmstrip'.




Figure 337: Selecting images to move to image viewer or bottom filmstrip

10. To reorder the filmstrip, drag one or more images to the new location.




Figure 338: Finished filmstrip

11. When finished with the bottom strip, click **Merge**.

The application generates an additional series with the selected images. The **Series Save** popup appears indicating if the merge operation was successful.



Figure 339: Merge Series Saved popup

12. Click the **OK** button.

13. Click **Close**.

The merged series will display in the **Patient Browser** and will be marked as a merged series under the Series Description column.

Using Show Info to view, update, and move a series

1. Click the **Patient Browser** tab.

2. Select the patient.

3. Select the series.

The **Show Info** button becomes active.

4. Click the **Show Info** button.

The **View/Update Information** dialog box appears.




Figure 340: View/Update Information dialog box

The following tabs appear:

Patient	Data about the patient, such as patient name, date of birth, and sex of patient.
Study	Data about the study, such as date, time, and referring physician.
Series	Data about the series, such as the position reference indicator, model, pixel-padding value, and series date and time.
Images	Data about the image, such as the instance number (sequential), exposure time.

5. Click the tab(s) to review and update the necessary information.
6. Double-click any editable field and make your change(s).
7. Click one of the following buttons

- Click the **Update** button to save your changes.

Note If information is invalid, you are prompted to correct the information and click **Update** again.

If the field cannot be edited, a prompt appears tell you the selected field is not editable.

- Click the **Move** button to show the **Move Series** popup.

Note The **Move Series** function is used when a scan has been acquired under the wrong patient file or to move a patient that was registered manually to the **Patient Registration** tab. Moving the patient to the **Patient Registration** tab will capture all the patient's information.

- Click the **Registration** button to move the patient into an existing patient or by manually creating a new patient using the **Patient Registration** process.
- Click the **Cancel** button to exit the **Move Series** popup.
- Click **Browser** to go to the **Patient Browser** tab and move the series.
- Click the **Move** button to confirm moving the series.
- Click the **Cancel** button to exit the **Patient Browser** tab and return to the **View/Update Information** popup.

8. Click the **Close** button to exit the **View/Update Information** popup.

Note An audit log of both old and new patient series, including the date and time of change and who performed it, is generated.

Modifying a series scanned under the wrong patient

If a patient has been scanned under the wrong identification, the series can be corrected.

1. Click the **Patient Browser** tab.
2. Select the series that was scanned with incorrect patient identification to modify the data.
3. Click the **Show Info** button.

The **View/Update Information** dialog box appears.



Figure 341: View/Update Information dialog box

4. Click the **Move** button.
- The **Move Series** popup appears, denoting where to retrieve patient information from.




Figure 342: Move Series popup

The following defines what each button performs:

Registration	If patient information is stored within hospital's HIS/RIS server, click the Registration button, which will open the Register Patient tab to let you choose patient/study information.
Browser	If patient information is stored within system's browser, click Browser button, which show the Patient Browser tab information to let you select a series with correct patient information.
Cancel	Returns you to the previous dialog box.

5. Perform one of the following:
 - If you clicked the **Registration** button in the previous step, go to the next step.
 - If you clicked the **Browser** button in the previous step, go to step 11.
6. Click the **Manual** button.
The **Exam Information** dialog box appears.
7. Enter the corrected data in any of the fields.
See "Registering the patient" on page 253.
8. Click the **Update** button to save the change(s).
9. Click the **Move** button.
A prompt appears to review changes made to the patient and/or series information for changes to take effect.
10. Click the **Ok** button and then the **Update** button.
The corrected patient and moved data will appear in the **Patient Browser**.
11. If you selected the **Browser** button, the **Patient Browser** tab is showing; select the correct patient and series.
12. Click the **Move** button.
A prompt appears to review changes made to the patient and/or series information for changes to take effect.
13. Click the **Ok** button.
14. Review the patient to ensure it is the proper one.
15. Click the **Update** button.

16. Click the **Cancel** button to return to the **Patient Browser**.

Loading a series into view

1. Click the **Patient Browser** tab.
2. Select the patient.
3. Select the series.
4. Click the **View Images** button or double-click the selected series. The **Viewing** tab opens, and the series appears for viewing and manipulating.

Appending a series

Note Regardless of how many series are appended, the series are listed chronologically. This tool can be used to put all images from a patient together on a CD or to PACS.

1. Click the **Patient Browser** tab.
2. Select the study to append.
3. Select the first series.
4. Right-click the mouse to select the second series. The **Append Images** appears on the floating menu.




Figure 343: Floating menu - Append Images

5. Click Append Images.
The **Please Wait** popup appears.



Figure 344: Please Wait popup

A new series is created with (Appended) at the end of the description.

PACS	Stored	Media	Series Date	Series Time	Modality	Series Description	Protocol	Pos
✓	✓	✗	04/27/2017	2:12:21 PM	CT	Helical Chest	Adult Chest Helical	HFS
✓	✓	✗	04/27/2017	2:12:21 PM	CT	Helical Chest (Appended)	Adult Chest Helical	HFS
✓	✓	✗	04/27/2017	2:02:45 PM	CT	Dose Report	Dose Report	HFS
✓	✓	✗	04/27/2017	2:02:44 PM	SR	Dose SR		

Figure 345: (Appended) series created

Chapter 11 Viewing Images

Viewing lets you see already-scanned images from previous examinations. To view images, select the patient in **Patient Browser** and then select the series to view. To open the image, click the **View Images** button or double-click the series.



Figure 346: Active Viewing tab

The following tables identify the tools in the **Viewing** tabs that let you manipulate images. Some image tools appear on specific viewing tabs, only. The view tabs are **2D**, **MPR**, **3D**, and **CTP** (if enabled).

Table 55: 2D, MPR, 3D, and CTP image tools









Image tools	Tool name	Action
Common tools		
	Clear Tool	Resets the tool to the default pointer device.
	Window Width/Center	Adjusts window width and center of image.
	Zoom	Magnifies the image.
	Pan	Adjusts image on X or Y axis.
	Invert	Inverts black to white and white to black.
	Capture	Saves a screen capture of a selected viewport.
	Capture All	Saves screen captures of all visible viewports.
	Reset	Reverts all images back to their original mode.
2D and CTP tools		





























Image tools	Tool name	Action
	Region of Interest (ROI)	Defines a circular ROI and displays the ROI information.
	Arrow	Draws an arrow on the image.
2D, CTP, and MPR tools		
	Line	Draws a line on the image and is used for measurement.
	Angle	Draws an angle on the image and displays the angle information.
2D only tools		
	Add Annotation	Create text box for annotation.
	Rotate	Rotates images.
	Reverse Image Stack	Reverses the order in which images display.
	Flip Vertically	Flips images up or down.
	Flip Horizontally	Flips images right or left.
	Cine Reverse	Cine backward through the image.
	Cine Forward	Cine forward through the image.
	Stop	Stop the cine loop.
MPR only tools		
	Tilt	When selected a White 'steering' wheel allows you to correct a rotated image.
3D only tool		

Image tools	Tool name	Action
	Color Preset	Dropdown menu allows you to select from multiple color options.
	Render Mode	Dropdown menu allows you to display images in Color, MIP, or Grayscale.
	Orientation	Dropdown menu that allows you to select from multiple orientation options.
	Rotate	Rotates the 3D image.
	Undo	Reverses the most recent action taken.
	Redo	Restores the most recent Undo action taken.
CTP only tools		
	Perfusion Artery/Vein Selection	Select to place the arterial and venous ROIs on the images.
	Calculate CBF, CBV, MTT Map	Select to calculate the CT Perfusion maps.
	Clear Perfusion Map	Cancels the calculations and returns to Calculation mode.
	Show Artery/Vein Flow Graph	Displays the Arterial Venous Flow graph. 
	Peak Image	Displays the image that has the highest HU value based on the arterial ROI placement.

Using keyboard shortcuts

Keys are a quick way to navigate around. The table below provides keyboard shortcuts you can use to manipulate images in the **Viewing** tab.

Table 56: Arrow key navigation

Arrow keys	Action
	To scroll through images.
	To adjust the window center.
	To quickly scroll through images.

Setting window width and center

Note Any modifications you make are not saved to the image.

1. Select a patient from **Patient Browser**, select the series to view.
2. To open the image, click the **View Images** button or double-click the series.
The **Viewing** tab is enabled and the **2D** viewer opens.
3. The following options allow you to adjust the window width and center of the image:
 - Click the Window Width/Center icon in the Tools menu, then while holding the left mouse button down drag up/down to modify Window Center and right/left to modify Window Width.
 - To adjust with a preset, click the **Windowing** dropdown and select a preset.




Figure 347: Windowing preset dropdown list

- Type values in the **Width** and **Center** text boxes and click the **Apply** button.




Figure 348: Windowing Width and Center text boxes, and the Apply button

- Right click over an image and use the **Activate Window Tool** option then while holding the left mouse button down drag up/down to modify Window Center and right/left to modify Window Width.




Figure 349: Right click menu

Viewing images in 2D

2D lets you view scanned images in a **2-Dimensional** space. Standard **2D** mode is used when *only* one dataset is loaded. The default layout is a 2 x 2 grid.

The **Viewing** tab and **2D** viewer opens when you select a dataset from the **Patient Browser**.

1. Select a patient from **Patient Browser** and select the series to view.
2. To open the image, click the **View Images** button or double-click the series.

The **Viewing** tab is enabled and the **2D** tab is opened.



Figure 350: 2D tools

3. Use any of the image tools to manipulate your images.
4. Click the **Reset** button to reset images back to the original setting(s). You cannot undo this action.

Comparing images

You can compare images in two different ways:

Note You can compare two series from the same patient or two series from different patients.

Using the floating menu to compare images

1. Select the patient in **Patient Browser**.
2. Select the first series from the series window.
3. Right-click and click **Mark for Compare** from the menu.



Figure 351: Floating menu - Mark For Compare

4. Select the second series or a series from a different patient.
5. Right-click and then click **Compare With Selected Series** from the floating menu.



Figure 352: Floating menu - Compare with Selected Series

Both series are loaded into **Viewing** to compare.




Figure 353: Compared series

6. Click the **Link** button to link both images together to view.



Figure 354: Link button

The **Unlink** button replaces the **Link** button.

7. Click the **Reset** button to reset images back to the original settings.

Using the Compare button to compare two images

1. Select the patient in **Patient Browser**.

2. Select the first series.
3. Press and hold the **Ctrl** key.
4. Select the second series.
Both series are highlighted.
5. Click the **Compare** button.



Figure 355: Using the Compare button

Comparing a scout and a scan

1. Select a patient from **Patient Browser**.
2. Select a scout from the series window.
3. Press and hold the **Ctrl** key on the keyboard.




Figure 356: Scout and scan selected to compare

4. Select the scan from the series window.
Both images are highlighted.
5. Click the **Compare** button.
The scout and the scan will appear on screen at the same time. A green localizer line appears on the scout.

6. Compare the scout to the scan.




Figure 357: Comparing a scout (1) and a scan (2)

7. To remove the scout, right-click on the scout and click **Hide Scout Viewpoint**.

To return the scout to view, right-click in the viewing window and click **Show Scout Viewport**.

Using the ROI

1. Select a patient from **Patient Browser** and select the series to view.
2. To open the image, click the **View Images** button or double-click the series.
The **Viewing** tab is enabled and the **2D** tab is opened.
3. Click the **ROI** tool.
4. Move the mouse pointer to the image where you want the **ROI** located.
5. Click the left-mouse button and drag the **ROI** diameter to the required size. To lock the **ROI** in place, click the left-mouse a second time.
 - To change the location of the ROI or the details of the ROI, click the **ROI** or measurements you wish to move. The ROI and its measurements will turn yellow, and the pointer becomes a hand. Click and hold the mouse button on either the ROI or its measurements and drag to a new location. Click anywhere outside the ROI to freeze it in the new location.
 - When you move the **ROI** to a different location the measurements of the ROI are automatically updated based on the new location.

- To remove the **ROI**, left click anywhere on the **ROI**, right click to see the floating menu, and click **Delete Annotation**, or click on the **ROI** and when it turns yellow, press **Delete** on the keyboard.




Figure 358: ROI

Using Layout and Rotate in 2D view

Layout lets you alter the number of images presented on the **Viewing** tab. **Rotate** lets you turn the images. Select a patient from **Patient Browser**, select the series to view.

- Select a patient from **Patient Browser** and select the series to view.
- To open the image, click the **View Images** button or double-click the series.
The **Viewing** tab is enabled and the **2D** tab is opened.
- To adjust the layout of the viewing area, click the **Rows** and/or **Columns** dropdowns to select the number of rows or columns you want to show.




Figure 359: Layout (viewing tools)

- To rotate the image, click the **Rotate** dropdown and select the number of degrees to rotate the images.




Figure 360: Rotate dropdown

5. Click the **Rotate** button to see the images turn to the new angle.




Figure 361: Rotate (viewing tools)

Applying a grid to your images in 2D

1. Select a patient from **Patient Browser** and select the series to view.
2. To open the images, click the **View Images** button or double-click the series.
The **Viewing** tab is enabled and the **2D** tab is opened.
3. Change the size of the grid in the **Grid (mm)** text box.




Figure 362: Grid (mm)

4. Click the **Apply** button to apply a grid over the image.
5. Perform one of the following to remove the grid:

- Click the **Reset** button.
- Enter 0 in the **Grid (mm)** text box and click the **Apply** button.

Viewing images in MPR

Multi-Planar Reformation (MPR) allows images to be created from the original **Axial** plane into **Coronal**, **Sagittal** or **Transverse (Axial)** planes. **MPR** is fast, uses all the attenuation values in the dataset, and can be easily performed on the workstation. **MPR** however, provides on a **two-dimensional (2D)** display of the image data.

Viewer layout is 2 x 2 as seen below.

1. Select a patient from **Patient Browser** and select the series to view.
2. To open the images, click the **View Images** button or double-click on the series.
The **Viewing** tab is enabled and the **2D** tab is opened.
3. Click the **MPR** tab.
The **MPR** screen appears.



Figure 363: MPR tools

4. Select the image reformat at the bottom of the screen.



Figure 364: Image reformat selections

5. Use any of the image tools to manipulate the images.
6. The tilt tool can be used to modify the rotation of the images.
7. Adjust the image angle by moving the circle.




Figure 365: Tilt tool

8. Click the reset button to reset the images back to the original settings.
You cannot undo this action.

Understanding and using slab

Through the reformation process, axial images are stacked creating a volume, or slab, which can be assessed in different planes. The thickness and spacing of each slab can be varied to meet the needs of the viewer. The reformations can be displayed in an average, maximum or minimum projection.

MPR's should be created using 1.25mm slices with a spacing of 0.625mm whenever possible.

Creating the slab

1. Select a patient from **Patient Browser** and select the series to view.
2. To open the images, click the **View Images** button or double-click the series.
The **Viewing** tab is enabled and the 2D tab is opened.
3. Click the **MPR** tab.
The **MPR** screen appears.
4. Click the **Sagittal**, **Coronal**, or **Transverse** plane to create your slab.




Figure 366: Image formats

5. The **Secondary Series** option is enabled.
6. Select **Enable Slab**.




Figure 367: Enable Slab option

The **Enable Slab** option is inactive if no **MPR** view is selected.



Figure 368: Enable Slab option under Secondary Series

7. Set the **Cyan** lines to determine the beginning and end of the slab.



Figure 369: Cyan Line and cyan circle to drag for FOV

8. Use the **Cyan** circle to drag the planned slab if required.
9. Define the **Slab Thickness** and **Slab Spacing** in the text boxes.



Figure 370: Slab Thickness and Slab Spacing text boxes

You can also use the **yellow** squares found on the slab thickness display to modify the **Slab Thickness**.



Figure 371: Small yellow boxes to manually drag for desired slab thickness

10. Click the **Slab Rendering Options** dropdown to select the appropriate option.




Figure 372: Slab Rendering Options dropdown

The following options are available in MPR Slab mode:

Slab Thickness	The thickness of the MPR slab.
Slab Spacing	The space between the start of one slab and the next.
Slab Rendering Options	Where you define the pixel values that will be displayed in each slab: options include, Average, Maximum Intensity and Minimum Intensity.

Maximum Intensity	The highest pixel values for all slices within the slab is displayed.
Minimum Intensity	The lowest pixel values for all slices within the slab is displayed.
Average	The pixel values of all slices within the slab are combined and the average value for each pixel is displayed.
Series Description	Text field for naming the series of images created when clicking the Generate button.
Yellow lines	Define the slab thickness. The boxes on the lines allow you to adjust the thickness using the mouse.
Cyan lines	Define the slab FOV and dictate the range of the new series to be generated. The cyan lines are adjustable by clicking and dragging on the lines themselves; both lines are moved by clicking and dragging the central circle marker.
Red, blue, and green lines	Define the cross sections of the anatomy being viewed.
Generate	Generates a new series with the name given in the Series Description field, based on the selected MPR view pane.

11. Select the **Tilt Tool** to correct any rotation on the image.




Figure 373: Tilt tool

12. Use the mouse pointer to move the white **Tilt** circle.

Note The circle does not represent the Field of View




Figure 374: Tilt white circle

13. Enter the slab name in the **Series Description** text box.
14. The slab can be previewed in the bottom right viewport.

Note Make sure **Zoom** is at 100% or below.

Ensure all expected anatomy is included when previewing the created MPR.

15. When you are ready to save, click the **Generate** button.
A **Saving Series** popup appears.

When the series is complete, the **Capture Complete** pop-up appears.



Figure 375: Capture Complete popup message – Series saved

16. Click **OK** to close the **Capture Complete** popup message.
17. The new **MPR** images appear in the Patient Browser with the description in the Series Description text box.

CT	50cm Helical	Adult Abdomen Helical
CT	Coronal Abdomen	Adult Abdomen Helical
CT	Lateral Scout	Adult Abdomen Helical

Figure 376: MPR images in Patient Browser

18. To create additional MPR's, select the Reset button in MPR mode, select the MPR view you want to create and perform the steps above to create the new view.

Viewing images in 3D

In **3D** viewing, a 3-Dimensional image is created by stacking all the images of a scan on top of one another to create a 3D-volume. The initial display shows the **3D** volume and a box appears around it. **MPR** planes also appear.

1. Select a patient from **Patient Browser** and select the series to view.
2. To open the images, click the **View Images** button or double-click the series.
The **Viewing** tab is enabled and the 2D tab is opened.
3. Click the **3D** tab.




Figure 377: 3D tools

4. To rotate the image up to 360°, click **Rotate** and move the image with the mouse pointer to the rotation of choice.
5. You can change the **Color Presets** from the dropdown menu:




Figure 378: 3D Color Presets

6. You can change the **Render Mode** from the dropdown menu:




Figure 379: 3D Render modes

7. Click the **Orientation** drop-down box to assign an orientation:



Figure 380: 3D Orientation options

8. Click the **Reset** button to reset images back to the original settings.
You cannot undo this action.

Chapter 12 Post Reconstruction

The system stores multiple patient series of raw data to allow post reconstruction of images. **Post Reconstruction** allows reconstructing of the acquired data using different algorithms, slice thicknesses, or use of image enhancement algorithms, such as **Metal Artifact Reduction**, **Noise Reduction** and **Windmill Correction**.

Reconstruction Overview

Metal artifact reduction

Streak artifacts are often seen around metal leads, prostheses, applicators, bone, or metal screws. Numerous factors can contribute to these streaks including under-sampling, photon starvation, patient motion, beam hardening, and scatter. You can use **Metal Artifact Reduction (MAR)** to reduce these streaks. **MAR** removes the metal from the image to reconstruct the soft tissue only; then it adds it back, to reduce the artifacts. This is currently used only on **Axial** scans.

Noise reduction

Noise appears as grain on the image and is caused by a low signal to noise ratio. This occurs more commonly when a thin-slice thickness is used. It can also occur when the radiation dose is insufficient to penetrate the anatomy being scanned.

Note Noise reduction applies to post-processing filters that reduce the amount of noise in the images. In clinical practice, using noise reduction may allow for a reduction in CT patient-dose depending on the clinical task, patient size, anatomical location, and clinical practice. Consult with the site's radiologist and physicist to determine the appropriate dose to obtain diagnostic image quality for a particular clinical task.

Windmill Correction











The Windmill Correction reduces artifacts that are common in **Helical** scans due to the nature of the cone-beam reconstruction.



Figure 381: Active Post Reconstruction tab

The tools available to **Post Reconstruction** are identified in the table below.

Table 57: Reconstruction tools

Image tools	Tool name	Action
	Load Images	Loads images from selected series into viewing.
	Stop	Cancels the current, post-reconstruction request. All images are generated until you click the Stop button.
	FOV	Adjusts the FOV prior to reconstruction.
	Clear Tool	Resets tool to default pointer device.
	Window Width/Center	Adjusts the width and level of selected image.
	Zoom	Magnifies the image.
	Pan	Adjusts the image on X or Y axis.
	Reset	Resets the display to default viewer settings.
	Start	Begins your Post Reconstruction .
	Resend	Sends the last acquired scan from the recon workstation to the Patient Browser .

Performing Post Reconstruction

The following figure identifies parts of **Post Reconstruction**:




Figure 382: Post Reconstruction areas

1. Click the **Post Reconstruction** tab.
2. Select a study in the **Exam Table**.
When you select a study, all the scanned series for that study appear in the **Series Table**.
3. Select the series to reconstruct.




Figure 383: Post Reconstruction study and series tables

4. Click the **Load** button.
The scan or series will load into the viewer. The series protocol table and the **Reconstruction Parameters** and **Options** tabs are active. The scan will appear in the **Image Viewing Area**.
5. View the study in the Image Viewing Area.




Figure 384: Post Reconstruction viewing image area

6. To modify the **FOV**, perform the following:

- Click the **FOV** tool, click and drag the mouse to form a square on the image in the **Image Viewing Area**.
 - The size of the square appears in the Width box in the **Reconstruction Parameters and Options** tab.
- Click the circle in the middle of the FOV square and drag to move the **FOV**.
- Click one of the two drag boxes on the corners to adjust the size of box.
 - The **Width** dimension, and **X/Y** coordinates are adjusted as the size changes.
 - The **FOV** size cannot exceed the range of 50 - 600mm square.




Figure 385: FOV resizing boxes

Note You can also enter a number in the **Width** or X and Y box to define a specific **FOV**.

Alternatively, click the **Use FOV Max** option to use the maximum **FOV**.

- Click the **Sharpness** dropdown to select a reconstruction algorithm from the **Reconstruction Parameters** tab.




Figure 386: Sharpness on the Reconstruction Parameters tab




Figure 387: Reconstruction Parameters Sharpness dropdown

- Click the **Slice Thickness/Spacing** dropdown to select.




Figure 388: Reconstruction Parameters Slice Thickness/Spacing dropdown

The slice thickness and spacing options available are determined by the type of scan that was acquired (Axial vs. Helical).

- The **# of Expected Images** text box shows the calculated number of images that will be reconstructed based on the parameters used for the reconstruction.




Figure 389: # of Expected Images

10. Click the **Options** tab. The following are only available during **Post Reconstruction**.

11. Perform the appropriate action:

- If desired, select **Noise Reduction** for an **Axial** or **Helical** scan.



Figure 390: Noise Reduction on the Options tab for a Helical scan

- If desired, select **Perform Windmill Correction** for a **Helical** scan.



Figure 391: Perform Windmill Correction and/or Noise Reduction on the Options tab for a Helical scan

- If desired, select **Metal Artifact Removal** for an **Axial** scan.




Figure 392: Metal artifact removal

12. Click the **Start** button to generate a new dataset.

When you click the **Start** button, the reconstructed images appear in the viewing pane.




Figure 393: Please wait while the system performs data reconstruction message

13. When the reconstruction is complete, the images appear in **Patient Browser**.

Resending images from the scanner to the workstation

Pressing the Resend button lets you send the last acquired scan from the recon computer to the **Patient Browser**. This may be necessary when you have wireless interruptions and/or workstation shuts down unexpectedly.




Figure 394: Resend button

Chapter 13 Accessories and Options

In this chapter you will learn how to convert a bed, stretcher, or any type of adjustable surface into a scanning platform using the Universal Transfer Board.

To request the catalog(s) to reference product descriptions/details and part numbers for the available accessories/options that are used with the BodyTom Elite, see “Contact information” on page 25:

When using a fixed scanner, the table moves from one portion of anatomy to another while the gantry remains stationary. With the BodyTom Elite, an in-place scanning platform remains stationary while the gantry or scanner translates from one point to the other to cover the anatomy.

The universal transfer board can be used for most beds or stretchers. It is placed under the patient and secured to the bed or stretcher with straps.



WARNING NeuroLogica Corp. recommends that the weight of the patient being positioned on the scan board does not exceed the bed manufacturer’s safe, recommended, operating patient load. Realizing patient safety is of the utmost importance, it is recommended that safe judgment be exercised at all times when it comes to the clinical care of patients. There are a number of varying factors, such as the condition of the bed being used, unique patient anatomy, as well as the proper scan board and positioning of the patient, per NeuroLogica Corp.’s clinical training guidelines and product labeling. If any excessive wear or damage is noticed to any scan board, do not use it for a patient scan; contact a qualified service technician to assess, repair, and/or replace the device.

Using the Universal Transfer Board

The universal transfer board is a carbon-fiber, radiolucent board that is designed to work with any ICU bed or stretcher. The carbon-fiber board comes with a 0.5-inch-thick headboard and 2-inch x 5-foot straps to strap the board to the ICU bed or stretcher.

You can use the universal transfer board on any bed, table, or stretcher. Because you can attach the universal transfer board to almost any type of surface, it is used anywhere throughout the hospital including the ICU, OR, and ER. The universal transfer board is placed on the mattress and secured with a strap or placed directly on a surgical table under the cushions. The patient lies on the board with the patient's head in the head holder. The BodyTom Elite is moved into position and the scan is performed.

The universal transfer board is always used with mattress stiffeners.

The mattress stiffeners provide a solid surface at the head of the bed to prevent the mattress from sagging when a scan is performed. There are usually four mattress stiffeners stored with the BodyTom Elite for easy transport.

The universal transfer board is used for adults.

Note The universal transfer board is an optional accessory that does not come with the system.

Table 58: Universal Transfer Board weight-bearing restrictions

The weight limit of the Universal Transfer board is equal to the weight limit of the patient bed. The weight limit on the portion of the Universal Transfer board that supports the patients head is 7.5 kg / 17lbs. The universal transfer board is used to support and scan the patients head, *only*.

See also "Parts that potentially come into contact with the patient" on page 88.



WARNING The weight limit for the superior portion of the Universal Transfer Board is 7.5kg or 17lbs.




Figure 395: Universal transfer board and stiffeners




Figure 396: Four types of mattress stiffeners

Note: Tipping of the board is a major concern. The universal transfer board **MUST** be securely fastened to the surface prior to placing the patient on the board.

1. Obey all warning labels when using the scan board.
2. Select the appropriate mattress stiffener for the mattress size and insert.



Figure 397: Mattress stiffener in place

3. The universal transfer board requires mattress stiffeners that provide a solid surface at the head of the bed to prevent the mattress from sagging with the weight of the patient when a scan is performed.
4. With the proper mattress stiffener properly inserted, apply the universal transfer board on top.
5. Position the board in accordance with the yellow, safety-warning stickers to avoid a tipping hazard. Do not extend the board beyond the mattress for proper placement.



Figure 398: Universal transfer board properly positioned on the bed on a mattress stiffener

6. When the board is properly positioned on the bed, secure it by using the safety strap.
7. The safety strap must be attached to the board, passed completely under the bed, and secured on the other side.



Figure 399: Universal transfer board with safety strap installed

8. When the universal transfer board is securely fastened to the bed, transfer the patient to the board, and secure the upper strap to the patient and the scan board.
9. When the patient is positioned and securely strapped in, position the scanner over the patient.
10. Initiate the scan.

Chapter 14 Cleaning and Storing the System and Workstation Specifications

Be familiar with this section before using the cleaning or storing the system.

Cleaning the scanner and workstation

When the system is between uses, NeuroLogica recommends keeping it clean as described below. This will help remove body fluids to prevent a health risk and damage to internal parts.

Note NeuroLogica recommends a solution of $\geq 99\%$ pure Isopropyl Alcohol (IPA) to sufficiently clean the equipment.



WARNING Do not use flammable or potentially explosive disinfecting sprays, since resultant vapor could ignite, causing personal injury and/or damage to the equipment.



WARNING To prevent short-circuiting or possible electrical shock, do not spray cleaning agents or spill liquid cleaning agents directly onto the machine.



WARNING Always electrically isolate this equipment from the main electrical supply before cleaning and disinfecting it to prevent short-circuiting or possible electrical shock.



CAUTION The unit surfaces may be cleaned with a soft cloth and the recommended solution or a similar mild non-abrasive cleaning solution. General purpose liquid disinfectant may also be used as necessary. Apply the cleaning solution to the cloth, not directly to the unit.



WARNING Before cleaning the workstation (drive system), be sure to disconnect the workstation from the wall outlet (power source). Failure to do so could result in electrical shock and cause severe injury to you and/or damage to electrical components.



CAUTION Do not allow electrical components to become wet. For eye and hand protection, it is important to wear safety glasses and rubber gloves, respectively.



CAUTION Do not spray cleaning agents directly on the system. Spray a clean cloth with the solution and then wipe down the scanner and workstation.

Cleaning the outside of the scanner and workstation

1. Prepare detergent/disinfectant (regulated by EPA as hospital disinfectant) solution according to instructions on label for correct usage.
 - Use a basin or spray bottle (with product label).
 - Use a pump (usually on detergent/disinfectant containers) to dispense the concentrate in the basin or spray bottle, then fill with correct amount of tap water.
 - If using a spray bottle, empty and rinse out after use.

Note The stability of the solution is unknown after 24 hours; therefore, a fresh preparation of cleaning solution **must** be prepared for each day of cleaning.

2. Use general purpose germicidal cleaner on the external covers and rails. Do not use cleaners on the screens.
3. Use swabs moistened with cleaning solution, clean and remove any dust, soil, dried contrast media, or foreign matter; allow all components to air dry.
4. Wipe down and clean the frame of device and allow to air dry; return to its storage area.

Note Wash (at 25°C) with neutral detergents, *only*; softening agents **are not allowed**.

The following recommended products are registered by the EPA as hospital disinfectant; these solutions are quaternary ammonium compounds and are used in environmental sanitation of non-critical surfaces:

- TB Quat™ is a cleaning solution manufactured by ABC Compounding Co.
- Wex-cide™ is a disinfectant manufactured by Wexford Labs, Inc., product number Wexcide128.

Maintenance of the workstation



WARNING Maintenance checks and all service must be performed by service personnel trained by NeuroLogica Corp. See “Contact information” on page 25.

Storing the system

Storing the scanner and workstation

Store the scanner in a dry, well ventilated, climate-controlled area. You can use the key to lock the scanner when not in use. See “Identifying operator control panel buttons” on page 80 to locate the lock.



CAUTION When the scanner is not in use and stored, it must be plugged in a 120V or 250V outlet to charge the batteries.



CAUTION When the workstation is not in use and stored, it must be plugged in a 120V (or other compatible) outlet to charge the batteries.

Store the scanner on its centipedes or castor wheels (feet).

Note If the floor surface is soft (spongy) store the system on its centipedes to disperse the weight of the system evenly.

1. After transporting the scanner to an acceptable storage location, you can either store the system in **Transport** mode (on its caster wheels) or **Scan** mode (on its centipedes).



Figure 400: BodyTom Elite castor wheels

2. Turn off the scanner and workstation.
The system is now ready to be stored.

See “Powering the workstation” on page 104.

Note It is recommended practice that the scanner is plugged in and turned on even when it is not in use.

Storing the QA phantom

Store the phantom in a secure location with easy access for the daily QA procedure.

Workstation specifications














Table 59: Workstation specifications






Phase	Single	
Voltage Range	100-240VAC \pm 10%	
Factory Outlet Recommendations	NEMA 5-20R	NEMA 6-30R
Frequency	50 or 60Hz	
Battery Capacity	Fully Charged/ 12 hrs. (Typical)	
Typical Usage	110-120 VAC 60 Hz	230-240VAC 50 Hz
Wiring	125V, 2 Pole, 3 Wire Grounding	250V, 2 Pole, 3 Wire Grounding
Battery Operating Voltage	51.8VDC	
Overall width	41in. (104cm)	
Overall height	79in. (199cm)	
Overall length	101in. (256.5cm)	
Weight EST	3510 lbs. (1592kg)	
Battery power (2) 12 VDC (lithium polymer)	800W	
Max programmed speed fwd.	1.6 MPH	
Max recharge time	~ 8hrs.	
Max continuous operation	8hrs.	
Locking and unlocking cycles	20	
Hrs transport over floors	2hrs.	
Hrs system locked no external power	2hrs.	
Max slope holding angle with scanner	7° C (44.6° F)	
Max doorway threshold	1in. (2.54cm)	

Max elevator threshold	¾in. (1.905cm)
Height to locking adapter	8.59in. - 8.69in.
Min/Max storage temperature	-25° C to 70° C (-13° F to 158° F)
Min/Max operating relative humidity	
Min/Max storage relative humidity	20% to 85% (non-condensing)
Min/Max ambient operating temperature	15° C to 35° C (59° F/95° F)

Understanding the symbols and product-marking plate

Table 60: Symbols and product-marking plate

Symbol	Description
	Indicates alternating current.
	Indicates protective earth (ground).
	Indicates a caution: consult accompanying documents.
	Indicates a caution: risk of electrical shock.
	Indicates electrostatic sensitive devices.
	Indicates a warning: high temperature.
	Indicates temperature limits.
	Indicates mechanical deactivation device.
	Indicates a radiation precaution; may be affected by radiation from other sources; may produce interference that affects other equipment.
	Indicates manufactured by.
	Indicates toe/foot crush hazard; positioned in the front of the unit adjacent to the safety bumper and is a warning for population when lowering the unit.
	Indicates a coil power cord.
	Indicates a chain hazard could cause severe personal injury.

Symbol	Description
	Indicates to keep away from rain for packaging.
	Indicates a humidity limit for packaging.
	Indicates a warning: battery charging.
	Indicates fuse usage.
	Indicates accompanying operating instructions in the user manual must be followed to safely operate equipment.

Note Disregarding information on safety is considered ***abnormal use***.

Locating the product-marking plate on the workstation




Figure 401: Product-marking plate on side of the workstation

Listing of replacement parts for workstation

To ensure proper compliance requirements of replacement parts, (for example, cables and accessories), parts must be purchased through NeuroLogica Corp.



WARNING Using other manufacturer cables and accessories may affect EMC performance. Unauthorized use of these items will void warranty and may cause harm to patient, others and/or equipment.

Product Safety and Electromagnetic Comparability

Tested by: Intertek Testing Services NA, Inc., 70 Codman Hill Road,
Boxborough, MA 01719

Appendix A Glossary

A

Algorithm	Mathematical filter applied to raw data during CT image reconstruction to remove blurring artifact inherent to back-projection. Also referred to as a kernel.
Annotation	User comments or text added to an image.
Anterior	Front of the patient's body
Application Entity (AE)	An end point of a DICOM information exchange, including the DICOM network or media interface software; that is, the software that sends or receives DICOM information objects or messages. A single device can have multiple AEs.
Attenuation	The reduction in intensity of a radiation beam as it passes through a substance.
Automatic Exposure Control (AEC)	Software used to adjust or modulate the mA throughout an acquisition to reduce patient radiation dose to a minimum.
Axial scan mode	Data acquisition while the scanner remains stationary. The scanner position may be incremented between exposures to collect data over a longer Z axis range. Also referred to as step and shoot.

B

Bolus Tracking	Monitors flow of contrast media in vessel and triggers scan at optimal timing. This is a scanner feature to automatically initiate a prescribed Axial, Helical or Dynamic scan when a threshold level of contrast enhancement is reached at a specified region of interest.
----------------	---

C

Collimation	Restricts x-ray to only the selected anatomy, minimizing dose to patient and reducing scatter.
-------------	--

Computed Tomography Angiography (CTA)	A test that uses x-rays to provide detailed pictures of the heart and the blood vessels that go to the heart, lung, brain, kidneys, head, neck, legs, and arms. A CT angiogram can show narrowed or blocked areas of a blood vessel.
Computed Tomography Dose Index (CTDI)	An approximate measure of the radiation dose received in a single CT section or slice.
Computed Tomography Dose Index Volume (CTDI _{vol})	Represents the dose for a specific scan protocol, which considers gaps and overlaps between the radiation-dose profiles from consecutive rotations of the x-ray source. It is not patient dose. The CT dose index volume is noted as CTDI _{vol} . It is reported in units of mGy. The CTDI _{vol} is based on measurements made by the manufacturer in a factory setting. The CTDI _{vol} is calculated differently for both the Axial and the Helical mode: For Axial scan mode: $CTDI_{vol} = [(N \times T)/I] \times CTDI_w$ For Helical scan mode: $CTDI_{vol} = 1/pitch \times CTDI_w$ See also dose and patient dose.
Computed Tomography Dose Index (CTDI _w) weighted average	The measure of ionizing radiation exposure per slice of data acquisition. CTDI represents the integrated dose along the Z axis from one axial CT scan (one rotation of the x-ray tube). The CT Dose Index is noted as CTDI _w .
Computed Tomography (CT) number	Relative value assigned to each pixel to quantify the attenuation occurring in each voxel in comparison with the attenuation of water. The calculated CT number for a given pixel is given in Hounsfield units (HU).
Computed Tomography Perfusion (CTP)	Evaluates cerebral perfusion or level of blood flow in the brain by monitoring the initial passing of iodinated contrast media through the vasculature of the brain.
Contrast media	Used to improve sensitivity and specificity of clinical diagnoses.

Contrast resolution	The ability of a CT system to detect an object with a small difference in linear attenuation coefficient from the surrounding tissue. Also referred to as low-contrast detectability or sensitivity.
---------------------	--

D

Digital Imaging Communication in Medicine (DICOM)	Digital Imaging and Communications in Medicine, or DICOM, is a standard that helps people doing work in the field of radiology. The DICOM standard is designed to promote communication and integration between a variety of radiology imaging systems and equipment used in filmless radiology.
Digital tilt	The ability to correct the image post acquisition and correct positional inaccuracies prior to sending to PACS.
Dose	The generic term that refers to the CTDI _{vol} , the standardized parameter to measure scanner radiation output – or the amount of amount of ionizing radiation absorbed by patient per unit mass.
Dose Length Product (DLP)	The measurement of dose for an entire series of CT images. DLP is equal to the calculated dose per section multiplied by the length of a CT acquisition along the Z axis.
Dynamic Host Control Protocol (DHCP)	A standardized network protocol used on Internet Protocol (IP) networks. The DHCP is controlled by a DHCP server that dynamically distributes network configuration parameters, such as IP addresses, for interfaces and services.
Dynamic scan mode (multiple detector widths)	Data acquisition at multiple time points over the same anatomic location(s).

E

Electromagnetic Compatibility (EMC)	The branch of electrical sciences that studies the unintentional generation, propagation, and reception of electromagnetic energy with reference to the unwanted effects (Electromagnetic interference (EMI)) that such energy may induce.
Electromagnetic Interference (EMI)	A disturbance generated by an external source that affects an electrical circuit by electromagnetic induction, electrostatic coupling, or conduction. The disturbance may degrade the performance of the circuit or even stop it from functioning.

F

Field of View (FOV)	The diameter of the acquired attenuation data displayed across the image matrix.
---------------------	--

H

Helical scan mode	A CT acquisition whereby an x-ray acquisition whereby the x-ray tube and scanner move continuously during scanning, yielding a data set in the form of a helix. Also referred to as spiral.
Hospital Information System/Radiology Information Systems (HIS/RIS)	A Radiology Information System (RIS) is the core system for the electronic management of imaging departments. The major functions of the RIS can include patient scheduling, resource management, examination performance tracking, examination interpretation, results distribution, and procedure billing. RIS complements Hospital information systems (HIS) and Picture Archiving and Communication System (PACS), and is critical to efficient workflow to radiology practices.
Hounsfield Unit (HU)	The unit of the CT number scale assigned to each pixel to quantify relative attenuation.

I

Interscan delay time	Minimum amount of time that must transpire between end of one scan and initiation of next scan. Interscan delay times include idle time between scans to allow tube cooling.
Iterative Bone Correction (IBC)	A feature build into the reconstruction software, which performs a correction on every single Axial image the scanner produces, including both primary series from a scan as well as secondary reconstruction images. Current IBC settings were chosen to provide optimal correction for standard medical imaging; however, the setting can be customized as needed.

K

Kernel	A mathematical filter applied to raw data during CT image reconstruction to remove blurring artifact inherent to back-projection. Also referred to as an algorithm.
--------	---

M

mAs	Tube current-time product: The product of tube current and exposure time per rotation, expressed in units of milliamperere seconds (mAs).
Matrix	Two-dimensional (2D) grid numbers arranged in rows and columns.
Maximum Intensity Projection (MIP)	The multiplanar reformation technique that displays only the maximum pixel value along a ray traced through the object to the viewers assumed perspective in front of the scanner display screen.
Mean Transit Time (MTT)	A common measurement during CTP studies of the brain. Refers to the average transit time, in seconds, needed for blood to pass through a given region of brain tissue.

milli amperage (mA)	Tube current: the number of electrons accelerated across an x-ray tube per unit time, expressed in units of milliampere (mA).
Modality Performed Procedure Step (MPPS)	A mechanism for modalities to pass information about the imaging performed back to the HIS/RIS or PACS.
Modality worklist manager	Scheduled (but not yet scanned) patient list.
motion artifact	Voluntary and involuntary patient motion during CT scan, appearing as a streak artifact on image; ghosting or blurring of image.
Multi-Planar Reformation (MPR)	The process of displaying CT images in a different orientation from the one used in the original reconstruction. Allows for reformation of images in planes that would otherwise be difficult or impossible to acquire with CT. Requires only image data. Raw data is not utilized.

N

Noise	Random statistical variations in the signal. Can be quantum noise, electronic noise due to lost signal, or artifact noise. Manifests itself as overall graininess of the reconstructed image.
-------	---

P

Partial volume artifact	Occurs when an object is only partly positioned within a voxel or is much smaller than the overall voxel volume. The object's attenuation is not accurately represented by the pixel value. Overlapping reconstructions further reduce partial volume artifacts.
Patient Browser, local database	Where the already-scanned patient list is stored.
Patient coordinates	References are as follows: <ul style="list-style-type: none"> • X left to right. • Y anterior to posterior. • Z head to feet.
Patient dose	The absorbed dose to a patient. See also CTDI _{vol} .

Peak kiloVoltage (kV)	The penetrating power of the photons coming from the x-ray tube.
Picture Archive and Communications Systems (PACS)	Stores medical information, including 2D images, and 3D medical images. All modern PACS setups will work with DICOM.
Pitch	In Helical mode, refers to the speed of the scanner movement over the table as the scanner rotates.
Pixel	A single, picture element of image matrix.
Post reconstruction	Prescribing the reconstruction parameters after scan acquisition.
Projection	View of anatomical cross-section from a particular vantage point.
Prone	Patient lying on stomach.
Protocol	Prescribes the acquisition and reconstruction parameters to be used for a scan.

Q

Quality Assurance (QA)	Procedure of performing periodic specified tests or measurements to assure that a set quality level, as specified by system manufacturer, has not been compromised.
------------------------	---

R

Radiation Safety Officer (RSO)	The person within an organization responsible for the safe use of radiation and radioactive materials as well as regulatory compliance.
Radio Frequency Interference (RFI)	Also called Electromagnetic Interference (EMI), is an unwanted disturbance that affects an electrical circuit due to electromagnetic radiation emitted from an external source.
Raw data	A transmission measurement obtained by the detectors used to mathematically reconstruct the CT image.
Reconstruction filter	Used to ensure accurate anatomical reconstruction. Also allows for either spatial resolution or low-contrast-resolution enhancement.

Region Of Interest (ROI)	Provides a quantitative analysis of the Hounsfield values of a specific anatomic area. A graphic outline in the shape of a circle is placed over an area on the image. Software calculates the average CT number in HU within the ROI.
Resolution	A scan time, per slice, in Axial mode, only.
Retrospective reconstruction	Reconstruction performed after the initial prospective reconstruction. Multiple retrospective reconstructions of raw data are possible, with changes to display FOV, kernel, slice thickness etc.

S

Scan delay	The time between the initiation of contrast agent administration and CT data acquisition. The chosen scan delay determines the phase of contrast enhancement for a given CT acquisition.
Scan protocol	A list of scanner-load parameters used to perform an x-ray exposure.
Scan types	Axial, Helical, Dynamic, Reference, and Scout.
Scout	Digital survey radiograph acquired by the CT system for the purpose of prescribing the cross-sectional acquisition. Like a conventional radiograph, the scout is produced by translating the scanner over the patient without tube or detector rotation. Also referred to as topogram or scanogram.
Series	A set of images acquired in a scan.
Size Specific Dose Estimate (SSDE)	Not dose to any specific organ but rather the mean dose in the center of the scanned volume. That is, SSDE is not the exact patient dose, as factors such as scan length and patient composition can differ from the assumptions used to calculate SSDE, for example conversion factors based on patient size provided to estimated patient dose for a patient of a particular size.
Slice spacing (Spacing)	The distance between the center of one CT slice and the center of the next slice.

Slice thickness	The dimension of a constructed CT slice along the longitudinal direction of acquisition (Z axis).
Spatial resolution	The ability of a CT imaging system to display fine details, separately. Given in units of line pairs per centimeter (lp/cm).
Supine	Lying on back.

T

Temporal resolution	The ability of a CT system to freeze motion and provide an image – free of blurring.
Test Bolus	Scan mode used to measure the contrast transit time using a small injection of contrast media.
Threshold	The CT number (Hounsfield Unit (HU)) where Bolus Tracking tool will trigger the system to begin the scan.
Time Attenuation Curve (TAC)	A graph of the contrast enhancement versus time. TAC is used to determine blood flow rate in seconds for contrast timing.
Time delay	Monitoring delay: Time from injection to the start of monitoring scans.
Transverse plane	Perpendicular to direction of Z axis.

V

Volume Rendering (VR) image or object	A 3D modeling technique that utilizes the entire acquired dataset but adjusts the opacity of the voxels included in the 3D image according to their tissue characteristics.
Voxel	Abbreviation of volume element. Refers to the volume of tissue represented by a pixel in the matrix used to display the CT image.

W







Window Level (WL)	The pixel value given in Hounsfield Units (HU) at the center of the window width. Window Level controls the brightness (density) of the CT image.
-------------------	---



Window Width (WW)	The range of pixel values assigned a shade of gray in the displayed CT image. Window width controls the contrast of the CT image.
-------------------	---


Appendix B Listing of All Buttons, Tools, and Icons

Status bar icons

Table 61: Status bar icons












Status bar icon	Status bar icon name	Status description
	X-ray status	Identifies x-ray as on or off. The icon changes from a gray/black icon when x-ray is off to a rotating yellow/black icon when x-ray is on.
	System state	Identifies the system's current state. The orb changes color depending on the state the system is in. See Table 25 on page 118 for a list of the different orb colors and system states they identify.
	Scanner position	Identifies the system's current position relative to its zero reference.
	System E-STOP status	Identifies when E-STOP is engaged. The icon will flash when E-STOP is pressed.
	System tube heat status	Indicates the current X-Ray tube heat status. The values are color coded as follows: Blue 0% - 25% Yellow 26% - 50% Orange 51% - 75% Red 76% - 100%
	Scanner battery capacity status	Indicates the remaining scanner battery percentage available. The capacity values are color coded as follows: Green 100% - 51% Yellow 50% - 25% Red 24% - 0%



Status bar icon	Status bar icon name	Status description						
	System air freshness status	<p>Indicates the air freshness status; it is recommended that an air calibration be performed:</p> <ul style="list-style-type: none">• Every eight (8) hours.• When the air freshness status falls below 50%.• If the scanner is moved to an area with a dramatic change in humidity and/or temperature. <p>The calibration status values are color coded as follows:</p> <table><tr><td>Green</td><td>100% - 51%</td></tr><tr><td>Yellow</td><td>50% - 25%</td></tr><tr><td>Orange</td><td>24% - 0%</td></tr></table> <p>After calibration it returns to 100%.</p>	Green	100% - 51%	Yellow	50% - 25%	Orange	24% - 0%
Green	100% - 51%							
Yellow	50% - 25%							
Orange	24% - 0%							
	Workstation battery capacity status	<p>Indicates the remaining workstation battery capacity available. The capacity values are color coded as follows:</p> <table><tr><td>Green</td><td>100% - 21%</td></tr><tr><td>Yellow</td><td>20% - 11%</td></tr><tr><td>Red</td><td>10% - 0%</td></tr></table> <p>You will be prompted to plug the workstation into an outlet to charge if the battery capacity is low; a scan cannot complete when the battery capacity is 10% or lower.</p> <p>When the workstation reaches the red capacity range, the system will shut down. A message informs you that the system will shut down due to low battery.</p> <p>The lightning bolt icon signifies that the workstation is currently charging and goes away when unplugged.</p>	Green	100% - 21%	Yellow	20% - 11%	Red	10% - 0%
Green	100% - 21%							
Yellow	20% - 11%							
Red	10% - 0%							

Status bar icon	Status bar icon name	Status description						
	Image storage space status	<p>Indicates the available disk space for image storage. The available space values are color coded as follows:</p> <table><tr><td>Green</td><td>100% - 51%</td></tr><tr><td>Yellow</td><td>50% - 20%</td></tr><tr><td>Red</td><td>19% - 0%</td></tr></table>	Green	100% - 51%	Yellow	50% - 20%	Red	19% - 0%
Green	100% - 51%							
Yellow	50% - 20%							
Red	19% - 0%							

System state orbs









Table 62: System state orbs












Orb	Color	State
	Dark gray	The system is in an unknown state.
	Light gray	The system is powering up or down.
	Dark purple	The system is busy.
	Purple	The system is completing air calibration.
	Light purple	The system is archiving.
	Blue	The system is idle.
	Green	The system is ready to perform a scan.
	Light yellow	The system is planning.
	Dark yellow	The system is preparing.
	Light orange	The system is reconstructing.
	Dark orange	The system is scanning.














Orb	Color	State
	Pink	The system is not ready.
	Red	The system is in fault.













Workstation buttons














Table 63: BodyTom Elite workstation buttons


Workstation button	Action
	In Protocol Manager – adds a new protocol from the list.
	In Patient Browser – selects the archive destination for selected information.
	In Patient Browser – used to begin a protocol.
	In Protocol Manager – used to create a new protocol from a previously saved protocol.
	In Patient Browser – generates the dose for the selected patient.
	In Patient Registration – cancels the current query. In Patient Browser – cancels any series being imported.
	In Acquisition – cancels the entire exam being performed.
	In Acquisition – cancels the current scan within a protocol.
	In System Configuration – clears information in fields.

Workstation button	Action
	In Protocol Manager – closes popup. In Store/Print Queue – closes the Store/Print Queue popup. In System Configuration – closes System or User Configuration dialog boxes.
	In Patient Browser – allows you to select multiple series of patient images to compare in Viewing .
	In Acquisition – authorizes the scanner to move to the next step.
	In Patient Registration – deletes patient(s) from the Stored Results list. In Patient Browser – deletes selected exam information from Patient Browser . In Store/Print Queue – deletes a series to be stored or a series that failed to store. In Protocol Manager – deletes a saved protocol. In System Configuration – clears a saved dose setting to remove the restriction.
	In System Configuration – deletes saved dose settings to remove all restrictions.
	In Store/Print Queue – when you select one or more series, displays an explanation of why a series failed to store.
	Move selected item down the list.
	In Patient Registration – used to modify protocols.
	In System Configuration under DICOM Servers – echoes the selected server to test the connection.
	In Protocol Manager – exports protocols to a media device.
	In Acquisition – provides extension of current protocol.

Workstation button	Action
	In Acquisition – completes the examination. Completes all protocols, builds Dose SR and images, and directs user to Patient Browser .
	In Viewing - generates a new series with the Series Description field information – based on the selected MPR .
	In Patient Browser – imports the exam information from PACS or Media . In Protocol Manager – imports previously exported protocols to the workstation.
	In Patient Registration – manually enters a new patient and, when completed, takes you to the Acquisition tab to acquire the data for a scan.
	In Patient Browser – used to select the destination for patient data to media.
	In Patient Browser – combines two different image sets.
	In Patient Browser – minimizes the Import for Media popup.
	In Protocol Manager – used to create a new protocol.
	To accept selections you make.
	In Patient Browser – used to select the destination for patient data to PACS .
	In Patient Browser – used to select only patient information for a patient when using Register feature.
	In Patient Browser – used to select all exam information for a patient when using Register feature.
	In Store/Print Queue – when you select one or more series, temporarily stops the series from being stored. This is a toggle button with the Resume button.










Workstation button	Action
	In Acquisition – pauses entire protocol.
	In Acquisition – pauses current scan within a protocol.
	In System Configuration under Audio Configuration – used to play audio files.
	For Quality Assurance – used to prepare workstation to run a Quality Assurance test.
	In Acquisition – selects an existing protocol for the current study.
	In Patient Registration – searches the HIS/RIS server for scheduled patients.
	In Patient Registration – registers the selected patient and then takes you to the Acquisition tab to acquire the data for an examination (scan). In Patient Browser – opens the Create New Study dialog box and moves selected patient to Acquisition .
	In Patient Registration – used to take patient information stored in HIS/RIS server to Patient Registration to choose the patient and study.
	In Acquisition – repeats the last scan that was performed.
	In Post Reconstruction , sends the last acquired scan from the recon workstation to the Patient Browser.
	In Viewing – reverts all images back to original state. In System Configuration – resets information back to default or clears information in fields.
	In Store/Print Queue – when you select one or more series, continues to store previously paused series. This is a toggle button with the Pause button.









Workstation button	Action
	In Store/Print Queue – when you select one of more series, tries to archive the selections.
	In System Configuration – saves updated information.
	In Patient Registration – searches queried patient entries for specific information.
	In Patient Browser – shows patient, study, series, and image information; used to modify series scanned under a wrong patient.
	For Daily Calibration – begins the daily (air) calibration.
	In Acquisition – begins any post-reconstructions that were defined during the protocol setup. In Post Recons – begins a manual reconstruction
	In System Configuration under Audio Configuration – stops audio files from playing.
	In Patient Registration – selects patient(s) from query results and moves them into the Stored Results list.
	In System Configuration under User Accounts – used by administrators to unlock a user's account.
	Move selected item up the list.
	In Protocol Manager – updates information on an existing protocol. In System Configuration – updates information.
	In Acquisition – prompts the application to send the selected protocol to the scanner and verify that the scanner has tube and battery capacity to perform the protocol.
	In Patient Registration – shows selected patient details. In System Configuration – shows information.













Workstation button	Action
	In Viewing – to load and views images.

Viewing tools

Table 64: Viewing tools

Tool	Tool name	Action
	Angle	2D, CTP, and Viewing tool – draws an angle on the image and displays the angle information.
	Arrow	2D, CTP, and Viewing tool – draws an arrow on the image, which can be repositioned.
	Calculate CBF, CBV, MTT Map	CTP only tool – calculates the Cerebral Blood Flow (CBF), Cerebral Blood Volume (CBV) and Mean Transit Time (MTT) maps.
	Clear Perfusion Map	CTP only tool – cancels the calculations and returns to Calculation mode.
	Capture	Common tool - saves a screen capture of selected viewport.
	Capture all Viewports	Common tool - saves screen captures of all visible viewports.
	Cine Backward	2D only tool – cines backward through the image.
	Cine Forward	2D only tool – cines forward through the image.
	Clear Active	Common tool (Acquisition, Post Reconstruction, Viewing) - resets the tool to the default pointer device.

Tool	Tool name	Action
	Change the Window Width Level	Common tool (Acquisition, Post Reconstruction, Viewing) – click and move pointer over image. Left click and hold down the mouse button and drag in chosen direction to adjust image width and level. Width and level values appear in the Width/Level status display. A pre-defined width level setting can also be selected. Select the preset from the dropdown list below the WL Preset button. Width and level presets can also be saved or deleted.
	Field Of View	Post Reconstruction tool – adjusts the Field Of View (FOV) prior to reconstruction.
	Flip Horizontal	2D only and Viewing tool – flips images right or left.
	Flip Vertical	2D only and Viewing tool – flips images up or down.
	Invert	Common tool (Viewing) - inverts black to white and white to black.
	Load	Viewing tool – loads images from selected series into viewing.
	Measure (Line)	2D, CTP, and Viewing tool – draws a line on the image and displays length information.
	Pan	Common tool (Acquisition, Post Reconstruction, Viewing) - click and move pointer over image. Left click and hold down the mouse button and drag the image in the chosen direction. Release mouse button to position image in new location.

Tool	Tool name	Action
	Peak Image	CTP only tool – displays the image that has the most visible contrast (based on arterial ROI placement).
	Perfusion Artery/Vein Selection	CTP only tool – selects the artery and vein to be used for performing perfusion calculations.
	Redo	Axial and Helical Viewing tool – restores the last text editing or resizing and positioning of controls – if no other action occurred since last time the Undo button was clicked.
	Redo	Viewing tool – performs the most recent action, again. The button is disabled if the application cannot redo the application.
	Region of Interest (ROI)	2D, CTP, and Viewing tool – defines a circular ROI and displays the ROI information (5mm diameter by default).
	Reset	Post Reconstruction and Acquisition tool – resets the display to default viewer settings.
	Reverse Image Stack	2D only and Viewing tool – reverses the order in which images display.
	Rotate	MPR only and Viewing tool – rotates the image.
	Scan Region Re-Draw	Acquisition tool – if scout lines and the scan region is deactivated, allows you to reactivate.
	Show Artery/Vein Flow Graph	CTP only tool – displays the Arterial Venous Flow graph. 
	Stop	Post Reconstruction tool – cancels the current, post-reconstruction request. All images are generated until you click the Stop button.






















Tool	Tool name	Action
	Stop Cine	2D only tool – stops the cine forward and backward.
	Text (Annotation)	2D only and Viewing tool – creates text box for annotation.
	Tilt	MPR only tool – corrects a rotated image.
	Toggle Layout	Acquisition tool – changes the layout to 2x2. Repeat process to return to 1x1.
	Toggle Scouts	Acquisition tool – removes scouts from Acquisition.
	Undo	Axial and Helical Viewing tool – reverses the most recent action taken (a successful copy, cut, delete, undo or paste action).
	Undo	Viewing tool – removes the most recent action performed on image. The workstation remembers the last five adjustments made. The tool is disabled if the workstation cannot redo the adjustment.
	Windowing	Common tool - adjusts the width and level of the selected image.
	Zoom	Common tool (Acquisition, Post Reconstruction, Viewing) - click and move the pointer over the image. Left-click the mouse and hold down the left-mouse button and move in upward direction to zoom in (enlarge) and downward to zoom out (shrink).

Table 65: Pendant buttons

Pendant	Button	Description	Action
		POWER	Illuminates when power is supplied to pendant.
		LASER	Turns on all three positional lasers. While the lasers are on, the scanner spins for the internal laser to be seen within the scanner opening.
		GO TO SCAN PLANE	Moves the scanner forward approximately 30cm. This is the distance between the internal and external lasers.
		ZERO REFERENCE	Sets the scanner to zero before starting a scout or a scan.
		MOVE BACKWARD (slow)	Pressing and holding moves the scanner backward 10mm per second.
		MOVE FORWARD (slow)	Pressing and holding moves the scanner forward 10mm per second.
		MOVE BACKWARD (fast)	Pressing and holding moves the scanner backward 60mm per second.
		MOVE FORWARD (fast)	Pressing and holding moves the scanner forward 60mm per second.
		SET MEMORY	Allows the user to program Scan and Rest positions for the scanner.

		SCAN POSITION	Moves the scanner to the Scan Position saved using the Set Memory feature.
		REST POSITION	Moves the scanner to the Rest Position saved using the Set Memory feature.

Appendix C Sample of Reference Protocols Provided

Table 66: Sample of BodyTom Elite adult protocols and important estimates

Protocol Name	Type	kV	mA	Slice Thickness /Spacing	Sharpness	Resolution	Coverage	CTDI _{vol} (mGy)	DLP ⁵ (mGy.cm)
Adult Head Axial	Axial	120	200	5.0 x 5.0	Soft Tissue	1 Sec.	250	47.32	1183
C-Spine Helical	Helical	120	250	1.25 x 1.25	Soft Tissue - Abdomen	Pitch = 0.8	400	24.16	966.4
Adult Chest Helical	Helical	120	150	1.25 x 1.25	Bone	Pitch =0.8	450	14.49	652.05
Adult Abdomen Helical	Helical	120	250	2.5 x 2.5	Soft Tissue - Abdomen	Pitch = 0.8	500	24.16	1208

Table 67: Sample of BodyTom Elite pediatric protocols and important estimates

Protocol Name	Type	kV	mA	Slice Thickness /Spacing	Sharpness	Resolution	Coverage	CTDI _{vol} (mGy)	DLP ⁶ (mGy.cm)
Pediatric Head Axial	Axial	100	175	5.0 x 5.0	Soft Tissue	1 Sec.	200	32.02	640.4

⁵ DLP is based on length from coverage column⁶ DLP is based on length from coverage column

Appendix D Automatic Exposure Control

1 Introduction:

Automatic Exposure Control (AEC) is a feature which allows the exposure to automatically be modified based on the attenuation of the scanned object. The main objective of AEC is to optimize the x-ray current based on prior knowledge of the scanned objects profile. AEC is used to optimize patient exposure while attempting to maintain acceptable diagnostic quality of the reconstructed images.

AEC uses image noise to optimize the scan current. The image noise on CT scanners can be traced to two sources: **Electronic Noise** and **Quantum Noise**. Electronic Noise is generated by the electronic components of the Data Acquisition System (DAS). **Quantum Noise** is related to x-ray generation. Currently **Quantum Noise** is the major component of noise on CT images, the contribution of **Electronic Noise** has become less significant since the early days of CT scanners.

2 Image Noise:

2.1 Electronic Noise:

The DAS is composed of the crystals, the photodiodes, the Analog to Digital Convertors (ADC) and other electronic components known as “converter cards” since they convert x-rays into a quantifiable current. Thermal Noise is the most common source of electronic noise in the CT system. As the scanners internal temperature increases thermal noise becomes the dominant component of the DAS’s electronic noise. Imperfections in the semiconductor chips used in the DAS also contribute to the Electronic Noise. However, with the advance in semiconductor crystals this has become less relevant. Currently Electronic Noise has no significant impact on image quality when using proper scan parameters, i.e., scan voltage, current and exposure.

2.2 Quantum Noise:

The generation of x-ray photons can be described by a Poisson random process. Poisson random processes are used to describe event generation over a fixed time interval. A Poisson random process is used to describe Queues in general. In a queue the number of new arrivals to the queue over a fixed time interval follow a Poisson distribution.

The **Quantum Noise** is related to the standard deviation of the Poisson distribution which inversely proportional to the square roots of the number of events:

$$\sigma_I \propto \frac{1}{\sqrt{N}}$$

Where σ is the image noise and N is the number of detected photons. The above equation can help relate the image noise to the scan current, I , the scan time, t , and the slice thickness, Sw , since the number of photons is proportional to either one:

$$N \propto It \text{ and } N \propto Sw$$

The image noise of a given scanner can be written as a function of the scan parameters:

$$\sigma = \frac{K}{\sqrt{I \cdot t \cdot Sw}}$$

To reduce the image noise, we can either increase the scan current, the scan time or the slice width, or any combination thereof. K is constant that is dependent on the image reconstruction process. It follows, for two different scan currents I_1 and I_2 over the scan times t_1 and t_2 and using the same slice thickness, the image noises for the two scans can be related using the following equation:

$$\frac{\sigma_2}{\sigma_1} = \sqrt{\frac{I_1 \cdot t_1}{I_2 \cdot t_2}} \text{ or } \sigma_2 = \sqrt{\frac{I_1 \cdot t_1}{I_2 \cdot t_2}} \sigma_1$$

The same relation exists between the image noise and the slice thickness

$$\frac{\sigma_2}{\sigma_1} = \sqrt{\frac{Sw_1}{Sw_2}} \text{ or } \sigma_2 = \sqrt{\frac{Sw_1}{Sw_2}} \sigma_1$$

As a result, doubling the slice thickness can reduce the image noise by almost 40%. Figure 1 shows the Noise as a function of the scan current (left) and the inverse of the scan current (right). The scan done for the same slice thickness of 2.5 mm and same scan time of 1 second.




Figure 1: The Noise in a 20 cm water phantom as a function of the scan current I (left) and the inverse of square root of the current $1.0/\sqrt{I}$.

3 AEC working instructions:

3.1 AEC input parameters:

AEC requires three different input parameters: Minimum mA, Maximum mA, and the Noise Level.

Minimum mA: The Minimum mA is set to prevent an unacceptable amount of noise in the reconstructed images.

Maximum mA: The Maximum mA value is used to prevent over-exposure of the patient being scanned. This is typically set to reduce the chances of over exposure which may lead to radiation sickness. However, in the case of the BodyTom Elite, the maximum scan current is set to protect the x-ray tube.

The noise level: The first step of AEC is to build an equivalent stack of cylindrical water phantoms or Water Equivalent Diameter (WED). The WED is created based on the measured attenuation from the scout. The noise level represents the noise in each section of the stack of water equivalent phantoms. The selected noise level should be within an acceptable range, and it should be dictated by the scan protocol.

One of the key features of the scanners AEC is the ability the user has to re-adjust the noise level based on the estimated mA until the desired mA profile is attained. The feature works as follows:

1. A protocol is created with the AEC feature enabled, and a Minimum mA, Maximum mA, and Noise Level are defined.
2. The above-mentioned protocol is used to create a Scout

image.

3. The user can toggle the AEC graph to view the scan current profile. The profile will be overlayed on top of the scout.
4. If the scan current profile is acceptable, then the user will initiate the diagnostic scan as desired.
5. If the scan profile is not acceptable, the user can **adjust** the noise level and **recalculate** the scan current.
6. The process can be repeated as many times as needed until the user is satisfied with the current profile. The AEC tool will allow the user to view the scan current before initiating the actual scan.

3.2 The scan parameters:

The scan protocol parameters are not needed for AEC however they do affect the current estimation:

1. **kV**: the scan kV is used to select the appropriate noise table used for estimating the scan mA.
2. **Slice Thickness**: The noise is measured at a slice thickness of 5.0 mm; however, the scan protocol slice thickness can be any of the allowable thickness values. The selected slice thickness is then used to adjust the noise table using the equation in section 2.2. The entire noise table will be multiplied by the square root of the slice thickness ratio. The multiplication factor is:

$$\alpha = \sqrt{\frac{5.0}{\text{Selected Sw in mm}}}$$

3. **The reconstruction kernel**: The noise image depends on the reconstruction kernel. AEC is limited to SoftTissue and PostFossa Kernels. AEC will be disabled if the user selects a different reconstruction kernel.

3.3 Notes

When AEC is selected the user should be aware of:

1. **Patient Positioning**: The patient should be properly positioned as close as possible to the scanner iso-center. Failure to do so can lead to an over-estimate of the scan current leading to an increase in patient dose.

2. **Presence of metal implants:** AEC should not be used if the patient has metal implants in the region to be scanned.
3. **The measured noise:** the final noise in the image depends on the size of the scanned patient. AEC assumes that the patient is cylindrical, as such the measured noise level could be different then the selected noise level.
4. **Anatomical features:** AEC should be used when the region to be scanned includes significant differences in attenuation, such as the chest and abdomen. Anatomical regions with slight differences in attenuation like the head, will not benefit from AEC use.

3.4 Sample protocols:

Below are some suggested protocols. The noise levels depend on the size and weight of the patient. The noise levels in the table below are for illustration purposes. The site physicist and CT manager should dictate the final noise levels.

	Noise Level	Minimum mA	Maximum MA	Slice thickness
Chest scanning	15	50	250	1.25
Chest/Abdomen	13	50	280	2.5

4 AEC algorithm description:

AEC uses the measured attenuation of the scanned object and the selected noise level to estimate the scan current at each planned scan location during the scan. The mA is typically estimated using different water phantom diameters.

The BodyTom Elite uses Z-modulated AEC where each planned scan location is modeled using a cylindrical water phantom, or WED. The WED is calculated using the scout profile. Once the WED is estimated a specific mA value is assigned to each planned scan location. The Flowchart (Figure 2) below describes the basic steps for using AEC:

1. Select the appropriate AEC parameters to be used, those values are:
 - a. Minimum mA
 - b. Maximum mA
 - c. Noise Level

2. Acquire a Scout using the same kV that will be used for the Axial or Helical acquisition.
 - a. An AP or Lateral Scout can be used; however, AP scouts are preferred.
3. For each planned slice, the system calculates the WED and assigns a specific mA value to that location. Figure 3 shows the WED of the equivalent water phantoms as calculated based on the scout image (left). The estimated diameters (right) shows that the water portion of the phantom match the true diameter of the phantom.
4. Using the measured noise in different diameters water phantoms at different mA levels, find the mA that generates the selected noise level.
5. Adjust the mA based on the scan mode. Figure 4 shows the estimated mA as well as the adjusted mA for the axial scan mode.




Figure 2: The AEC flowchart



Figure 3: the AP scout of a 20 cm water phantom and the estimated radius based on the scout.




Figure 4: The estimated current using the equivalent diameter (left) and the scan current for axial scan mode (right).

5 The Noise measurements:


The noise tables used for predicting the scan current is measured for each scan voltage. For each scan voltage the noise is measured using a set of predefined scan currents of 50, 100, 150, 200, 250 and 300 mA. The noise is measured using different water phantoms. In our case we have used 150, 200 and 300-mm water phantoms. Table 68 shows a sample of the noise table at 120 kV.

Table 68: The measured noise at 120 kV

Phantom Diameter (mm)	Scan Currents (mA)					
	50	100	150	200	250	300
100	7.5868	5.4078	4.4009	3.8519	3.4769	3.1679
150	10.76	7.6690	6.2493	5.4768	4.8741	4.4806
200	17.41	12.3721	10.1349	8.8315	7.8699	7.1964
250	28.400	20.2000	16.4000	14.3000	12.7000	11.6000
300	49.952	35.6299	28.7523	24.9660	22.3726	20.2553


For example, if we desired to determine scan current in a 175mm water phantom that has 10 HU noise at 120 kV using the above table. The Noise is modeled using a 2-dimensional function of the current and the scan current. The scan current is then extracted from the 2D model. The noise is measured using a nominal slice width, typically 5.0 mm. The noise table will then be normalized based on the selected slice thickness of the scan protocol.

Appendix E Rotating Anode X-Ray Tube



GS-3073

Rotating Anode X-Ray Tube Tubes
Radiogènes à Anode Tournante
Drehanoden - Röntgenröhre
Tubos de Rayos-X con Anodo Giratorio



Note: Document originally drafted in the English language.
Note : Document à l'origine rédigé dans l'anglais.
Anmerkung: Dokument ursprünglich gezeichnet in der englischen Sprache.
Nota: Documento elaborado originalmente en la lengua inglesa.

Product Description The GS-3073 is a 5.5" (140 mm) 150 kV, 2.5 MJ (3.5 MHU) maximum anode heat content, rotating anode insert. This insert is specifically designed for CT Scanners. The insert features a 7° tungsten-rhenium facing on molybdenum with a graphite backed target and is available with the following nominal focal spot:	Description du Produit Le tube GS-3073, est une tube à anode tournante de plateau 140 mm, (5,5 pouces), 150 kV, d'une capacité thermique de 2,5 MJ (3,5 MUC). Il est à spécialement conçu pour une utilisation avec les scanners CT. Le pente de l'anode en molybdène traitée, tungstène, rhénium, recouverte de graphite, est de 7°. La dimension des foyers est de:	Produktbeschreibung Die GS-3073 ist eine 140 mm (5.5") Doppelfokus Drehanoden-Röntgenröhre, mit einer Anoden Wärmespeicherkapazität von 2,5 MJ (3,5 MHU) und einer max. Spannungsfestigkeit von 150 kV. Die Röntgenröhre wurde für den Einsatz an CT Scannern entwickelt. Der rückseitig graphitbeschichtete Wolfram Rhenium-Molybdän Anodenteller besitzt einen Winkel von 7°. Folgende Brennfleckkombination ist lieferbar:	Descripción del Producto El GS-3073 es un tubo de ánodo giratorio de 140 mm (5.5"), 150 kV, 2.5 MJ (3.5 kUC), la cual es el máximo almacenaje termal del ánodo, es diseñado específicamente para uso en CT scanners. El blanco emisor es una combinación de tungsteno, renio y molibdeno con grafito en la parte posterior con un rayo central de 7 grados. Disponible con las siguientes combinaciones de marcas focales:
0.7 x 0.8 1.2 x 1.4 IEC 60336 Loading Factor for slit focal: Small - 120 kV, 100 mA Large - 120 kV, 200 mA Maximum Anode Cooling Rate: 8,750 W (12,250 HU/sec) Maximum continuous anode heat dissipation: 3,400 W (4,760 HU/sec) Nominal Anode Input Power: Small - 23 kW IEC 60613 Large - 42 kW IEC 60613 Reference Axis: Perpendicular to port face. This insert is intended for use in a Varex Imaging B-240H housing.	0.7 x 0.8 1.2 x 1.4 CEI 60336 Facteur de charge pour foyer à fente: Petit - 120 kV, 100 mA Grand - 120 kV, 200 mA Taux maximum de refroidissement de l'anode: 8,750 W (12,250 UC/sec) Description calorifique maximum de l'anode (en continu): 3,400 W (4,760 UC/sec) Puissance Nominale de l'anode: Petit - 23 kW CEI 60613 Grand - 42 kW CEI 60613 Référence axe: Perpendiculaire à la face de sortie. Ce tube est essentiellement destiné à être employé dans les gaines Varex Imaging des séries B-240H.	0.7 x 0.8 1.2 x 1.4 IEC 60336 Ladefaktor: Klein - 120 kV, 100 mA Gross - 120 kV, 200 mA Nennleistung der Anode: 8,750 W (12,250 HU/sek) Maximale kontinuierliche Wärmeableitung des Anodentellers: 3,400 W (4,760 HU/sek) Nominale Anoden Eingangsleistung: Klein - 23 kW IEC 60613 Gross - 42 kW IEC 60613 Referenz Achsen: Senkrecht zum Strahlenaustrittsfenster. Die Röntgenröhre ist für den Einbau in die Varex Imaging Strahlerhaube B-240H vorgesehen.	0.7 x 0.8 1.2 x 1.4 IEC 60336 Carga Eléctrica Para la Abertura Focal: Pequeño - 120 kV, 100 mA Grande - 120 kV, 200 mA Medida Máxima del Enfriamiento del Anodo: 8,750 W (12,250 HU/seg) Máxima disipación termal continuo del Anodo: 3,400 W (4,760 HU/seg) El Poder de Penetración para el Anodo Nominal: Pequeño - 23 kW IEC 60613 Grande - 42 kW IEC 60613 Referencia de axes: Perpendicular a la abertura facial. Este tubo es diseñado, para uso en los encajes Varex Imaging de la serie B-240H.

133595-000 Rev A 01/17

Copyright © 2017, Varex Imaging Corporation. All Rights Reserved.



GS-3073

Volumetric / Helical Scan Ratings IEC 60613
 Tableaux des Caractéristiques Nominale de Balayage Volumétrique/Hélicoïdal CEI 60613
 Volumen-/Spiralbelichtungs-Leistungdiagramme IEC 60613
 Volumétrico/Clasificación Gráfica del Escán/Helicoideo IEC 60613

3Ø 60 Hz

0.7 x 0.8 Focal Spot 7 Degrees
 0.7 x 0.8 Dimension Focale 7 Degrés
 0.7 x 0.8 Brennfleck 7 Grad
 0.7 x 0.8 De Marcas Focales 7 Grados

Volume Scan Time (Seconds)	MAXIMUM ALLOWED TUBE CURRENT (mA) AS A FUNCTION OF THE FOLLOWING STARTING HEAT STORAGE AND TUBE VOLTAGES								
	Starting H.S. = 16%			Starting H.S. = 33%			Starting H.S. = 50%		
	120 kV	130 kV	140 kV	120 kV	130 kV	140 kV	120 kV	130 kV	140 kV
1	125	100	100	125	100	100	125	100	100
2	125	100	100	125	100	100	125	100	100
4	125	100	100	125	100	100	125	100	100
10	125	100	100	125	100	100	125	100	100
20	125	100	100	125	100	100	125	100	100
30	125	100	100	125	100	100	125	100	100
40	125	100	100	125	100	100	125	100	100
50	125	100	100	125	100	100	125	100	100
60	125	100	100	125	100	100	125	100	100
70	125	100	100	125	100	100	100 (a)	100 (a)	100 (a)

3Ø 180 Hz

0.7 x 0.8 Focal Spot 7 Degrees
 0.7 x 0.8 Dimension Focale 7 Degrés
 0.7 x 0.8 Brennfleck 7 Grad
 0.7 x 0.8 De Marcas Focales 7 Grados

Volume Scan Time (Seconds)	MAXIMUM ALLOWED TUBE CURRENT (mA) AS A FUNCTION OF THE FOLLOWING STARTING HEAT STORAGE AND TUBE VOLTAGES								
	Starting H.S. = 16%			Starting H.S. = 33%			Starting H.S. = 50%		
	120 kV	130 kV	140 kV	120 kV	130 kV	140 kV	120 kV	130 kV	140 kV
1	175	175	150	175	175	150	175	175	150
2	175	175	150	175	175	150	175	175	150
4	175	175	150	175	175	150	175	175	150
10	175	175	150	175	175	150	175	175	150
20	175	175	150	175	175	150	175	175	150
30	175	175	150	175	175	150	175	175	150
40	175	175	150	175	175	150	175	175	150
50	175	175	150	175	175	150	150 (a)	150 (a)	125 (a)
60	175	175	150	175	175	150	125 (a)	125 (a)	100 (a)
70	175	175	150	150 (a)	150 (a)	125 (a)	100 (a)	100 (a)	100 (a)

3Ø 60 Hz

1.2 x 1.4 Focal Spot 7 Degrees
 1.2 x 1.4 Dimension Focale 7 Degrés
 1.2 x 1.4 Brennfleck 7 Grad
 1.2 x 1.4 De Marcas Focales 7 Grados

Volume Scan Time (Seconds)	MAXIMUM ALLOWED TUBE CURRENT (mA) AS A FUNCTION OF THE FOLLOWING STARTING HEAT STORAGE AND TUBE VOLTAGES								
	Starting H.S. = 16%			Starting H.S. = 33%			Starting H.S. = 50%		
	120 kV	130 kV	140 kV	120 kV	130 kV	140 kV	120 kV	130 kV	140 kV
1	250	250	225	250	250	225	250	250	225
2	250	250	225	250	250	225	250	250	225
4	250	250	225	250	250	225	250	250	225
10	250	250	225	250	250	225	250	250	225
20	250	250	225	250	250	225	250	225	225
30	250	250	225	250	250	225	225	225	200
40	250 (b)	225 (b)	200 (b)	250 (b)	225 (b)	200 (b)	175 (a)	175 (a)	150 (a)
50	225	225	200	225 (a)	200 (a)	175 (a)	150 (a)	150 (a)	125 (a)
60	225	200	200	175 (a)	175 (a)	150 (a)	125 (a)	125 (a)	100 (a)
70	200 (a)	200 (a)	175 (a)	150 (a)	150 (a)	125 (a)	100 (a)	100 (a)	100 (a)

3Ø 180 Hz

1.2 x 1.4 Focal Spot 7 Degrees
 1.2 x 1.4 Dimension Focale 7 Degrés
 1.2 x 1.4 Brennfleck 7 Grad
 1.2 x 1.4 De Marcas Focales 7 Grados

Volume Scan Time (Seconds)	MAXIMUM ALLOWED TUBE CURRENT (mA) AS A FUNCTION OF THE FOLLOWING STARTING HEAT STORAGE AND TUBE VOLTAGES								
	Starting H.S. = 16%			Starting H.S. = 33%			Starting H.S. = 50%		
	120 kV	130 kV	140 kV	120 kV	130 kV	140 kV	120 kV	130 kV	140 kV
1	375	350	325	375	350	325	375	350	325
2	375	350	325	375	350	325	375	350	325
4	375	350	325	375	350	325	375	350	325
10	350 (b)	300 (b)	300 (b)	350 (b)	300 (b)	300 (b)	325 (b)	300 (b)	300 (b)
20	350 (b)	300 (b)	300 (b)	350 (b)	300 (b)	300 (b)	325	300	275
30	300 (b)	275 (b)	250 (b)	300 (b)	275 (b)	250 (b)	250 (a)	225 (a)	200 (a)
40	250 (b)	225 (b)	200 (b)	250 (b)	225 (b)	200 (b)	175 (a)	175 (a)	150 (a)
50	250 (b)	225 (b)	200 (b)	225 (a)	200 (a)	175 (a)	150 (a)	150 (a)	125 (a)
60	250 (b)	225 (b)	200 (b)	175 (a)	175 (a)	150 (a)	125 (a)	125 (a)	100 (a)
70	200 (a)	200 (a)	175 (a)	150 (a)	150 (a)	125 (a)	100 (a)	100 (a)	100 (a)

Note:

1. Limits are based on maximum track rating except for the following codes:
 a - Limited by available heat storage.
 b - Limited by window heating.
 c - Limited by filament emission.
 2. H.S. = Heat Storage
 kV = Tube Voltage

Remarque:

1. Les limites sont fonction de l'indice maximal de surface de l'anode, sauf pour les codes suivants:
 a - Limite par le stockage thermique disponible.
 b - Limite par le chauffage de la fenêtre.
 c - Limite par le rayonnement des filaments.
 2. H.S. = Stockage Thermique
 kV = Tube Voltage

Anmerkungen:

1. Grenzwerte basieren auf der maximalen Anodenoberflächenleistung mit Ausnahme der folgenden Codes:
 a - Durch verfügbare Wärmekapazität begrenzt.
 b - Durch Öffnungserwärmung begrenzt.
 c - Durch Glühfadenemission begrenzt.
 2. H.S. = Wärmekapazität
 kV = Röhre Spannung

Nota:

1. La clasificación de la marca máxima son limitadas, excepto por los siguientes códigos:
 a - Limitado por el almacenaje de calor disponible.
 b - Limitado por el calor de conducción de la ventana.
 c - Limitado por la emisión del filamento.
 2. H.S. = Almacenaje de calor
 kV = Tubo Voltage

Note:

Rating charts reflect maximum tube performance. Tube operation is ultimately limited by system software.

Remarque:

Abaque de caractéristiques représentent des valeurs maximales. L'utilisation du tube est finalement limitée par le logiciel du système.

Anmerkungen:

Die Leistungsdiagramme reflektieren die maximale Röhrenleistung. Der Röhrenbetrieb ist ultimativ zu begrenzen durch die Systemkontrollsoftware.

Nota:

El máximo poder del tubo es reflejada en el clasificación diagrama. La operación del tubo es ultimamente limitada por el control del sistema programado.



Copyright © 2017, Varex Imaging Corporation. All Rights Reserved.



GS-3073

3 Ø ==

Cathode Emission Characteristics Charts IEC 60613
 Caractéristiques d'Émission du Filament CEI 60613
 Kathoden - Emissionskennlinien IEC 60613
 Características de Emisión del Catodo IEC 60613



Copyright © 2017, Varex Imaging Corporation. All Rights Reserved.



GS-3073

B-240H Housing

Le Gaine B-240H

Das B-240H Gehäuse

Encaje de B-240H


Maximum Peak Voltage	150 kV
Anode to Ground	75 kV
Cathode to Ground	75 kV
Maximum X-ray Tube Assembly Heat Content	1.5 MJ (2.0 MHU)
Maximum Continuous Heat Dissipation (Includes stator heat)	3.7 kW (5.18 kHU/sec)
Maximum Heat Exchanger Dissipation	5.0 kW (7.0 kHU/sec)
Focal Point Position (Central Ray) Within 1mm (X, Y Direction from the center of radiation port.)	
X-Ray Tube Assembly	
Permanent filtration	1.0 mm Al IEC 60522
Loading Factors for Leakage Radiation	150 kV, 23 mA
Federal Standard High Voltage Cable	72
Ambient Air Temperature Limits for Operation	5°C to 40°C
Temperature Limits for Storage and Transport	-20°C to +75°C
Humidity	+10% to +90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Weight: Housing	33.1 kg (73 lbs)
Heat Exchanger	15.4 kg (34 lbs)
IEC Classification	Class 1
Safety Devices: (Internal) Thermal Switch	
Normally Closed Contact	Opening at 85°C ±4°C Closes at 74°C ±3°C
Flow Switch - Normally Open contact	Contacts close with adequate oil flow.
Heat Exchanger Control (Internal) Thermal Switch	
Normally closed contact	Opening at 70°C ±4°C Closes at 59°C ±3°C
Filament Frequency Limits	50 HZ - 25 kHz
Power Supply	DC

Voltage Maximum	150 kV
Tension Anode - Terre	75 kV
Tension Cathode - Terre	75 kV
Capacité Thermique Maximale de L'Ensemble Tube/Gaine ..	1.5 MJ (2.0 MUC)
Dissipation thermique continue de la gaine (Inclut la chaleur statorique)	3.7 kW (5.18 kUC/sec)
Dissipation Maximale de l'échangeur de chaleur	5.0 kW (7.0 kUC/sec)
Position du foyer (rayon central) à 1mm près (Coordonnées X, Y par rapport au centre du port de rayonnement.)	
Ensemble Radiogène	
Filtre non amovible	1,0 mm Al CEI 60522
Facteur de Charge Poru Rayonnement de fuite	150 kV, 23 mA
Embout de Cables au Standard Federal	72
Température Ambiante Pendant L'usage	5°C à 40°C
Limites de Température Pour le Transport et Pour L'Emmasinage	
Humidité	-20°C à +75°C
Humidité	+10% à +90%
Limites de pression atmosphérique	70 kPa à 106 kPa
Poids: Gaine	33,1 kg (73 lbs)
Échangeur de Chaleur	15,4 kg (34 lbs)
Classification CEI	Classe 1
Dispositifs de Sécurité: Switch Thermique	
Normalement Fermé	Ouverture à 85°C ±4°C Se ferme à 74°C ±3°C
Contacteur de débit - Contact Normalement Ouverture	Contacts fermés en présence d'un débit d'huile adéquat.
Contrôle d'échangeur de chaleur (Interne) Switch Thermique	
Normalement Fermé	Ouverture à 70°C ±4°C Se ferme à 59°C ±3°C
Limites de Fréquence des Filaments	50 HZ - 25 kHz
Alimentation Demandée	Courant Continu

Maximale Spannungsfestigkeit	150 kV
Anode gegen Erde	75 kV
Kathode gegen Erde	75 kV
Maximale Wärmespeicherkapazität des Strahlergehäuse ..	1.5 MJ (2.0 MHU)
Maximale kontinuierliche Wärmeableitung des Strahlergehäuse (einschließlich Statorwärme)	3.7 kW (5.18 kHU/sek)
Maximale Wärmeaustauscher - Verlustleistung	5.0 kW (7.0 kHU/sek)
Brennfleckposition (Zentralstrahl) Innerhalb 1mm. (X-, Y-Achse von der Mitte des Strahleneintrittsfensters)	
Röntgenstrahlers	
Eigenfilterwert	1.0 mm Al IEC 60522
Ladefaktoren für Leckstrahlungsmessung	150 kV, 23 mA
Federal Standard Hochspannungsbuchsen	72
Umgebungstemperaturgrenzen für den Betrieb	5°C zu 40°C
Temperaturgrenzen für Aufbewahrung und Transport	-20°C zu +75°C
Feuchtigkeit	+10% zu +90%
Luftdruck	70 kPa zu 106 kPa
Gewicht: Gehäuse	33.1 kg (73 lbs)
Wärmetauscher	15.4 kg (34 lbs)
IEC Klassifizierung	Klasse 1
Sicherheitseinrichtungen: Thermoschalter	
normalerweise geschlossen Verbindung	Offen bei 85°C ±4°C Schließt an 74°C ±3°C
Strömungsschalter - Kontakte normalerweise offen	Kontakte schließen sich bei ausreichendem Ölfluß.
Wärmetauscher-Steuerung(Interne) Thermoschalter	
normalerweise geschlossen Verbindung	Offen bei 70°C ±4°C Schließt an 59°C ±3°C
Heizfaden Frequenzgrenze	50 HZ - 25 kHz
Netzanschluß	DC

Voltage Maximo Elevado	150 kV
Anodo a Tierra	75 kV
Catodo a Tierra	75 kV
Maximo Calor Contenido Ensamblaje del Tubo de Rayos X ..	1.5 MJ (2.0 MHU)
Difusion del calor continuo del encaje (Incluye el calor de la bovina)	3.7 kW (5.18 kHU/seg)
Disipación maxima del radiador	5.0 kW (7.0 kHU/seg)
Posición de la marca focal (Rayo Central) Dentro de 1mm. (La Dirección axial X, Y se refiere del centro de la Radiación Portal.)	
Ensamblaje de Tubo de Rayos X	
Filtración Permanente	1.0 mm Al IEC 60522
Especificaciones de Encaje para la fuga de Radiación	150 kV, 23 mA
Cable de Receptáculos Comun Federal	72
Temperatura Limitada de Operación	5°C a 40°C
Temperatura Limitada de Almacen y Transporte	-20°C a +75°C
Humedad	+10% a +90%
Limites de la presión atmosférica	70 kPa a 106 kPa
Peso: Encaje	33.1 kg (73 lbs)
Radiador	15.4 kg (34 lbs)
IEC Clasificación	Clase 1
Aparatos de Seguridad: Interruptor Termal	
Normalmente Cerrado	Apertura en 85°C ±4°C Se cierra en 74°C ±3°C
Interruptor de Flujo - Normalmente los contactos setan abiertos	Contactos cerrado con a decuado flujo de aceite.
Control del Radiador (Interno) Interruptor Termal	
Normalmente Cerrado	Apertura en 70°C ±4°C Se cierra en 59°C ±3°C
Limites de la frecuencia del filamento	50 HZ - 25 kHz
Suministrador-de-Poder	Corriente Directa

Copyright © 2017, Varex Imaging Corporation. All Rights Reserved.






GS-3073

Tube Housing Assembly Heating and Cooling IEC 60613
 Échauffement et Refroidissement de l'Ensemble CEI 60613
 Röhregehäusebaugruppe Aufheizung und Abkühlung IEC 60613
 Enfriamiento y Calentamiento del Encaje Asamblado IEC 60613

X-RAY TUBE ASSEMBLY HEATING AND COOLING CURVE



NOTE:
 1. Heat inputs into housing include tube power, filament power, and stator power.
 2. Heating curves based on no restrictions of natural convection around tube housing assembly.
 3. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

Remarque:
 1. L'apport calorifique dans la gaine inclut la puissance du tube, du filament et du stator.
 2. Courbes d'échauffement basées sur une circulation d'air naturelle sans entrave autour de l'ensemble gaine-tube.
 3. Les abaques d'échauffement et de refroidissement représentent des valeurs maximales. L'utilisation du tube est finalement limitée par le logiciel du système.

Anmerkungen:
 1. Der Erwärmungskurven berücksichtigen die Verlustleistung aus der Anode, der Kathode und des Stators.
 2. Die Heizkurven basieren auf keinerlei Einschränkung der natürlichen Konvektion in der Umgebung der Strahlerhaube.
 3. Die Angaben stellen die höchstzulässigen Betriebswerte dar. Der technische Betrieb muß im Rahmen der Belastungs- und Abkühlkennlinien erfolgen.

NOTA:
 1. La energía del encaje incluye el poder del tubo, el poder del filamento y el poder de la bobina.
 2. Las curvas de calentamiento no son afectadas por el calor natural creado en la parte exterior del encaje.
 3. El máximo poder del tubo es reflejada en el diagrama de enfriamiento y calentamiento del tubo es ultimamente limitada por el control del sistema programado.


Copyright © 2017, Varex Imaging Corporation. All Rights Reserved.



GS-3073

Stator - Wiring Diagram
Stator - Schéma de Câblage
Stator - Drahtfarbentabelle
Bovina - Diagramas

Stator Ratings and Characteristics
Spécificités et Caractéristiques du Stator
Statornennleistungen und Merkmale
Características y Clarificación de la Bovina



Wire Color Couleurs des Branchements Kabelfarben Cable de Color	Description Description Beschreibung Descripción
1 Red	Phase Shift
Rouge	Stator de Changement de Phase
Rot	Veränderliche Statorphase
Rojo	Cambio de Fase del Estator
2 Black	Phase
Noir	Phase
Schwarz	Phase
Negro	Fase
3 White	Common
blanc	Neutre
Weiss	Neutral
blanco	Común
4/5 Green	Housing Ground
Vert	Masse de la cage
Grün	Masse des Gehäuse
Verde	Encaje a Tierra
6 Orange	Thermal Switch
Orange	Switch Thermique
Orange	Thermoschalter
Anaranjado	Interruptor termal
7 Red	Thermal Switch
Rouge	Switch Thermique
Rot	Thermoschalter
Rojo	Interruptor termal
8 White	Thermal Switch
blanc	Switch Thermique
Weiss	Thermoschalter
blanco	Interruptor termal
9 Black	Thermal Switch
Noir	Switch Thermique
Schwarz	Thermoschalter
Negro	Interruptor termal

Stator Drive Frequency Fréquence d'alimentation du stator Statorantriebs Frequenz Frecuencia de la Impulsión del estator	RPM
50 Hz	2900 - 3000
60 Hz	3400 - 3600
150 Hz	8500 - 9000
180 Hz	9500 - 10,800


Stator Type: "R"	
Stator Coil Resistance:	
Black to White	14.0 Ohms ±15%
Green to White	46.0 Ohms ±15%
Starter Voltage:	
50/60 Hz	200 VAC
150/180 Hz	470 VAC
Time to Full Speed:	
50/60 Hz	0 - 3000 RPM 9.6 Sec.
150/180 Hz	0 - 8000 RPM 9.6 Sec.
X-Ray Tube Assembly:	
GS-3073/B-240H	IEC 60601-2-28

Genre Stator: "R"	
Résistance de la bobine du stator:	
(résistance ohmique)	
Noir - Blanc	14,0 Ohms ±15%
Vert - Blanc	46,0 Ohms ±15%
Tension de démarrage:	
50/60 Hz	200 alternatif au démarrage
150/180 Hz	470 alternatif au démarrage
Temps pour atteindre la vitesse maximum:	
50/60 Hz de 0 à 3000 trs./min	9,6 Sec.
150/180 Hz de 0 à 8000 trs./min	9,6 Sec.
Ensemble radiogène:	
GS-3073/B-240H	CEI 60601-2-28

Statortyp: "R"	
Stator - Spulenwiderstand	
Schwarz - Weiss	14,0 Ohms ±15%
Grün - Weiss	46,0 Ohms ±15%
Spannungen:	
50/60 Hz	200 VAC
150/180 Hz	470 VAC
Hochlaufzeit:	
50/60 Hz	0 - 3000 U/min 9,6 Sek.
150/180 Hz	0 - 8000 U/min 9,6 Sek.
Röntgenstrahler:	
GS-3073/B-240H	IEC 60601-2-28

Tipo de la Bovina: "R"	
Resistencia del Rollo de la Bovina:	
Negro a Blanco	14,0 Ohms ±15%
Verde a Blanco	46,0 Ohms ±15%
Voltage de la Obtenida:	
50/60 Hz	200 VAC
150/180 Hz	470 VAC
Tiempo Para la Velocidad Maxima:	
50/60 Hz	0 - 3000 RPM 9,6 Segundo
150/180 Hz	0 - 8000 RPM 9,6 Segundo
Ensamblaje de Tubos de Rayos X:	
GS-3073/B-240H	IEC 60601-2-28

Copyright © 2017, Varex Imaging Corporation. All Rights Reserved.



Appendix F Error Code

Table 69: Error code list

Error code	Error code description	Popup description	Cause
0	ABORT_EXAM	Exam has been aborted. Please try again. If problem persists, contact customer service with code:	Generic message whenever a scan/exam has been terminated abnormally. User would not typically see this because a more specific error should be posted.
1	CAN_DEVICE_DISCO NNECT_CENT1	Communications fault. Please contact customer service with code:	Cannot communicate with centipede 1 device.
2	CAN_DEVICE_DISCO NNECT_CENT2	Communications fault. Please contact customer service with code:	Cannot communicate with centipede 2 device.
3	CAN_DEVICE_DISCO NNECT_ROTATE	Communications fault. Please contact customer service with code:	Cannot communicate with rotate device.
4	CAN_DEVICE_DISCO NNECT_BIB	Communications fault. Please contact customer service with code:	Cannot communicate with BIB device.
5	CAN_DEVICE_DISCO NNECT_OIB1	Communications fault. Please contact customer service with code:	Cannot communicate with OIB1 device. Note that this alone will not cause a Fault state but Start and Cancel buttons on one side of scanner will not operate.
6	CAN_DEVICE_DISCO NNECT_OIB2	Communications fault. Please contact customer service with code:	Cannot communicate with OIB2 device. Note that this alone will not cause a Fault state but Start and Cancel buttons on one side of scanner will not operate.

Error code	Error code description	Popup description	Cause
7	CAN_DEVICE_DISCONNECTION_POWER	Communications fault. Please contact customer service with code:	Cannot communicate with CCB device.
8	CAN_DEVICE_DISCONNECTION_TRANS	Communications fault. Please contact customer service with code:	Cannot communicate with Transport device.
9	CAN_DEVICE_DISCONNECTION_DCB	Communications fault. Please contact customer service with code:	Cannot communicate with DCB device.
10	CAN_DEVICE_DISCONNECTION_HVG	Communications fault. Please contact customer service with code:	Cannot communicate with HVG device.
11	HVG_LATCH_ERROR_ENABLE	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	There was a problem with x-ray Enable signal.
12	HVG_LATCH_ERROR_INTERLOCK	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	An interlock condition was asserted.
13	HVG_LATCH_ERROR_110_TIMER	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	The 110% timer has expired and forced x-rays off.
14	HVG_LATCH_ERROR_XRT_THERM_SW	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	X-ray tube thermal switch asserted.
15	HVG_LATCH_ERROR_HE_FLOW_SW	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Heat exchanger flow switch asserted.
16	HVG_LATCH_ERROR_WDT	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	The watchdog timer has timed out and forced x-rays off (WD timer is controlled by the DCB firmware), this error would be unusual.

Error code	Error code description	Popup description	Cause
17	HVG_LATCH_ERROR_ARC_FAULT	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Arcs occurred.
18	HVG_LATCH_ERROR_HVG_FAULT	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	A generic HVG fault condition occurred (look at HVG_ERROR_code).
19	HVG_LATCH_ERROR_STARTER	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Anode (starter) device reported and error.
20	HVG_LATCH_ERROR_DAS_OVER_RANGE	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	DAS (detector) data values are out of range.
21	HVG_ERROR_MA_REGULATION	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	There was a problem with mA regulation.
22	HVG_ERROR_KV_REGULATION	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	There was a problem with kV regulation.
23	HVG_ERROR_ANODE_STARTER	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Anode (starter) device reported and error.
24	HVG_ERROR_INVERTER_TEMP	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Inverter over temp condition.
25	HVG_ERROR_UNCOMMAND_EXP	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	An exposure happened without being commanded.
26	HVG_ERROR_ANODE_OVER_VOLTAGE	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Anode over voltage condition.

Error code	Error code description	Popup description	Cause
27	HVG_ERROR_CATHODE_OVER_VOLTAGE	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Cathode over voltage condition.
28	HVG_ERROR_ANODE_OVER_CURRENT	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Anode over current condition.
29	HVG_ERROR_CATHODE_OVER_CURRENT	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Cathode over current condition.
30	HVG_ERROR_FILAMENT_OVER_CURRENT	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Filament over current condition.
31	HVG_ERROR_ARC_DETECTED	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Arcs occurred.
32	HVG_ERROR_CURRENT_RETURN_WIRE_DISCONNECT	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	Current return wire disconnect.
33	HVG_ERROR_MAOVER_PROG	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	The actual mA was higher than programmed.
34	HVG_ERROR_KV_OVER_PROG	High voltage fault. Please retry task. If problem persists, please contact customer service with code:	The actual kV was higher than programmed.
35	HVG_ERROR_FILAMENT_REGULATION	High voltage failure, retry protocol, if problem persists contact customer service with code	There was a problem with filament regulation.
36	DCB_ERROR_BAD_DETECTORS	Data acquisition fault. Please contact customer service with code:	DCB reported a bad detector condition.

Error code	Error code description	Popup description	Cause
40	RECON_PROTOCOL_REJECTED_INVALID_PROTOCOL_TYPE_OR_USAGE_PARAMETER	Recon Protocol Rejected due to invalid Type or Usage parameter	Recon rejected protocol because of invalid parameters.
41	RECON_PROTOCOL_REJECTED_INVALID_PROTOCOL	Recon Protocol Rejected due to invalid Protocol	Recon rejected protocol because of invalid parameters.
42	RECON_PROTOCOL_REJECTED_INVALID_PROTOCOL_IMAGE_COORDINATES	Recon Protocol Rejected due to invalid Image Coordinates	Recon rejected protocol because of invalid parameters.
43	RECON_PROTOCOL_REJECTED_INVALID_PROTOCOL_ROI_COORDINATES	Recon Protocol Rejected due to invalid ROI Coordinates	Recon rejected protocol because of invalid parameters.
44	RECON_PROTOCOL_REJECTED_HELICAL_QA_FAILED	Recon Protocol Rejected due to Helical QA Failure	Recon rejected protocol because of invalid parameters.
45	RECON_PROTOCOL_REJECTED_RECON_BUSY	Recon Protocol Rejected due to Recon Busy	Recon is in an invalid state to perform a protocol.
46	RECON_PROTOCOL_REJECTED_SERIAL_LINK_DISCONNECT_OCCURRED	Recon Protocol Rejected due to Serial Link Disconnect	Recon rejected protocol because the serial link is not connected.
47	RECON_PROTOCOL_REJECTED_INSUFFICIENT_MEMORY	Recon Protocol Rejected due to Insufficient Memory	Recon rejected protocol because of insufficient memory necessary to perform requested protocol.
48	RECON_PROTOCOL_REJECTED_INVALID_PROTOCOL_SLICE_COORDINATES	Recon Protocol Rejected due to invalid Slice Coordinates	Recon rejected protocol because of invalid parameters.
50	RECON_DATA_CORRUPTED	System has encountered an unexpected error. Please contact customer service with code:	The Recon app reported data corruption (view headers not correct).

Error code	Error code description	Popup description	Cause
51	RECON_OFFSET_CAL_FAILED	System failed to properly perform protocol calibration. Please try again. If problem persists, please contact customer service with code:	The Recon app reported Offset Cal failure.
52	RECON_AIR_CAL_FAILED	System failed to properly calibration system. Please try calibration again. If problem persists, please contact customer service with code:	The Recon app reported Air Cal failure.
53	RECON_OFFSET_CAL_TIMEOUT	System failed to properly perform protocol calibration. Please try again. If problem persists, please contact customer service with code:	Scanner control app timed out waiting for an offset cal to complete.
54	RECON_OFFSET_PROTOCOL_TIMEOUT	System failed to properly perform protocol calibration. Please try again. If problem persists, please contact customer service with code:	Scanner control app timed out waiting for Recon app to acknowledge a protocol request (offset or image).
55	RECON_PROTOCOL_TIMEOUT	System is unable to perform protocol. Please try again. If problem persists, please contact customer service with code:	Scanner control app timed out waiting for a scan to complete.
56	RECON_PROTOCOL_REJECTED	System is unable to perform protocol. Please try again. If problem persists, please contact customer service with code:	The Recon app rejected a scan protocol request.

Error code	Error code description	Popup description	Cause
57	RECON_POST_NO_SCAN_INFO	System is unable to perform post reconstruction. Please try again. If problem persists, please contact customer service with code:	The image data was not found for a Post Recon or Resend Images command (where image data could be scan_info file, .dcm files, .cor files, or entire directory).
58	RECON_POST_RECON_FAILED	System is unable to perform post reconstruction. Please try again. If problem persists, please contact customer service with code:	A Post Recon or Resend Images command could not be run due to scanner state.
60	PREPARE_FAILED	The scanner encountered a fault preparing for scan. Please contact customer service with code:	A generic Prepare error occurred. This should be accompanied by a more specific error. See scanner log for detail.
61	DISK_PREPARE_ERROR	Communications fault. Please contact customer service with code:	Disk subsystem reported an error during prepare.
62	CENTIPEDE_MOVE_TIMEOUT	The scanner failed to move to the correct location. Please check for obstructions or debris on the floor that could impede the movement and try again. If problem persists, please contact customer service with code:	Scanner Control app timed out waiting for scanner to reach requested position; OR scanner did not reach required velocity for a scout or helical scan.
63	ANODE_STARTUP_TIMEOUT	X-ray power sequence fault. Please contact customer service with code:	Scanner Control app timed out waiting for anode device to report anode rotation is at speed.

Error code	Error code description	Popup description	Cause
64	DISK_NOT_IN_TICK_MODE	Rotational speed fault. Please contact customer service with code:	Scanner Control app timed out waiting for disk to get to tick mode (during prepare for scans that specify disk rotation).
65	COLLIMATOR_POSITION_TIMEOUT	Failure to position collimator. Please contact customer service with code:	Scanner Control app timed out waiting for collimator to get to requested position.
66	DCB_READY_TIMEOUT	DCB communication fault. Please contact customer service with code:	Scanner Control app timed out waiting for DCB to report "ready" state.
67	HEAT_EXCHANGER_ERROR	X-ray cooling fault. Please contact customer service with code:	Heat exchanger did not come on during prepare.
68	FILAMENT_ERROR	X-ray filament fault. Please contact customer service with code:	Scanner Control app timed out waiting for DCB to report "filament on"; OR for filament monitor to report > 10.
69	SERIAL_LINK_NOT_UP	Communications fault. Please contact customer service with code:	Scanner Control app timed out waiting for serial link to come up; OR serial link was lost before start Acq.
70	ROTATE_COMMAND_FAILED	Rotation communication error. Please contact customer service with code:	Unused error code.
71	TRANSLATE_COMMAND_FAILED	Translate fault. Please contact customer service with code:	Move command to centipede device failed.
72	COULD_NOT_START_SSP	Software initialization fault. Please contact customer service with code:	SSP software (Scanner State and Position GUI, the GUIs that are on the scanner displays) did not start up.

Error code	Error code description	Popup description	Cause
73	ROTATE_TO_ANGLE_FAILED	Rotational fault. Please contact customer service with code:	Scanner Control app timed out waiting for "Rotate to Angle" operation.
74	COULD_NOT_CLEAR_ROTATE_FAULT	Rotation fault. Please contact customer service with code:	The rotate device has reported a fault, and the fault cannot be cleared.
75	ENCODER_CONSISTENCY_FAULT	Encoder consistency fault. Please contact customer service with code:	Unused error code.
76	INTERLOCK_TEST_FAILED	Interlock Test Failed. Please ensure that e-stop is not engaged. If e-stop is not engaged, please contact customer service with code:	The interlock circuit either remained continuous after a Software Interlock was applied; OR the interlock circuit was broken when it was expected to be continuous.
80	DISK_SPACE_LOW	System disk space low, please contact customer service with code:	Unused error code.
81	DISK_SPACE_QUERY_FAILED	System disk space low, please contact customer service with code:	Unused error code.
82	NOT_ALL_BASE_DEVICES_PRESENT	Base communications fault. Please contact customer service with code:	Not all devices were found on Base CAN bus (or a device has become disconnected). This should be accompanied by another error code describing which device(s) disconnected. You can look at Scanner Control GUI for a status of which devices are present.

Error code	Error code description	Popup description	Cause
83	NOT_ALL_DISK_DEVICES_PRESENT	Disk communications fault. Please contact customer service with code:	Not all devices were found on Disk CAN bus (or a device has become disconnected). This should be accompanied by another error code describing which device(s) disconnected. You can look at Scanner Control GUI for a status of which devices are present.
90	CCB_BATTERY_OPERATIONAL	Battery system fault. Please contact customer service with code:	CCB device reported a (not) Operational Alarm (scanning not possible).
91	CCB_BATTERY_INTERLOCK	Battery system communication error. Please contact customer service with code:	CCB device reported an Interlock Alarm (scanning not possible).
92	CCB_BATTERY_MAIN_BREAKER	Circuit breaker has been tripped. Please reset and contact customer service with code:	Battery main breaker turned off. Scanner is running on wall power only. (Scanning not possible.)
93	CCB_LOW_BATTERY_ALARM	Low battery condition. Please charge system as soon as possible.	CCB device reported a Low Battery Alarm (scanning not possible).
94	CCB_DEAD_BATTERY_ALARM	Dead battery condition. The system is shutting down. Please charge system and report condition to customer service with code:	CCB device reported a Dead Battery Alarm (auto-shutdown of scanner is imminent) (scanning DEFINITELY not possible).
95	CCB_HIGH_BATTERY_ALARM	High battery condition has occurred. Battery charging has been disabled.	CCB device reported a High Battery Alarm

Error code	Error code description	Popup description	Cause
96	CCB_OVERCGARGED_BATTERY_ALARM	Please power down and unplug system and contact customer service immediately with code:	CCB device reported an Overcharged Battery Alarm. This is beyond a High Battery warning and is serious. The system should be turned off and unplugged immediately.
97	CCB_BATTERY_HIGH_TEMP_ALARM	Battery system fault. Please contact customer service with code:	CCB device reported a High Temp Alarm (scanning not possible).
98	CCB_BATTERY_OVER_TEMP_ALARM	Over temperature battery condition. Please power down and unplug system and contact customer service immediately with code:	CCB device reported an Over Temp Alarm. This is beyond a High Temp warning and is serious. The system should be turned off and unplugged immediately.
99	CCB_BATTERY_MEASUREMENT_ERROR	Battery system fault. Please contact customer service with code:	CCB device reported a Measurement Error Alarm (scanning not possible)
100	CCB_BATTERY_IMBALANCE_WARNING	Battery system fault. Please contact customer service with code:	One or more battery voltage levels are not the same as the others. No action required. System will try to correct itself.
110	RECON_PROTOCOL_REJECTED_INVALID_PROTOCOL_USAGE_PARAMETER	Reconstruction fault. Please contact customer service with code:	Recon received invalid usage parameter in protocol.
111	RECON_PROTOCOL_REJECTED_INVALID_NUMBER_OF_VIEWS_PARAMETER	Reconstruction fault. Please contact customer service with code:	Recon received invalid number of views in protocol.
112	RECON_PROTOCOL_REJECTED_RUN_DMA_SETUP	Reconstruction fault. Please contact customer service with code:	Recon failed to initialize DMA in preparation for scan.

Error code	Error code description	Popup description	Cause
113	RECON_PROTOCOL_REJECTED_UNDEFINED_USAGE	Reconstruction fault. Please contact customer service with code:	Recon received undefined usage parameter in protocol.
114	RECON_PROTOCOL_REJECTED_INVALID_RAW_DATA_REPLAY	Reconstruction fault. Please contact customer service with code:	Unused error code.
115	RECON_PROTOCOL_REJECTED_FILES_RETRIEVE_FAILED	Reconstruction fault. Please contact customer service with code:	Unused error code.
116	RECON_PROTOCOL_REJECTED_INVALID_PARAMETER_STRUCTURE_SIZE	Reconstruction fault. Please contact customer service with code:	Recon received incorrect structure size.
117	RECON_PROTOCOL_REJECTED_INVALID_PROTOCOL_FLASH_IO_CMD	Reconstruction fault. Please contact customer service with code:	Unused error code.
118	RECON_PROTOCOL_REJECTED_PREPARE_AND_PRIME_POST_RECON_USAGE	Reconstruction fault. Please contact customer service with code:	Unused error code.
119	RECON_PROTOCOL_REJECTED_INVALID_POST_RECON_STATE	Reconstruction fault. Please contact customer service with code:	Recon received a protocol while still processing previous protocol/post recon.
120	RECON_PROTOCOL_REJECTED_INVALID_MESSAGE_BODY_LENGTH	Reconstruction fault. Please contact customer service with code:	Recon received incorrect structure size.
121	RECON_PROTOCOL_REJECTED_INVALID_RELOAD_PARAMETER_FILES	Reconstruction fault. Please contact customer service with code:	Unused error code.
122	RECON_PROTOCOL_REJECTED_INVALID_UNSUPPORTED_COMMAND	Reconstruction fault. Please contact customer service with code:	Recon was sent an unsupported command from scanner control.
123	RECON_PROTOCOL_REJECTED_INVALID_HELICAL_FILTER_KERNEL_TYPE	Reconstruction fault. Please contact customer service with code:	Recon was sent an invalid helical filter kernel from workstation/scanner control.

Error code	Error code description	Popup description	Cause
124	RECON_PROTOCOL_REJECTED_INVALID_NUMBER_OF_HELICAL_IMAGES_FOR_WINDMILL	Reconstruction fault. Please contact customer service with code:	Recon received an invalid number of images for windmill.
125	RECON_PROTOCOL_REJECTED_GPU_FAILED_TO_START	Reconstruction fault. Please contact customer service with code:	Recon failed to initialize GPU during preparation for scan.
126	RECON_PROTOCOL_REJECTED_RECON_BUSY	Reconstruction fault. Please contact customer service with code:	Recon received a protocol while still processing previous protocol/post recon.
130	RECON_AIR_CAL_FAILED_NON_AIR_IMAGE	Reconstruction fault. Please contact customer service with code:	Air image above threshold for air calibration.
131	RECON_AIR_CAL_FAILED_SEND_EVENT	Reconstruction fault. Please contact customer service with code:	Air calibration failed to be performed.
132	RECON_AIR_CAL_IMAGE_EXCEEDS_THRESHOLD	Reconstruction fault. Please contact customer service with code:	Air image above threshold for air calibration.
133	RECON_AIR_CAL_FAILED_NO_VIEW_DATA	Reconstruction fault. Please contact customer service with code:	No view data received during an air calibration.
134	RECON_AIR_CAL_FAILED_CORRUPTED_VIEW_DATA	Reconstruction fault. Please contact customer service with code:	Corrupted views received during air calibration.
135	RECON_OFFSET_CAL_FAILED_SEND_EVENT	Reconstruction fault. Please contact customer service with code:	Offset calibration failed to be performed.
136	RECON_OFFSET_CAL_FAILED_NO_VIEW_DATA	Reconstruction fault. Please contact customer service with code:	No view data received during an offset calibration.
137	RECON_OFFSET_CAL_FAILED_CORRUPTED_VIEW_DATA	Reconstruction fault. Please contact customer service with code:	Corrupted views received during offset calibration.
138	RECON_OFFSET_CAL_FAILED_BAD_REFERENCE	Reconstruction fault. Please contact customer service with code:	Offset calibration failed due to bad reference.

Error code	Error code description	Popup description	Cause
139	RECON_OFFSET_CAL_80_PERCENT_BAD_REFERENCE	Please contact Customer Service immediately and run a Quality Assurance (QA) Phantom test to verify image quality. Error code:	80% of reference detector values are above the acceptable threshold during an offset cal.
140	UPS_LOW_BATTERY_ALARM	Workstation low battery condition. Please charge system as soon as possible.	UPS device reported a Low Battery Alarm.
141	UPS_DEAD_BATTERY_ALARM	Workstation dead battery condition. The Workstation is shutting down. Please charge cart and report condition to customer service with code:	UPS device reported a Dead Battery Alarm.
142	UPS_HIGH_BATTERY_ALARM	Workstation high battery condition has occurred. Cart battery charging has been disabled.	UPS device reported a High Battery Alarm.
143	UPS_OVERCHARGED_BATTERY_ALARM	Please power down and unplug workstation cart and contact customer service immediately with code:	UPS device reported an Overcharged Battery Alarm. This is beyond a High Battery warning and is serious. The cart should be turned off and unplugged immediately.
144	UPS_BATTERY_HIGH_TEMP_ALARM	Workstation battery system fault. Please contact customer service with code:	UPS device reported a High Temp Alarm.
145	UPS_BATTERY_OVER_TEMP_ALARM	Workstation over temperature battery condition. Please power down workstation, unplug cart and contact customer service immediately with code:	UPS device reported an Over Temp Alarm. This is beyond a High Temp warning and is serious. The cart should be turned off and unplugged immediately.
146	UPS_BATTERY_MEASUREMENT_ERROR	Workstation battery system fault. Please contact customer service with code:	UPS device reported a Measurement Error Alarm.

Error code	Error code description	Popup description	Cause
147	UPS_BATTERY_IMBALANCE_WARNING	Workstation battery system fault. Please contact customer service with code:	One or more battery voltage levels are not the same as the others. No action required. UPS will try to correct itself.
148	UPS_CHARGER_FAULT	Workstation battery system fault. Please contact customer service with code:	UPS device reported a Charger fault.

Appendix G Revision History

Table 70: Workstation (only) revision history before combined with BodyTom Elite User Manual (below)

Revision	ECO number	Effective date	Author	Changes
00	ECO-002019	None supplied	Keith Almeida	Releasing 1-00132-060 User Manual, WS CART NL4000
01	ECO-002266	None supplied	Keith Almeida	Added “Installing a UPS” instructions
02	ECO-002735	None supplied	Rich DeSalvo	Removed “Cart Assembly” instructions Updated TOC, Figures, Tables Updated cart info
N/A	N/A	N/A See Table 71, below.	Cynthia Crow	Consolidated the <i>BodyTom Workstation User Manual</i> with the <i>BodyTom User Manual</i> . See Table 71, below.

Table 71: Revision history

Revision	ECO number	Effective date	Author	Changes
00	ECO-001673	2012/01/27	Christofer Krueger/Mike Limoli/Ibrahim Bechwati	New Release
01	ECO-001877	2012/11/25	Christofer Krueger	Correct minor edits uncovered during translations' review and pilot installations. Also needed to add 3 rd edition and HC requirements
02	ECO-002003	2013/12/19	Christofer Krueger	Update to meet IEC 60601-1 3 rd ed, 60825-1, and CFDA Requirements. Error codes, additional warnings, focal spot, CTDI ₁₀₀ results, scatter values, product label location, etc.
03	ECO-001992	2014/04/03	Keith Almeida	Updated to software release 01.06
04	ECO-002326	2015/01/28	Keith Almeida	Releasing NL4000 Software Version 01.07
05	ECO-002441	2015/04/22	Keith Almeida	Releasing NL4000 Software Version 01.07.01
06	ECO-002869	2016/04/18	Keith Almeida	Update for consistent formatting and Samsung branding. Added 1.07 Error codes Reduced use of pronouns Inclusion of Information Bulletins CSB-00006 and CSB-00007.
07	ECO-002915	2016/07/13	Keith Almeida	Addition of Slab, Kernel, and Floor Flatness Certification corresponding with version 1.08 of the software.
08	ECO-002961	2016/10/08	Ross Caisse	Removing references to bolus tracking for software release 1.08.01.
09	ECO-003015	2016/11/18	Ross Caisse	Adding references to bolus tracking for software release 1.08.02.

Revision	ECO number	Effective date	Author	Changes
10	ECO-003600	2017/09/14	Cynthia Crow	Changed branding from BodyTom to BodyTom Elite. Applied new template/styles, step-by-step procedures, and workflow for improved audience experience. Consolidated the previous (Revision 02) <i>BodyTom Workstation User Manual</i> (1-00132-060) with the <i>BodyTom User Manual</i> , Revision 10 (1-NL4000-060).
11	ECO-003712	2018/04/13	Derek Decoux	Updated the tube heat capacity icon description with correct ranges for different colors. Added Essential Performance factors to the preface. Added a Functional Earth Symbol (Table 55)
12	ECO-004280	2018/11/01	Christofer Krueger	Adding a note above Table 13 on EMC Characteristics. Modified Table 14 to detail the correct standard parameters. Changed title to include “NL4000 BodyTom”. Added “BodyTom” to the “About this Manual” section. Table 4 ~ 0.25Ω to 0.105Ω.
13	ECO-004754	2020/03/03	Christofer Krueger	Added details to comply with IEC 60825-1, 3 rd ed. Table 2 and Laser Safety section of Chapter 1. EU MDR updates: Added Serious Injury, Electronic Copy and Clinical Benefit notes to the preface

Revision	ECO number	Effective date	Author	Changes
14	ECO-004919	2020/06/17	Christofer Krueger	Added Laser warning statement from the FDA. Updated for new marking plate. Removed broken bookmarks. Corrected the workstation temp range. Modified an incomplete sentence in table 58. Updated glossary to remove the word "withring". Bug 5304. Added new statement and Figure 203 on Geometric efficiency for being less than 70%.
15	ECO-005032	2020/09/015	Stephen Dunn	Updated Geometric Efficiency wording, added ACR information
16	ECO-005301	2021/02/05	Gina Cunsolo	Updated to remove reference to BodyTom. Added symbols to Table 2 to reflect all symbols on the label.
17	ECO-005794	2022/02/23	Karen Reed/Keith Kaser/Stephen Dunn	Added statement to contact CS if there is an adverse event under Contact Info Sec. Added reference for eIFU in "About this user manual" Section Updated Symbols Table for clerical purposes. Updated Table 9 to reference both BodyTom and BodyTom Elite PNs. Updated WS Cart weight in Table 11 for correct weight. Removed old references of standards already mentioned under Compliance Statement. Added Table 42 (p. 244) and 43 (p. 245).

				<p>Updated images to remove older logos. Removed Floor calibration. Remove Surgical clamp details. Added Interventional scanning information. Updated Cleaning instructions to remove references from Chapter One and move to Cleaning Chapter 14. Updated Pendant information. Modified multiple bulleted tables to make them actual tables. Modified Dose Alert information. Added detailed AEC information to Appendix D. Updated entire document Font to Calibri 12. Modified Appendix C to include only a sample of the reference protocols. Modified all line spacing to 1.0. Updated International Customer Service contacts to make them more user friendly. Made many changes to wording and grammar throughout document. Updated all images to remove red callouts and arrows and replaced with gold. Removed small filament information from Focal Spot section of Chapter One. Removed E-Stop instructions from Chapter 1 and moved to Chapter 3. Added scanner Battery Capacity figure to Chapter 1. Added X-Ray tube capacity information to Chapter 1. Removed Floor flatness certification information from Chapter 1. Removed multiple icons/orbs from text of document. Removed redundant 'Saving a Preset' under Applying Windowing Presets in Chapter 5. Removed AEC information from 'Creating a new protocol' section</p>
--	--	--	--	---

Revision	ECO number	Effective date	Author	Changes
				in Chapter 6 since it is covered in detail in 'Scanning with special features' of Chapter 9. Added 'Editing and existing protocol' to Chapter 6. Updated 'Starting Quality Assurance' section of Chapter 7 to include better instructions for this process. Modified all contents of Chapter 10 to match appearance of items on the GUI. Updated all sections of Chapter 11 'Viewing Images' to make them clearer. Removed redundant 'Applying kernels in Post Reconstruction' section of Chapter 12. Updated 'Using Universal Transfer Board' in Chapter 13 to include more information. Updated Footers of pages 2-24. Added information to Site Specifications related to floor supporting product weight. Updated workstation weight. Updated Figure 5. Updated Figure 2, 12, 33, 38, 40 and 399 to remove BodyTom reference.
18	ECO-006198		Keith A. Kaser	Corrected multiple spacing issues. Added Brain CT Perfusion scanning instructions under the "Scanning with Special Features" section and moved Viewing Images in the CTP Viewer from the Viewing Images section to the new CT Perfusion section mentioned above.
19	ECO-006502		Keith A. Kaser	Updated 'Intended use of system' section. Replaced Figure 232 to change highlight section from red to gold. Replaced Figure 345 to change arrow from red to gold.

Revision	ECO number	Effective date	Author	Changes
				Updated Laser Section – changed numbers from Red to Green and Blue to match Figure 5 colors. Replaced Figure 343 to change highlight section from red to gold. Replaced Figure 392 to change highlight section from red to gold. Replaced Figure 398 to change arrows from red to gold. Updated Table 21 and 65 to correct direction of travel for Pendant buttons.

Copyright © applied and printed in 2022